

Datasheet for ABIN1048350

anti-Bestrophin 2 antibody (Internal Region)

2 Images



Overview

| Overview | |
|-----------------------|---|
| Quantity: | 50 μg |
| Target: | Bestrophin 2 (BEST2) |
| Binding Specificity: | Internal Region |
| Reactivity: | Human, Monkey, Dog |
| Host: | Rabbit |
| Clonality: | Polyclonal |
| Conjugate: | This Bestrophin 2 antibody is un-conjugated |
| Application: | Immunohistochemistry (IHC), Immunohistochemistry (Paraffin-embedded Sections) (IHC (p)) |
| Product Details | |
| Brand: | IHC-plus™ |
| Immunogen: | Synthetic 15 amino acid peptide from internal region of human BEST2 / Bestrophin-2. Percent identity with other species by BLAST analysis: Human, Gorilla, Monkey, Marmoset, Dog (100%), Mouse, Hamster, Panda, Bovine (93%), Opossum (87%), Elephant, Pig (80%). |
| | Type of Immunogen: Synthetic peptide |
| Specificity: | Human BEST2 / Bestrophin-2. BLAST analysis of the peptide immunogen showed no homology |
| | with other human proteins, except BEST1 (47 %). |
| Predicted Reactivity: | with other human proteins, except BEST1 (47 %). Percent identity with other species by BLAST analysis: Human, Gorilla, Monkey, Marmoset, Dog (100%) Mouse, Hamster, Panda, Bovine (93%) Opossum (87%) Elephant, Pig (80%). |

Target Details

| Target: | Bestrophin 2 (BEST2) |
|-------------------|---|
| Alternative Name: | BEST2 / Bestrophin-2 (BEST2 Products) |
| Background: | Name/Gene ID: BEST2 |
| | Subfamily: Bestrophin |
| | Family: Ion channel |
| | |
| | Synonyms: BEST2, Bestrophin-2, Bestrophin 2, Hbest2, VMD2-like gene 1, VMD2L1 |
| Gene ID: | 54831 |

Application Details

| Assay Procedure: | The IHC-pro Immunohistochemistry Protocol |
|--------------------|---|
| Comment: | Target Species of Antibody: Human |
| Application Notes: | Approved: IHC, IHC-P (2.5 - 5 μg/mL) |

Tissue Preparation

Formalin fixation and embedding in paraffin wax

Tissue Sectioning

Make 4-µm sections and place on pre-cleaned and charged microscope slides.

Heat in a tissue-drying oven for 45 minutes at 60°C

Deparaffinization

Wash slides in 3 changes of xylene – 5 minutes each at room temperature.

Rehydration

Wash slides in 3 changes of 100% alcohol – 3 minutes each at room temperature.

Wash slides in 2 changes of 95% alcohol – 3 minutes each at room temperature.

Wash slides in 1 change of 80% alcohol – 3 minutes at room temperature.

Rinse slides in gentle running distilled water – 5 minutes at room temperature.

Antigen retrieval

Steam slides in 0.01 M sodium citrate buffer, pH 6.0 at 99-100°C - 20 minutes

Remove from heat and let stand at room temperature in buffer - 20 minutes

Rinse in 1X TBS with Tween (TBST) – 1 minute at room temperature.

Immunostaining

Do not allow tissues to dry at any time during the staining procedure.

Apply a universal protein block – 20 minutes at room temperature.

Drain protein block from slides, apply diluted primary antibody – 45 minutes at room temperature.

Rinse slides in 1X TBST - 1 minute at room temperature.

Apply a biotinylated secondary antibody (specific to the host of the primary antibody) - 30 minutes at room temperature.

Rinse slides 1X TBST – 1 minute at room temperature.

Apply alkaline phosphatase streptavidin – 30 minutes at room temperature.

Rinse slides in 1X TBST - 1 minute at room temperature.

Apply alkaline phosphatase chromogen substrate - 30 minutes at room temperature.

Wash slides in distilled water – 1 minute at room temperature.

Dehydrate

This method should only be used if the chromogen substrate is alcohol insoluble.

Wash slides in 2 changes of 80% alcohol – 1 minute each at room temperature.

Wash slides in 2 changes of 95% alcohol - 1 minute each at room temperature.

Wash slides in 3 changes of 100% alcohol – 1 minute each at room temperature.

Wash slides in 3 changes of xylene – 1 minute each at room temperature.

Apply coverslip

Restrictions:

For Research Use only

Handling


| Format: | Liquid |
|--------------------|--|
| Concentration: | Lot specific |
| Buffer: | PBS, less than 0.1 % sodium azide. |
| Preservative: | Sodium azide |
| Precaution of Use: | This product contains Sodium azide: a POISONOUS AND HAZARDOUS SUBSTANCE which should be handled by trained staff only. |
| Storage: | 4 °C,-20 °C |
| Storage Comment: | Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for |

up to 1 month. Avoid freeze-thaw cycles.

Expiry Date:


12 months

Images



Immunohistochemistry

Image 1. Anti-BEST2 / Bestrophin-2 antibody ABIN1048350 IHC staining of human testis. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval.



Immunohistochemistry

Image 2. Anti-BEST2 / Bestrophin-2 antibody IHC of human Lymph Node, Non-Hodgkins Lymphoma. Immunohistochemistry of formalin-fixed, paraffinembedded tissue after heat-induced antigen retrieval.