antibodies.com

# Datasheet for ABIN1048454 anti-CNR2 antibody (Cytoplasmic Domain)

2 Images



### Overview

Quantity:	50 µg
Target:	CNR2
Binding Specificity:	Cytoplasmic Domain
Reactivity:	Human, Monkey, Dog
Host:	Rabbit
Clonality:	Polyclonal
Conjugate:	This CNR2 antibody is un-conjugated
Application:	Immunohistochemistry (IHC), Immunohistochemistry (Paraffin-embedded Sections) (IHC (p))

# Product Details

Brand:	IHC-plus™
Immunogen:	Synthetic 20 amino acid peptide from 3rd cytoplasmic domain of human CNR2 / CB2. Percent identity with other species by BLAST analysis: Human, Gorilla, Gibbon, Monkey, Marmoset, Dog (100%), Mouse, Rat, Elephant, Panda (95%), Hamster, Bat, Horse (90%), Rabbit (85%).
	Type of Immunogen: Synthetic peptide
Specificity:	Human CNR2 / CB2. BLAST analysis of the peptide immunogen showed no homology with other human proteins, except CNR1 (60 % ).
Predicted Reactivity:	Percent identity with other species by BLAST analysis: Human, Gorilla, Gibbon, Monkey, Marmoset, Dog (100%) Mouse, Rat, Elephant, Panda (95%) Hamster, Bat, Horse (90%) Rabbit (85%).

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/4 | Product datasheet for ABIN1048454 | 09/12/2023 | Copyright antibodies-online. All rights reserved.

# Product Details

### Purification: Immur

## Immunoaffinity purified

# Target Details

l'alget Details	
Target:	CNR2
Alternative Name:	CNR2 / CB2 (CNR2 Products)
Background:	Name/Gene ID: CNR2
	Subfamily: Cannabinoid
	Family: GPCR
	Synonyms: CNR2, Cannabinoid receptor 2, CB2, CB2 cannabinoid receptor, CX5, CB2R, HCB2
	Cannabinoid CB2 receptor, CB-2, CB2 receptor
Gene ID:	1269
Application Details	
Application Notes:	Approved: IHC, IHC-P (13 - 26 µg/mL)
Comment:	Target Species of Antibody: Human
Assay Procedure:	The IHC-pro Immunohistochemistry Protocol
	Tissue Preparation
	Formalin fixation and embedding in paraffin wax
	Tissue Sectioning
	Make 4-µm sections and place on pre-cleaned and charged microscope slides.
	Heat in a tissue-drying oven for 45 minutes at 60°C
	Deparaffinization
	Wash slides in 3 changes of xylene – 5 minutes each at room temperature.
	Rehydration
	Wash slides in 3 changes of 100% alcohol – 3 minutes each at room temperature.
	Wash slides in 2 changes of 95% alcohol – 3 minutes each at room temperature.
	Wash slides in 1 change of 80% alcohol – 3 minutes at room temperature.
	Rinse slides in gentle running distilled water – 5 minutes at room temperature.

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 2/4 | Product datasheet for ABIN1048454 | 09/12/2023 | Copyright antibodies-online. All rights reserved.

#### Antigen retrieval

Steam slides in 0.01 M sodium citrate buffer, pH 6.0 at 99-100°C - 20 minutes Remove from heat and let stand at room temperature in buffer - 20 minutes Rinse in 1X TBS with Tween (TBST) – 1 minute at room temperature.

#### Immunostaining

Do not allow tissues to dry at any time during the staining procedure. Apply a universal protein block – 20 minutes at room temperature. Drain protein block from slides, apply diluted primary antibody – 45 minutes at room temperature. Rinse slides in 1X TBST - 1 minute at room temperature. Apply a biotinylated secondary antibody (specific to the host of the primary antibody) - 30 minutes at room temperature. Rinse slides 1X TBST – 1 minute at room temperature. Apply alkaline phosphatase streptavidin – 30 minutes at room temperature. Rinse slides in 1X TBST - 1 minute at room temperature. Apply alkaline phosphatase streptavidin – 30 minutes at room temperature. May alkaline phosphatase chromogen substrate – 30 minutes at room temperature. Wash slides in distilled water – 1 minute at room temperature.

#### Dehydrate

This method should only be used if the chromogen substrate is alcohol insoluble. Wash slides in 2 changes of 80% alcohol – 1 minute each at room temperature. Wash slides in 2 changes of 95% alcohol – 1 minute each at room temperature. Wash slides in 3 changes of 100% alcohol – 1 minute each at room temperature. Wash slides in 3 changes of xylene – 1 minute each at room temperature. Apply coverslip

Restrictions:

For Research Use only

### Handling


Format:	Liquid
Concentration:	Lot specific
Buffer:	PBS, less than 0.1 % sodium azide.
Preservative:	Sodium azide

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 3/4 | Product datasheet for ABIN1048454 | 09/12/2023 | Copyright antibodies-online. All rights reserved.

Precaution of Use:	This product contains Sodium azide: a POISONOUS AND HAZARDOUS SUBSTANCE which should be handled by trained staff only.
Storage:	4 °C,-20 °C
Storage Comment:	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Expiry Date:	12 months


## Images

Handling



### Immunohistochemistry

**Image 1.** Anti-CNR2 / CB2 antibody ABIN1048454 IHC staining of human tonsil, germinal center and mantle zone. Immunohistochemistry of formalin-fixed, paraffinembedded tissue after heat-induced antigen retrieval.



### Immunohistochemistry

**Image 2.** Anti-CNR2 / CB2 antibody IHC of human Breast, Carcinoma. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval.

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 4/4 | Product datasheet for ABIN1048454 | 09/12/2023 | Copyright antibodies-online. All rights reserved.