
antibodies .- online.com







anti-CD6 antibody (APC)



Image



Overview

Quantity:	100 tests
Target:	CD6
Reactivity:	Human
Host:	Mouse
Clonality:	Monoclonal
Conjugate:	This CD6 antibody is conjugated to APC
Application:	Flow Cytometry (FACS), Immunofluorescence (IF)

Product Details

Clone:	MAE-1C10
Isotype:	IgG1
Characteristics:	Monoclonal Mouse Anti-Human CD6 APC is recommended for use in flow cytometry for identification of human mature T cells and a sub-population of B cells in peripheral blood.

Target Details

Target:	CD6
Alternative Name:	CD6 (CD6 Products)
Background:	The antibody reacts with a 120 kDa single chain membrane glycoprotein specific for the CD6 antigen. The antigen is found on >80% normal mature T cells and a sub-population of B cells in
	peripheral blood. The antigen also exhibits weak expression on cortical thymocytes and the majority of chronic B lymphocytic leukaemias

Application Details

It is recommended for use in flow cytometry. This reagent is effective for direct immunofluorescence staining of human tissue for flow cytometric analysis using 20 μ l/10^6 cells.
Allophycocyanin (Febico, Far East Bio-Tech Co.).
1. Transfer 100 μ l of anticoagulated (EDTA) blood to a 12 x 75 mm polystyrene test tube (10% cells). 2. Add 20 μ l of CD6 APC and mix gently with a vortex mixer. The 20 μ l is a guideline only the optimal volume should be determined by the individual laboratory. 3. The recommended negative control is a non-reactive APC-conjugated antibody of the same isotype. 4. Incubate in the dark at room temperature at 4°C for 30 minutes or at room temperature (20-25 °C) for 15 minutes. 5. Add 1,5 ml of Lysing Solution to each sample and mix gently with a vortex mixer. Incubate for 10 minutes at room temperature in the dark. 6. Centrifuge at 1000 x g for 5 minutes. Gently aspirate the supernatant and discard it leaving approximately 50 μ l of fluid. 7. Add 2 ml 0.01 mol/I PBS (It betters that it containing 2% bovine serum albumin) and resuspent the cells by using a vortex mixer. 8. Centrifuge at 1000 x g for 5 minutes. Gently aspirate the supernatant and discard it leaving approximately 50 μ l of fluid. 9. Resuspend pellet in an appropriate fluid for flow cytometry, e.g. 0.3 ml PBS. The PBS should contain 1% paraformaldehyde (fixative) if samples are not analysed the same day. 10. Analyse on a flow cytometer or store at 2-8 °C in the dark until analysis. Samples can be run up to 24 hours after lysis.
For Research Use only
Liquid
The conjugate is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 0,09% Sodium azide, pH 7.2.
Sodium azide
1. The device is not intended for clinical use including diagnosis, prognosis, and monitoring of disease state, and it must not be used in conjunction with patient records or treatment. 2. This product contains Sodium azide (NaN3), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, Sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. As with any product derive

from biological sources, proper handling procedures should be used.

Storage:

4°C

Images




Image 1.