antibodies -online.com





Rat anti-Mouse IgE (Heavy Chain) Antibody (HRP)



age



Publication



Go to Product page


Overview	
Quantity:	1 mL
Target:	IgE
Binding Specificity:	Heavy Chain
Reactivity:	Mouse
Host:	Rat
Clonality:	Monoclonal
Conjugate:	HRP
Application:	ELISA, Immunohistochemistry (IHC)
Product Details	
Clone:	23G3
Isotype:	lgG1
Specificity:	Reacts with the heavy chain of BALB/c mouse IgE as demonstrated by ELISA.
Purification:	Affinity chromatography using phosphopeptide, cross-adsorbed with the corresponding non-phosphopeptide
Target Details	
Target:	IgE
Abstract:	IgE Products

Application Details

ELISA: 1:1,000-1:2,000 Representative data are included in this product insert. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. Restrictions: For Research Use only Handling Format: Liquid Buffer: 1.0 mL of stock solution in 50 % glycerol / 50 % PBS, pH 7.4. Preservative: Without preservative Handling Advice: Avoid repeated freezing and thawing. Dilute only prior to immediate use. Centrifuge product if not completely clear after standing at room temperature. Do NDT add Sodium Azide! Use of Sodium Azide will inhibit enzyme activity of horseradish peroxidase. Storage: 4 *C/-20 *C Storage Comment: Store at 2-8 *C or long-term at -20 *C. Publications Product cited in: Welc, Flores, Wehling-Henricks, Ramos, Wang, Bertoni, Tidball: "Targeting a therapeutic LIF transgene to muscle via the immune system ameliorates muscular dystrophy." in: Nature communications, Vol. 10, Issue 1, pp. 2788, (2019) (PubMed). Finno, Gianino, Perumbakkarn, Williams, Bordbari, Gardner, Burns, Peng, Durward-Akhurst, Valberg: "A missense mutation in MYH1 is associated with susceptibility to immune-mediated	Application Notes:	Working Dilution:
Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. Restrictions: For Research Use only Handling Format: Liquid Buffer: 1.0 mL of stock solution in 50 % glycerol / 50 % PBS, pH 7.4. Preservative: Without preservative Handling Advice: Avoid repeated freezing and thawing. Dilute only prior to immediate use Centrifuge product if not completely clear after standing at room temperature. Do NOT add Sodium Azide! Use of Sodium Azide will inhibit enzyme activity of horseradish peroxidase. Storage: 4 °C/-20 °C Storage Comment: Store at 2-8 °C or long-term at -20 °C. Publications Product cited in: Welc, Flores, Wehling-Henricks, Ramos, Wang, Bertoni, Tidball: "Targeting a therapeutic LIF transgene to muscle via the immune system ameliorates muscular dystrophy." in: Nature communications, Vol. 10, Issue 1, pp. 2788, (2019) (PubMed).		ELISA: 1:1,000-1:2,000
Other applications have not been tested but use in such assays should not necessarily be excluded. Restrictions: For Research Use only Handling Format: Liquid Buffer: 1.0 mL of stock solution in 50 % glycerol / 50 % PBS, pH 7.4. Preservative: Without preservative Handling Advice: Avoid repeated freezing and thawing. Dilute only prior to immediate use Centrifuge product if not completely clear after standing at room temperature. Do NOT add Sodium Azidel Use of Sodium Azide will inhibit enzyme activity of horseradish peroxidase. Storage: 4 °C /- 20 °C Storage Comment: Store at 2-8 °C or long-term at -20 °C. Publications Product cited in: Welc, Flores, Wehling-Henricks, Ramos, Wang, Bertoni, Tidball: 'Targeting a therapeutic LIF transgene to muscle via the immune system ameliorates muscular dystrophy.' in: Nature communications, Vol. 10, Issue 1, pp. 2788, (2019) (PubMed). Finno, Gianino, Perumbakkam, Williams, Bordbari, Gardner, Burns, Peng, Durward-Akhurst,		Representative data are included in this product insert.
Other applications have not been tested but use in such assays should not necessarily be excluded. Restrictions: For Research Use only Handling Format: Liquid Buffer: 1.0 mL of stock solution in 50 % glycerol / 50 % PBS, pH 7.4. Preservative: Without preservative Handling Advice: Avoid repeated freezing and thawing. Dilute only prior to immediate use Centrifuge product if not completely clear after standing at room temperature. Do NOT add Sodium Azidel Use of Sodium Azide will inhibit enzyme activity of horseradish peroxidase. Storage: 4 °C /- 20 °C Storage Comment: Store at 2-8 °C or long-term at -20 °C. Publications Product cited in: Welc, Flores, Wehling-Henricks, Ramos, Wang, Bertoni, Tidball: 'Targeting a therapeutic LIF transgene to muscle via the immune system ameliorates muscular dystrophy.' in: Nature communications, Vol. 10, Issue 1, pp. 2788, (2019) (PubMed). Finno, Gianino, Perumbakkam, Williams, Bordbari, Gardner, Burns, Peng, Durward-Akhurst,		
excluded. Restrictions: For Research Use only Handling Format: Liquid Buffer: 1.0 mL of stock solution in 50 % glycerol / 50 % PBS, pH 7.4. Preservative: Without preservative Avoid repeated freezing and thawing. Dilute only prior to immediate use Centrifuge product if not completely clear after standing at room temperature. Do NOT add Sodium Azide! Use of Sodium Azide will inhibit enzyme activity of horseradish peroxidase. Storage: 4 °C/-20 °C Storage Comment: Store at 2-8 °C or long-term at -20 °C. Publications Product cited in: Welc, Flores, Wehling-Henricks, Ramos, Wang, Bertoni, Tidball: "Targeting a therapeutic LIF transgene to muscle via the immune system ameliorates muscular dystrophy." in: Nature communications, Vol. 10, Issue 1, pp. 2788, (2019) (PubMed). Finno, Gianino, Perumbakkam, Williams, Bordbari, Gardner, Burns, Peng, Durward-Akhurst,		
Restrictions: For Research Use only Handling Format: Liquid Buffer: 1.0 mL of stock solution in 50 % glycerol / 50 % PBS, pH 7.4. Preservative: Without preservative Avoid repeated freezing and thawing. Dilute only prior to immediate use Centrifuge product if not completely clear after standing at room temperature. Do NOT add Sodium Azide! Use of Sodium Azide will inhibit enzyme activity of horseradish peroxidase. Storage: 4 °C/-20 °C Storage Comment: Store at 2-8 °C or long-term at -20 °C. Publications Product cited in: Welc, Flores, Wehling-Henricks, Ramos, Wang, Bertoni, Tidball: "Targeting a therapeutic LIF transgene to muscle via the immune system ameliorates muscular dystrophy." in: Nature communications, Vol. 10, Issue 1, pp. 2788, (2019) (PubMed). Finno, Gianino, Perumbakkam, Williams, Bordbari, Gardner, Burns, Peng, Durward-Akhurst,		
Format: Liquid Buffer: 1.0 mL of stock solution in 50 % glycerol / 50 % PBS, pH 7.4. Preservative: Without preservative Handling Advice: Avoid repeated freezing and thawing. Dilute only prior to immediate use Centrifuge product if not completely clear after standing at room temperature. Do NOT add Sodium Azide! Use of Sodium Azide will inhibit enzyme activity of horseradish peroxidase. Storage: 4 °C/-20 °C Storage Comment: Store at 2-8 °C or long-term at -20 °C. Publications Product cited in: Welc, Flores, Wehling-Henricks, Ramos, Wang, Bertoni, Tidball: "Targeting a therapeutic LIF transgene to muscle via the immune system ameliorates muscular dystrophy." in: Nature communications, Vol. 10, Issue 1, pp. 2788, (2019) (PubMed). Finno, Gianino, Perumbakkam, Williams, Bordbari, Gardner, Burns, Peng, Durward-Akhurst,		excluded.
Format: Liquid Buffer: 1.0 mL of stock solution in 50 % glycerol / 50 % PBS, pH 7.4. Preservative: Without preservative Handling Advice: Avoid repeated freezing and thawing. Dilute only prior to immediate use Centrifuge product if not completely clear after standing at room temperature. Do NOT add Sodium Azide! Use of Sodium Azide will inhibit enzyme activity of horseradish peroxidase. Storage: 4 °C/-20 °C Storage Comment: Store at 2-8 °C or long-term at -20 °C. Publications Product cited in: Welc, Flores, Wehling-Henricks, Ramos, Wang, Bertoni, Tidball: "Targeting a therapeutic LIF transgene to muscle via the immune system ameliorates muscular dystrophy." in: Nature communications, Vol. 10, Issue 1, pp. 2788, (2019) (PubMed). Finno, Gianino, Perumbakkam, Williams, Bordbari, Gardner, Burns, Peng, Durward-Akhurst,	Restrictions:	For Research Use only
Buffer: 1.0 mL of stock solution in 50 % glycerol / 50 % PBS, pH 7.4. Preservative: Without preservative Handling Advice: Avoid repeated freezing and thawing. Dilute only prior to immediate use Centrifuge product if not completely clear after standing at room temperature. Do NOT add Sodium Azide! Use of Sodium Azide will inhibit enzyme activity of horseradish peroxidase. Storage: 4 °C/-20 °C Storage Comment: Store at 2-8 °C or long-term at -20 °C. Publications Product cited in: Welc, Flores, Wehling-Henricks, Ramos, Wang, Bertoni, Tidball: "Targeting a therapeutic LIF transgene to muscle via the immune system ameliorates muscular dystrophy." in: Nature communications, Vol. 10, Issue 1, pp. 2788, (2019) (PubMed). Finno, Gianino, Perumbakkam, Williams, Bordbari, Gardner, Burns, Peng, Durward-Akhurst,	Handling	
Preservative: Without preservative Handling Advice: Avoid repeated freezing and thawing. Dilute only prior to immediate use Centrifuge product if not completely clear after standing at room temperature. Do NOT add Sodium Azide! Use of Sodium Azide will inhibit enzyme activity of horseradish peroxidase. Storage: 4 °C/-20 °C Storage Comment: Store at 2-8 °C or long-term at -20 °C. Publications Product cited in: Welc, Flores, Wehling-Henricks, Ramos, Wang, Bertoni, Tidball: "Targeting a therapeutic LIF transgene to muscle via the immune system ameliorates muscular dystrophy." in: Nature communications, Vol. 10, Issue 1, pp. 2788, (2019) (PubMed). Finno, Gianino, Perumbakkam, Williams, Bordbari, Gardner, Burns, Peng, Durward-Akhurst,	Format:	Liquid
Handling Advice: Avoid repeated freezing and thawing. Dilute only prior to immediate use Centrifuge product if not completely clear after standing at room temperature. Do NOT add Sodium Azide! Use of Sodium Azide will inhibit enzyme activity of horseradish peroxidase. Storage: 4 °C/-20 °C Storage Comment: Store at 2-8 °C or long-term at -20 °C. Publications Product cited in: Welc, Flores, Wehling-Henricks, Ramos, Wang, Bertoni, Tidball: "Targeting a therapeutic LIF transgene to muscle via the immune system ameliorates muscular dystrophy." in: Nature communications, Vol. 10, Issue 1, pp. 2788, (2019) (PubMed). Finno, Gianino, Perumbakkam, Williams, Bordbari, Gardner, Burns, Peng, Durward-Akhurst,	Buffer:	1.0 mL of stock solution in 50 % glycerol / 50 % PBS, pH 7.4.
Centrifuge product if not completely clear after standing at room temperature. Do NOT add Sodium Azide! Use of Sodium Azide will inhibit enzyme activity of horseradish peroxidase. Storage: 4 °C/-20 °C Storage Comment: Store at 2-8 °C or long-term at -20 °C. Publications Product cited in: Welc, Flores, Wehling-Henricks, Ramos, Wang, Bertoni, Tidball: "Targeting a therapeutic LIF transgene to muscle via the immune system ameliorates muscular dystrophy." in: Nature communications, Vol. 10, Issue 1, pp. 2788, (2019) (PubMed). Finno, Gianino, Perumbakkam, Williams, Bordbari, Gardner, Burns, Peng, Durward-Akhurst,	Preservative:	Without preservative
Do NOT add Sodium Azide! Use of Sodium Azide will inhibit enzyme activity of horseradish peroxidase. Storage: 4 °C/-20 °C Storage Comment: Store at 2-8 °C or long-term at -20 °C. Publications Product cited in: Welc, Flores, Wehling-Henricks, Ramos, Wang, Bertoni, Tidball: "Targeting a therapeutic LIF transgene to muscle via the immune system ameliorates muscular dystrophy." in: Nature communications, Vol. 10, Issue 1, pp. 2788, (2019) (PubMed). Finno, Gianino, Perumbakkam, Williams, Bordbari, Gardner, Burns, Peng, Durward-Akhurst,	Handling Advice:	Avoid repeated freezing and thawing. Dilute only prior to immediate use
peroxidase. Storage: 4 °C/-20 °C Storage Comment: Store at 2-8 °C or long-term at -20 °C. Publications Product cited in: Welc, Flores, Wehling-Henricks, Ramos, Wang, Bertoni, Tidball: "Targeting a therapeutic LIF transgene to muscle via the immune system ameliorates muscular dystrophy." in: Nature communications, Vol. 10, Issue 1, pp. 2788, (2019) (PubMed). Finno, Gianino, Perumbakkam, Williams, Bordbari, Gardner, Burns, Peng, Durward-Akhurst,		Centrifuge product if not completely clear after standing at room temperature.
Storage: 4 °C/-20 °C Storage Comment: Store at 2-8 °C or long-term at -20 °C. Publications Product cited in: Welc, Flores, Wehling-Henricks, Ramos, Wang, Bertoni, Tidball: "Targeting a therapeutic LIF transgene to muscle via the immune system ameliorates muscular dystrophy." in: Nature communications, Vol. 10, Issue 1, pp. 2788, (2019) (PubMed). Finno, Gianino, Perumbakkam, Williams, Bordbari, Gardner, Burns, Peng, Durward-Akhurst,		Do NOT add Sodium Azide! Use of Sodium Azide will inhibit enzyme activity of horseradish
Storage Comment: Store at 2-8 °C or long-term at -20 °C. Publications Product cited in: Welc, Flores, Wehling-Henricks, Ramos, Wang, Bertoni, Tidball: "Targeting a therapeutic LIF transgene to muscle via the immune system ameliorates muscular dystrophy." in: Nature communications, Vol. 10, Issue 1, pp. 2788, (2019) (PubMed). Finno, Gianino, Perumbakkam, Williams, Bordbari, Gardner, Burns, Peng, Durward-Akhurst,		peroxidase.
Product cited in: Welc, Flores, Wehling-Henricks, Ramos, Wang, Bertoni, Tidball: "Targeting a therapeutic LIF transgene to muscle via the immune system ameliorates muscular dystrophy." in: Nature communications, Vol. 10, Issue 1, pp. 2788, (2019) (PubMed). Finno, Gianino, Perumbakkam, Williams, Bordbari, Gardner, Burns, Peng, Durward-Akhurst,	Storage:	4 °C/-20 °C
Product cited in: Welc, Flores, Wehling-Henricks, Ramos, Wang, Bertoni, Tidball: "Targeting a therapeutic LIF transgene to muscle via the immune system ameliorates muscular dystrophy." in: Nature communications , Vol. 10, Issue 1, pp. 2788, (2019) (PubMed). Finno, Gianino, Perumbakkam, Williams, Bordbari, Gardner, Burns, Peng, Durward-Akhurst,	Storage Comment:	Store at 2-8 °C or long-term at -20 °C.
transgene to muscle via the immune system ameliorates muscular dystrophy." in: Nature communications , Vol. 10, Issue 1, pp. 2788, (2019) (PubMed). Finno, Gianino, Perumbakkam, Williams, Bordbari, Gardner, Burns, Peng, Durward-Akhurst,	Publications	
communications, Vol. 10, Issue 1, pp. 2788, (2019) (PubMed). Finno, Gianino, Perumbakkam, Williams, Bordbari, Gardner, Burns, Peng, Durward-Akhurst,	Product cited in:	Welc, Flores, Wehling-Henricks, Ramos, Wang, Bertoni, Tidball: "Targeting a therapeutic LIF
Finno, Gianino, Perumbakkam, Williams, Bordbari, Gardner, Burns, Peng, Durward-Akhurst,		transgene to muscle via the immune system ameliorates muscular dystrophy." in: Nature
		communications , Vol. 10, Issue 1, pp. 2788, (2019) (PubMed).
Valberg: "A missense mutation in MYH1 is associated with susceptibility to immune-mediated		Finno, Gianino, Perumbakkam, Williams, Bordbari, Gardner, Burns, Peng, Durward-Akhurst,
		Valberg: "A missense mutation in MYH1 is associated with susceptibility to immune-mediated

myositis in Quarter Horses." in: Skeletal muscle, Vol. 8, Issue 1, pp. 7, (2018) (PubMed).

Huang, Ge, Izzi, Greenspan: "a3 Chains of type V collagen regulate breast tumour growth via glypican-1." in: Nature communications, Vol. 8, pp. 14351, (2018) (PubMed).



Immunohistochemistry

Image 1. Paraffin embedded mouse small intestine section was stained with Rat Anti-Mouse IgE-UNLB followed by an HRP conjugated secondary antibody, DAB, and hematoxylin.