ANTIBODIES ONLINE

Datasheet for ABIN1379973 IL18BP ELISA Kit

Image

Go to Product page

Overview

Quantity:	96 tests
Target:	IL18BP
Reactivity:	Human
Method Type:	Sandwich ELISA
Detection Range:	94-6000 pg/mL
Minimum Detection Limit:	94 pg/mL
Application:	ELISA

Product Details

Purpose:

The OmniKine? Human IL-18 BPa ELISA Kit contains the components necessary for quantitative determination of natural or recombinant Human IL-18 BPa concentrations within any experimental sample including cell lysates, serum and plasma. This particular immunoassay utilizes the quantitative technique of a "Sandwich" Enzyme- Linked Immunosorbent Assay (ELISA) where the target protein (antigen) is bound in a "sandwich" format by the primary capture antibodies coated to each well-bottom and the secondary detection antibodies added subsequently by the investigator. The capture antibodies coated to the bottom of each well are specific for a particular epitope on Human IL-18 BPa while the useradded detection antibodies bind to epitopes on the captured target protein. Amid each step of the procedure, a series of wash steps must be performed to ensure the elimination of nonspecific binding between proteins to other proteins or to the solid phase. After incubation and "sandwiching" of the target antigen, a peroxidase enzyme is conjugated to the constant heavy chain of the secondary antibody (either covalently or via Avidin/Streptavidin-Biotin interactions), allowing for a colorimetric reaction to ensue upon substrate addition. When the substrate TMB

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/6 | Product datasheet for ABIN1379973 | 07/26/2024 | Copyright antibodies-online. All rights reserved. Brand:

(3, 3', 5, 5'-Tetramethylbenzidine) is added, the reaction catalyzed by peroxidase yields a blue color that is representative of the antigen concentration. Upon sufficient color development, the reaction can be terminated through addition of Stop Solution (2 N Sulfuric Acid) where the color of the solution will turn yellow. The absorbance of each well can then be read by a spectrophotometer, allowing for generation of a standard curve and subsequent determination of protein concentration. OmniKine™ Sample Type: Cell Lysate, Serum, Plasma Analytical Method: Quantitative **Detection Method:** Colorimetric Specificity: The Human IL-18 BPa ELISA Kit allows for the detection and quantification of endogenous levels of natural and/or recombinant Human IL-18 BPa proteins. Cross-Reactivity (Details): The Human IL-18 BPa ELISA is capable of recognizing both recombinant and naturally

produced Human IL-18 BPa proteins. The antigens listed below were tested at 50 ng/mL and did not exhibit significant cross reactivity or interference. Human: IL-18 BPc, IL-18BPd Characteristics: The Human IL-18 BPa ELISA Kit allows for the detection and quantification of endogenous levels of natural and/or recombinant Human IL-18 BPa proteins within the range of 94-6000 pg/mL.

Components:	•	Microstrips Coated w / Capture Antibody: 12 x 8-Well Microstrips
	•	Protein Standard: Lyophilized (100 ng), Red container
	•	Biotinylated Detection Antibody: Lyophilized, Yellow container

- 400x Streptavidin-HRP: 30 µL, Blue container
- Wash Buffer (10x): 50 mL, Clear containter
- · Assay Diluent: 50 mL, Clear container
 - Ready-to-Use Substrate: 12 mL, Brown container
 - Stop Solution: 12 mL, Clear container
 - Adhesive Plate Sealers: 4 Sheets
- Technical Manual 1 Manual

Material not included:	The following materials and/or equipment are NOT provided in this kit but are necessary to
	successfully conduct the experiment:
	Microplate reader able to measure absorbance at 450 nm (with correction wavelength set to
	540 nm or 570 nm)
	Micropipettes with capability of measuring volumes ranging from 1 μl to 1 mL

Deionized or sterile water

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 2/6 | Product datasheet for ABIN1379973 | 07/26/2024 | Copyright antibodies-online. All rights reserved.

Squirt bottle, manifold dispenser, multichannel pipette reservoir or automated microplate washer Graph paper or computer software capable of generating or displaying logarithmic functions Absorbent paper or vacuum aspirator Test tubes or microfuge tubes capable of storing ≥1 mL Bench top centrifuge (optional) Bench top vortex (optional)

Target Details

Target:	IL18BP
Alternative Name:	IL-18 BPa (IL18BP Products)
Background:	The protein encoded by IL-18BP functions as an inhibitor of the proinflammatory cytokine, IL18.
	It binds IL18, prevents the binding of IL18 to its receptor, and thus inhibits IL18-induced IFN-
	gamma production, resulting in reduced T-helper type 1 immune responses. This protein is
	constitutively expressed and secreted in mononuclear cells. Elevated level of this protein is
	detected in the intestinal tissues of patients with Crohn's disease. Alternatively spliced
	transcript variants encoding different isoforms have been described for this gene. Isoform A
	binds to IL-18 and inhibits its activity while IL-18BP functions as an inhibitor of the early TH1
	cytokine response and is strongly expressed in heart, lung, placenta and spleen. Source: Entrez
	Gene, Swiss-Prot

Application Details

Plate:	Pre-coated
Protocol:	This particular immunoassay utilizes the quantitative technique of a Sandwich Enzyme-Linked
	Immunosorbent Assay (ELISA) where the target protein (antigen) is bound in a sandwich
	format by the primary capture antibodies coated to each well-bottom and the secondary
	detection antibodies added subsequently by the investigator. The capture antibodies coated to
	the bottom of each well are specific for a particular epitope on the Human IL-18 BPa cytokine
	while the user-added detection antibodies bind to epitopes on the captured target protein. Amid
	each step of the procedure, a series of wash steps must be performed to ensure the elimination
	of non-specific binding between proteins to other proteins or to the solid phase. After

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 3/6 | Product datasheet for ABIN1379973 | 07/26/2024 | Copyright antibodies-online. All rights reserved.

Application Details

	incubation and sandwiching of the target antigen, a peroxidase enzyme is conjugated to the constant heavy chain of the secondary antibody (either covalently or via Avidin/Streptavidin-
	Biotin interactions), allowing for a colorimetric reaction to ensue upon substrate addition. When
	the substrate TMB (3, 3', 5, 5'-Tetramethylbenzidine) is added, the reaction catalyzed by
	peroxidase yields a blue color that is representative of the antigen concentration. Upon
	sufficient color development, the reaction can be terminated through addition of Stop Solution
	(2 N Sulfuric Acid) where the color of the solution will turn yellow. The absorbance of each well
	can then be read by a spectrophotometer, allowing for generation of a standard curve and
	subsequent determination of protein concentration.
Sample Preparation:	If samples are to be used within 24 hours, aliquot and store at 4 °C. If samples are to be used
	over a long period of time, aliquot and store between -20 °C and -80 °C, depending on the
	duration of storage.
	Note: Samples containing a visible precipitate or pellet must be clarified prior to use in the
	assay.
	Caution: Avoid repeated freeze/thaw cycles to prevent loss of biological activity of proteins in
	experimental samples.
	Cell Lysate and Supernatants:
	 Remove large cell components via centrifugation and perform the assay. Cell lysates and supernatants require a dilution using Assay Diluent. A serial dilution may be performed to determine a suitable dilution factor for the sample. For future use of the sample, follow the sample storage guidelines stated above. Serum:
	 Allow samples to clot in a serum separator tube (SST) for 30 minutes. After sufficient clotting, centrifuge at 1000 x g for 15 minutes and remove serum from SST in preparation for the assay. Serum samples require at least a 1:50 dilution using Assay Diluent. For future use of the sample, follow the storage guidelines above. Plasma:
	Use heparin, citrate or EDTA as an anticoagulant to gather plasma from original biological sample. After collection of the plasma, centrifuge for 15 minutes at 1000 x g. This step must be performed within 30 minutes of plasma collection. Plasma samples require at least a 1:50 dilution using Assay Diluent. Afterwards, perform the assay or for future use of the sample, follow the storage guidelines stated above.
Assay Procedure:	Note: If possible, all incubation steps should be performed on an orbital shaker to equilibrate
	solutions when added to the microplate wells. Also, all provided solutions should be at ambient
	temperature prior to use.
	Note: Avoid adding solutions into wells at an angle, always keep pipette tip perpendicular to

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 4/6 | Product datasheet for ABIN1379973 | 07/26/2024 | Copyright antibodies-online. All rights reserved.

Application Details

	Reconstitution of Provided Materials:
	1. Reconstitute the Biotin-Conjugated Detection Antibody in 67 μ L of ddH IO for a concentration of 180 μ g/ml.
	2. Reconstitute the Protein Standard in 100 μ L of ddH \mathbb{I} O for a concentration of 340 ng/ml.
	3. Dilute the 50 mL of 10x Wash Buffer in 450 mL of ddH20 for 500 mL of 1x Wash Buffer.
	Addition of Known Standard and Unknown Sample to Immunoassay:
	The OmniKine™ Human CD163 ELISA Kit allows for the detection and quantification of endogenous levels of natural and/or recombinant Human CD163 proteins
Calculation of Results:	Generation of Standard Curve and Interpretation of Data
	1. Average the duplicate or triplicate readings for each standard, control and sample and
	subtract the average zero standard optical density.
	2. Generate a standard curve by using Microsoft Excel or other computer software capable of
	establishing a 4- Parameter Logistic (4-PL) curve fit. If using Excel or an alternative graphing
	tool, plot the average optical density values in absorbance units (y-axis) against the known
	standard concentrations in pg/ml (x-axis). Note: Only use the values in which a noticeable
	gradient can be established. Afterwards, generate a best fit curve or trend-line through the
	plotted points via regression analysis.
Restrictions:	For Research Use only
Handling	
Handling Precaution of Use:	Reagents provided in this kit may be harmful if ingested, inhaled or absorbed through the skin.
	Reagents provided in this kit may be harmful if ingested, inhaled or absorbed through the skin. Please carefully review the MSDS for each reagent before conducting the experiment.
	Please carefully review the MSDS for each reagent before conducting the experiment.
	Please carefully review the MSDS for each reagent before conducting the experiment. Stop Solution contains 2 N Sulfuric Acid (H2SO4) and is an extremely corrosive agent. Please
	Please carefully review the MSDS for each reagent before conducting the experiment. Stop Solution contains 2 N Sulfuric Acid (H2SO4) and is an extremely corrosive agent. Please wear proper eye, hand and face protection when handling this material. When the experiment is
Precaution of Use:	Please carefully review the MSDS for each reagent before conducting the experiment. Stop Solution contains 2 N Sulfuric Acid (H2SO4) and is an extremely corrosive agent. Please wear proper eye, hand and face protection when handling this material. When the experiment is finished, be sure to rinse the plate with copious amounts of running water to dilute the Stop
Precaution of Use:	Please carefully review the MSDS for each reagent before conducting the experiment. Stop Solution contains 2 N Sulfuric Acid (H2SO4) and is an extremely corrosive agent. Please wear proper eye, hand and face protection when handling this material. When the experiment is finished, be sure to rinse the plate with copious amounts of running water to dilute the Stop Solution prior to disposing the plate.
Precaution of Use:	 Please carefully review the MSDS for each reagent before conducting the experiment. Stop Solution contains 2 N Sulfuric Acid (H2SO4) and is an extremely corrosive agent. Please wear proper eye, hand and face protection when handling this material. When the experiment is finished, be sure to rinse the plate with copious amounts of running water to dilute the Stop Solution prior to disposing the plate. This ELISA kit is intended for research purposes only, NOT diagnostic or clinical procedures of
Precaution of Use:	 Please carefully review the MSDS for each reagent before conducting the experiment. Stop Solution contains 2 N Sulfuric Acid (H2SO4) and is an extremely corrosive agent. Please wear proper eye, hand and face protection when handling this material. When the experiment is finished, be sure to rinse the plate with copious amounts of running water to dilute the Stop Solution prior to disposing the plate. This ELISA kit is intended for research purposes only, NOT diagnostic or clinical procedures of any kind.
Precaution of Use:	 Please carefully review the MSDS for each reagent before conducting the experiment. Stop Solution contains 2 N Sulfuric Acid (H2SO4) and is an extremely corrosive agent. Please wear proper eye, hand and face protection when handling this material. When the experiment is finished, be sure to rinse the plate with copious amounts of running water to dilute the Stop Solution prior to disposing the plate. This ELISA kit is intended for research purposes only, NOT diagnostic or clinical procedures of any kind. Materials included in this kit should NOT be used past the expiration date on the kit label.
	 Please carefully review the MSDS for each reagent before conducting the experiment. Stop Solution contains 2 N Sulfuric Acid (H2SO4) and is an extremely corrosive agent. Please wear proper eye, hand and face protection when handling this material. When the experiment is finished, be sure to rinse the plate with copious amounts of running water to dilute the Stop Solution prior to disposing the plate. This ELISA kit is intended for research purposes only, NOT diagnostic or clinical procedures of any kind. Materials included in this kit should NOT be used past the expiration date on the kit label. Reagents or substrates included in this kit should NOT be mixed or substituted with reagents or
Precaution of Use:	 Please carefully review the MSDS for each reagent before conducting the experiment. Stop Solution contains 2 N Sulfuric Acid (H2SO4) and is an extremely corrosive agent. Please wear proper eye, hand and face protection when handling this material. When the experiment is finished, be sure to rinse the plate with copious amounts of running water to dilute the Stop Solution prior to disposing the plate. This ELISA kit is intended for research purposes only, NOT diagnostic or clinical procedures of any kind. Materials included in this kit should NOT be used past the expiration date on the kit label. Reagents or substrates included in this kit should NOT be mixed or substituted with reagents or substrates from any other kits.

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 5/6 | Product datasheet for ABIN1379973 | 07/26/2024 | Copyright antibodies-online. All rights reserved.

	The assay is designed to eliminate interference and background by other cellular
	macromolecules or factors present within any biological samples. However, the possibility of
	background noise cannot be fully excluded until all factors have been tested using the assay kit
	Reagents provided in this kit may be harmful if ingested, inhaled or absorbed through the skin.
	Please carefully review the MSDS for each reagent before conducting the experiment.
	Stop Solution contains 2 N Sulfuric Acid (H2SO4) and is an extremely corrosive agent. Please
	wear proper eye, hand and face protection when handling this material. When the experiment is
	finished, be sure to rinse the plate with copious amounts of running water to dilute the Stop
	Solution prior to disposing the plate.
Storage:	4 °C
Storage Comment:	Note: If used frequently, reagents may be stored at 4 °C.
	Unopened Kits: Store at 4 °C for 6 months.
	• Microstrips Coated w/ Capture Antibody, 400x Streptavidin-HRP Wash Buffer (10x), Assay
	Diluent Ready-to-Use Substrate, Stop Solution: 6 Months at 4 °C
	Protein Standard, Biotinylated Detection Antibody: Lyophilized: 6 Months (if Reconstituted: 1

Protein Standard, Biotinylated Detection Antibody: Lyophilized: 6 Months (if Reconstituted: 1 Month) at 4 °C

Images



ELISA

Image 1. This is an example of what a typical standard curve will look like. You must make your own standard curve. Do not use this example as your own standard curve.

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 6/6 | Product datasheet for ABIN1379973 | 07/26/2024 | Copyright antibodies-online. All rights reserved.