



[Go to Product page](#)

Datasheet for ABIN1390460
anti-INSL3 antibody (AA 11-110) (FITC)

Overview

Quantity:	100 µL
Target:	INSL3
Binding Specificity:	AA 11-110
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Conjugate:	This INSL3 antibody is conjugated to FITC
Application:	Immunofluorescence (Cultured Cells) (IF (cc)), Immunofluorescence (Paraffin-embedded Sections) (IF (p))

Product Details

Immunogen:	KLH conjugated synthetic peptide derived from human INSL3
Isotype:	IgG
Predicted Reactivity:	Human
Purification:	Purified by Protein A.

Target Details

Target:	INSL3
Alternative Name:	INSL3 (INSL3 Products)
Background:	Synonyms: INSL 3, INSL-3, INSL3_HUMAN, Insulin like 3, Insulin-like 3 A chain, Ley IL, Ley-I-L,

Target Details

Leydig insulin like peptide, Leydig insulin-like peptide, Relaxin like factor, Relaxin-like factor, RLF, RLNL.

Background: This gene encodes a member of the insulin-like hormone superfamily. The encoded protein is mainly produced in gonadal tissues. Studies of the mouse counterpart suggest that this gene may be involved in the development of urogenital tract and female fertility. This protein may also act as a hormone to regulate growth and differentiation of gubernaculum, and thus mediating intra-abdominal testicular descent. Mutations in this gene may lead to cryptorchidism. Alternate splicing results in multiple transcript variants. [provided by RefSeq, May 2012].

Pathways: [Hormone Activity](#), [cAMP Metabolic Process](#)

Application Details

Application Notes: IF(IHC-P) 1:50-200
IF(IHC-F) 1:50-200
IF(ICC) 1:50-200

Restrictions: For Research Use only

Handling

Format: Liquid

Concentration: 1 µg/µL

Buffer: Aqueous buffered solution containing 0.01M TBS (pH 7.4) with 1 % BSA, 0.03 % Proclin300 and 50 % Glycerol.

Preservative: ProClin

Precaution of Use: This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be handled by trained staff only.

Storage: -20 °C

Storage Comment: Store at -20°C. Aliquot into multiple vials to avoid repeated freeze-thaw cycles.

Expiry Date: 12 months