antibodies .- online.com







anti-ANKS1B antibody (Alexa Fluor 350)



\sim	
()\/△	rview
\cup	1 410 44

Quantity:	100 μL
Target:	ANKS1B
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Conjugate:	This ANKS1B antibody is conjugated to Alexa Fluor 350
Application:	Western Blotting (WB), Immunofluorescence (Paraffin-embedded Sections) (IF (p))

Product Details

Immunogen:	KLH conjugated synthetic peptide derived from human ANKS1B/AIDA1
Isotype:	IgG
Cross-Reactivity:	Human, Mouse, Rat
Purification:	Purified by Protein A.

Target Details

Target:	ANKS1B
Alternative Name:	ANKS1B (ANKS1B Products)
Background:	Synonyms: AIDA 1, AIDA-1, Amyloid-beta protein intracellular domain-associated protein 1,
	Anks1b, Ankyrin repeat and sterile alpha mot domain-containing protein 1B, ANS1B_HUMAN,
	E2A-PBX1-associated protein, EB-1.
	Background: TheAmyloid protein precursor (AbPP) is a widely expressed transmembrane

Target Details

protein that is processed into the b-Amyloid (Ab) peptide, which accumulates in insoluble plaques in the brain of Alzheimer?s disease patients and AbPP intracellular domain (AID). AID may function as a pro-apoptotic peptide, a regulator of calcium homeostasis and a molecule involved in transcriptional regulation. The AID associated protein 1 (AIDA-1) is highly expressed in the brain and is regulated by AbPP. It interacts with AbPP to play a role in brain development. AIDA-1 also interacts with coilin in Cajal bodies to regulate pre-mRNA splicing.

Gene ID:

56899

UniProt:

Q7Z6G8

Application Details

Application Notes: IF(IHC-P) 1:50-200

Restrictions: For Research Use only

Handling

Format:	Liquid
Concentration:	1 μg/μL
Buffer:	Aqueous buffered solution containing 0.01M TBS (pH 7.4) with 1 % BSA, 0.03 % Proclin300 and 50 % Glycerol.
Preservative:	ProClin
Precaution of Use:	This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be handled by trained staff only.
Storage:	-20 °C
Storage Comment:	Store at -20°C. Aliquot into multiple vials to avoid repeated freeze-thaw cycles.
Expiry Date:	12 months