Datasheet for ABIN1401485

antibodies.com

## anti-SERPINB10 antibody (AbBy Fluor® 555)



Overview

Quantity:	100 μL
Target:	SERPINB10
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Conjugate:	This SERPINB10 antibody is conjugated to AbBy Fluor® 555
Application:	Western Blotting (WB), Immunofluorescence (Paraffin-embedded Sections) (IF (p))
Product Details	
Immunogen:	KLH conjugated synthetic peptide derived from human SERPINB10
lsotype:	lgG
Cross-Reactivity:	Human
Purification:	Purified by Protein A.
Target Details	
Target:	SERPINB10
Alternative Name:	SERPINB10 (SERPINB10 Products)
Background:	Synonyms: Bomapin, PI10, Protease inhibitor 10, Protease inhibitor 10 ovalbumin type bomapin,
	Serpin B10, Serpin peptidase inhibitor clade B ovalbumin member 10, SPB10_HUMAN.
	Background: SERPINB10 is a serine proteinase inhibitor of the ovalbumin like B clade of serpins.
	It was first discovered in human bone marrow in a search for serpins similar to SerpinB6.

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/2 | Product datasheet for ABIN1401485 | 06/07/2024 | Copyright antibodies-online. All rights reserved.

Target Details	
	SERPINB10 was identified in lymphocytes, and was elevated in chronic myeloid leukemia and chronic myelomonocytic leukemia cells. The monocytic cell lines THP-1 and AML-193 also showed elevated SERPINB10 expression.
Gene ID:	5273
Application Details	
Application Notes:	IF(IHC-P) 1:50-200
Restrictions:	For Research Use only
Handling	
Format:	Liquid
Concentration:	1 µg/µL
Concentration: Buffer:	•
	1 μg/μL Aqueous buffered solution containing 0.01M TBS ( pH 7.4) with 1 % BSA, 0.03 % Proclin300 and
Buffer:	1 μg/μL Aqueous buffered solution containing 0.01M TBS ( pH 7.4) with 1 % BSA, 0.03 % Proclin300 and 50 % Glycerol.
Buffer: Preservative:	1 μg/μL         Aqueous buffered solution containing 0.01M TBS ( pH 7.4) with 1 % BSA, 0.03 % Proclin300 and 50 % Glycerol.         ProClin         This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be
Buffer: Preservative: Precaution of Use:	<ul> <li>1 μg/μL</li> <li>Aqueous buffered solution containing 0.01M TBS ( pH 7.4) with 1 % BSA, 0.03 % Proclin300 and 50 % Glycerol.</li> <li>ProClin</li> <li>This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be handled by trained staff only.</li> </ul>