antibodies

Datasheet for ABIN1401486 anti-SERPINB10 antibody (Alexa Fluor 647)



| | Go to Froduct page |
|-------------------|--|
| Overview | |
| Quantity: | 100 µL |
| Target: | SERPINB10 |
| Reactivity: | Human |
| Host: | Rabbit |
| Clonality: | Polyclonal |
| Conjugate: | This SERPINB10 antibody is conjugated to Alexa Fluor 647 |
| Application: | Western Blotting (WB), Immunofluorescence (Paraffin-embedded Sections) (IF (p)) |
| Product Details | |
| Immunogen: | KLH conjugated synthetic peptide derived from human SERPINB10 |
| lsotype: | lgG |
| Cross-Reactivity: | Human |
| Purification: | Purified by Protein A. |
| Target Details | |
| Target: | SERPINB10 |
| Alternative Name: | SERPINB10 (SERPINB10 Products) |
| Background: | Synonyms: Bomapin, PI10, Protease inhibitor 10, Protease inhibitor 10 ovalbumin type bomapin, |
| | Serpin B10, Serpin peptidase inhibitor clade B ovalbumin member 10, SPB10_HUMAN. |
| | Background: SERPINB10 is a serine proteinase inhibitor of the ovalbumin like B clade of serpins. |

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/2 | Product datasheet for ABIN1401486 | 03/07/2024 | Copyright antibodies-online. All rights reserved.

It was first discovered in human bone marrow in a search for serpins similar to SerpinB6.

| Target Details | |
|--|---|
| | SERPINB10 was identified in lymphocytes, and was elevated in chronic myeloid leukemia and chronic myelomonocytic leukemia cells. The monocytic cell lines THP-1 and AML-193 also showed elevated SERPINB10 expression. |
| Gene ID: | 5273 |
| Application Details | |
| Application Notes: | IF(IHC-P) 1:50-200 |
| Restrictions: | For Research Use only |
| Handling | |
| Format: | Liquid |
| | |
| Concentration: | 1 µg/µL |
| Concentration: Buffer: | • |
| | 1 μg/μL Aqueous buffered solution containing 0.01M TBS (pH 7.4) with 1 % BSA, 0.03 % Proclin300 and |
| Buffer: | 1 μg/μL Aqueous buffered solution containing 0.01M TBS (pH 7.4) with 1 % BSA, 0.03 % Proclin300 and 50 % Glycerol. |
| Buffer: Preservative: | 1 μg/μL Aqueous buffered solution containing 0.01M TBS (pH 7.4) with 1 % BSA, 0.03 % Proclin300 and 50 % Glycerol. ProClin This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be |
| Buffer: Preservative: Precaution of Use: | 1 μg/μL Aqueous buffered solution containing 0.01M TBS (pH 7.4) with 1 % BSA, 0.03 % Proclin300 and 50 % Glycerol. ProClin This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be handled by trained staff only. |