.-online.com antibodies

## Datasheet for ABIN1419718 anti-RNF187 antibody (HRP)



Overview	
Quantity:	100 µL
Target:	RNF187
Reactivity:	Human, Rat, Mouse
Host:	Rabbit
Clonality:	Polyclonal
Conjugate:	This RNF187 antibody is conjugated to HRP
Application:	Western Blotting (WB), Immunohistochemistry (Paraffin-embedded Sections) (IHC (p))
Product Details	
Immunogen:	KLH conjugated synthetic peptide derived from human RNF187
Isotype:	lgG
Cross-Reactivity:	Human, Mouse, Rat
Purification:	Purified by Protein A.
Target Details	
Target:	RNF187
Alternative Name:	RNF187 (RNF187 Products)
Background:	Synonyms: E3 ubiquitin-protein ligase RNF187, Protein RNF187, RACO-1, RING domain AP1 coactivator 1, RING finger protein 187, RN187_HUMAN, RNF187.
	Background: E3 ubiguitin-protein ligase that acts as a coactivator of JUN-mediated gene

Background: E3 ubiquitin-protein ligase that acts as a coactivator of JUN-mediated gene

activation in response to growth factor signaling via the MAP3K1 pathway, independently from

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/2 | Product datasheet for ABIN1419718 | 03/07/2024 | Copyright antibodies-online. All rights reserved.

Target Details		
	MAPK8.	
Gene ID:	149603	
Application Details		
Application Notes:	WB 1:300-5000	
	IHC-P 1:200-400	
Restrictions:	For Research Use only	
Handling		
Format:	Liquid	
Concentration:	1 μg/μL	
Buffer:	Aqueous buffered solution containing 0.01M TBS ( pH 7.4) with 1 % BSA, 0.03 % Proclin300 and 50 % Glycerol.	
Preservative:	ProClin	
Precaution of Use:	This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be handled by trained staff only.	
Handling Advice:	Do NOT add Sodium Azide! Use of Sodium Azide will inhibit enzyme activity of horseradish peroxidase.	
Storage:	-20 °C	
Storage Comment:	Store at -20°C. Aliquot into multiple vials to avoid repeated freeze-thaw cycles.	
Expiry Date:	12 months	