

## Datasheet for ABIN1423372

# anti-KCNE1-Like antibody (HRP)



#### Overview

Quantity:	100 μL
Target:	KCNE1-Like (KCNE1L)
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Conjugate:	This KCNE1-Like antibody is conjugated to HRP
Application:	Western Blotting (WB), Immunohistochemistry (Paraffin-embedded Sections) (IHC (p))
Product Details	
Immunogen:	KLH conjugated synthetic peptide derived from human KCNE1L
Isotype:	IgG
Cross-Reactivity:	Human, Mouse, Rat
Purification:	Purified by Protein A.
Target Details	
Target:	KCNE1-Like (KCNE1L)
Alternative Name:	KCNE1L (KCNE1L Products)
Background:	Synonyms: AMME syndrome candidate gene 2 protein, AMMECR2 protein, Cardiac voltage gated potassium channel accessory subunit 5, KCNE1 like, KCNE5, Mink, MinK like protein, Potassium voltage gated channel subfamily E member 1 like protein, Potassium voltage gated

channel, Isk related family, member 1 like, Voltage gated potassium channel accessory subunit

5 KCE1I	HUMAN.
J. NULTE	I IUIVIAIN.

Background: KCNE1L belongs to the potassium channel KCNE family which represent the most complex class of voltage-gated ion channels from both functional and structural standpoints. Their diverse functions include regulating neurotransmitter release, heart rate, insulin secretion, neuronal excitability, epithelial electrolyte transport, smooth muscle contraction, and cell volume.

Gene ID:

23630

### **Application Details**

Application Notes:	WB 1:300-5000
	IHC-P 1:200-400
Restrictions:	For Research Use only

# Handling

Format:	Liquid
Concentration:	1 μg/μL
Buffer:	Aqueous buffered solution containing 0.01M TBS (pH 7.4) with 1 % BSA, 0.03 % Proclin300 and 50 % Glycerol.
Preservative:	ProClin
Precaution of Use:	This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be handled by trained staff only.
Handling Advice:	Do NOT add Sodium Azide! Use of Sodium Azide will inhibit enzyme activity of horseradish peroxidase.
Storage:	-20 °C
Storage Comment:	Store at -20°C. Aliquot into multiple vials to avoid repeated freeze-thaw cycles.
Expiry Date:	12 months