antibodies -online.com







anti-BSX antibody (AA 101-200) (Alexa Fluor 555)



Overview		
Quantity:	100 μL	
Target:	BSX	
Binding Specificity:	AA 101-200	
Reactivity:	Mouse, Rat	
Host:	Rabbit	
Clonality:	Polyclonal	
Conjugate:	This BSX antibody is conjugated to Alexa Fluor 555	
Application:	Western Blotting (WB), Flow Cytometry (FACS), Immunofluorescence (Cultured Cells) (IF (cc)), Immunofluorescence (Paraffin-embedded Sections) (IF (p))	
Product Details		
Immunogen:	KLH conjugated synthetic peptide derived from human BSX	
Isotype:	IgG	
Cross-Reactivity:	Mouse, Rat	
Predicted Reactivity:	Human,Dog,Cow,Sheep,Pig,Horse,Rabbit	
Purification:	Purified by Protein A.	
Target Details		
Target:	BSX	
Alternative Name:	BSX (BSX Products)	

Target Details

•				
Background:	Synonyms: Brain specic homeobox, Brain specic homeobox protein homolog, BSH, BSX1,			
	BSH_HUMAN.			
	Background: BSX(Brain-specific homeobox protein homolog) may play a role in the			
	determination and function of cell types in the brain. In brain, it is restricted to a few specific			
	developing brain structures such as pineal gland, telencephalic septum, hypothalamic pre-			
	mammillary body and arcuate nucleus. Bsx might be considered an important molecular			
	marker for early embryonic stages of epiphysis development. It belongs to the distal-less			
	homeobox family and contains 1 homeobox DNA-binding domain.			
Gene ID:	390259			
Pathways:	Feeding Behaviour			
Application Details				
Application Notes:	FCM 1:20-100			
	IF(IHC-P) 1:50-200			
	IF(IHC-F) 1:50-200			
	IF(ICC) 1:50-200			
Restrictions:	For Research Use only			
Handling				
Format:	Liquid			
Concentration:	1 μg/μL			
Buffer:	Aqueous buffered solution containing 0.01M TBS (pH 7.4) with 1 % BSA, 0.03 % Proclin300 and			
	50 % Glycerol.			
Preservative:	ProClin			
Precaution of Use:	This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be			
	handled by trained staff only.			
Storage:	-20 °C			
Storage Comment:	Store at -20°C. Aliquot into multiple vials to avoid repeated freeze-thaw cycles.			
Expiry Date:	12 months			