antibodies.com

Datasheet for ABIN1699202 anti-Destrin antibody (AA 75-165) (AbBy Fluor® 647)



Overview

Quantity:	100 µL	
Target:	Destrin (DSTN)	
Binding Specificity:	AA 75-165	
Reactivity:	Human	
Host:	Rabbit	
Clonality:	Polyclonal	
Conjugate:	This Destrin antibody is conjugated to AbBy Fluor® 647	
Application:	Western Blotting (WB), Immunofluorescence (Cultured Cells) (IF (cc)), Immunofluorescence (Paraffin-embedded Sections) (IF (p))	

Product Details

Immunogen:	KLH conjugated synthetic peptide derived from human Destrin	
Isotype:	lgG	
Predicted Reactivity:	Human,Mouse,Rat,Dog,Cow,Horse,Rabbit,Monkey	
Purification:	Purified by Protein A.	
Target Details		

Target:	Destrin (DSTN)	
Alternative Name:	Destrin (DSTN Products)	
Background:	Synonyms: 2610043P17Rik, ACTDP, Actin depolymerizing factor, Actin-depolymerizing factor,	

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/3 | Product datasheet for ABIN1699202 | 06/07/2024 | Copyright antibodies-online. All rights reserved.

ADF, AU042046, bA462D18.2, corn1, DEST_HUMAN, Destrin actin depolymerizing factor,
Destrin, DSN, Dstn, Sid 23, sid23p.

Background: Actin-depolymerizing factor (ADF), also known as destrin, is a member of the ADF/Cofilin/destrin superfamily that has the ability to rapidly depolymerize F-Actin in a stoichiometric manner. The Actin-depolymerizing activity of ADF is reversibly controlled by changes in KCl concentration but is insensitive to calcium concentration. ADF depolymerizes F-Actin by interacting directly with F-Actin protomers. ADF shares 71 % sequence homology with Cofilin, however the two proteins differ in their interaction with Actin. The difference in the function of ADF and Cofilin results from the subtle difference in their amino acid sequence rather than possible differences in posttranslational modifications. As a result of different cleavage sites on ADF and Cofilin, the proteins differ in their overall tertiary folds. Sensitivity to polyphosphoinositides may be a common feature in vitro among Actin-binding proteins such as ADF and Cofilin that can bind to G-Actin and regulate the state of Actin polymerization. ADF and Cofilin are Actin-depolymerizing proteins whose activities are possibly regulated by their phosphorylation/dephosphorylation.

Gene ID:

11034

Application Details

Application Notes: IF(IHC-P) 1:50-200 IF(IHC-F) 1:50-200 IF(ICC) 1:50-200 Restrictions: For Research Use only		
IF(ICC) 1:50-200		
Restrictions: For Research Use only		
Handling		
Format: Liquid		
Concentration: 1 µg/µL	1 μg/μL	
Buffer: Aqueous buffered solution containing 0.01M TBS (pH 7.4) with 1 % BSA	, 0.03 % Proclin300 and	
50 % Glycerol.		
Preservative: ProClin		
Precaution of Use: This product contains ProClin: a POISONOUS AND HAZARDOUS SUBST.	ANCE, which should be	
handled by trained staff only.		
Storage: -20 °C		

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 2/3 | Product datasheet for ABIN1699202 | 06/07/2024 | Copyright antibodies-online. All rights reserved.

lond	lina
land	
10110	

Storage Comment:

Store at -20°C. Aliquot into multiple vials to avoid repeated freeze-thaw cycles.

Expiry Date:

12 months

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 3/3 | Product datasheet for ABIN1699202 | 06/07/2024 | Copyright antibodies-online. All rights reserved.