



[Go to Product page](#)

Datasheet for ABIN1705697
anti-POLH antibody (AA 51-160) (Cy5.5)

Overview

Quantity:	100 µL
Target:	POLH
Binding Specificity:	AA 51-160
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Conjugate:	This POLH antibody is conjugated to Cy5.5
Application:	Western Blotting (WB), Immunofluorescence (Cultured Cells) (IF (cc)), Immunofluorescence (Paraffin-embedded Sections) (IF (p))

Product Details

Immunogen:	KLH conjugated synthetic peptide derived from human DNA polymerase eta
Isotype:	IgG
Predicted Reactivity:	Human,Mouse,Rat,Rabbit
Purification:	Purified by Protein A.

Target Details

Target:	POLH
Alternative Name:	DNA polymerase eta (POLH Products)
Background:	Synonyms: DNA polymerase eta, FLJ16395, FLJ21978, POLH, POLH_HUMAN, polymerase DNA

Target Details

directed eta, RAD30, RAD30 homolog A, RAD30A, Xeroderma pigmentosum variant type protein, XP V, XPV.

Background: DNA polymerase specifically involved in DNA repair. Plays an important role in translesion synthesis, where the normal high fidelity DNA polymerases cannot proceed and DNA synthesis stalls. Plays an important role in the repair of UV-induced pyrimidine dimers.

Depending on the context, it inserts the correct base, but causes frequent base transitions and transversions. May play a role in hypermutation at immunoglobulin genes. Forms a Schiff base with 5'-deoxyribose phosphate at abasic sites, but does not have lyase activity. Targets POLI to replication foci.

Gene ID: 5429

Pathways: [DNA Damage Repair](#)

Application Details

Application Notes: IF(IHC-P) 1:50-200
IF(IHC-F) 1:50-200
IF(ICC) 1:50-200

Restrictions: For Research Use only

Handling

Format: Liquid

Concentration: 1 µg/µL

Buffer: Aqueous buffered solution containing 0.01M TBS (pH 7.4) with 1 % BSA, 0.03 % Proclin300 and 50 % Glycerol.

Preservative: ProClin

Precaution of Use: This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be handled by trained staff only.

Storage: -20 °C

Storage Comment: Store at -20°C. Aliquot into multiple vials to avoid repeated freeze-thaw cycles.

Expiry Date: 12 months