# antibodies -online.com







## anti-Ephrin A5 antibody (AA 21-203)

3 Images



Go to Product page

$\sim$	
( )\/△	rview
$\cup$	1 410 44

Quantity:	100 μL
Target:	Ephrin A5 (EFNA5)
Binding Specificity:	AA 21-203
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Conjugate:	This Ephrin A5 antibody is un-conjugated
Application:	Western Blotting (WB), Immunohistochemistry (IHC), Immunocytochemistry (ICC)

## **Product Details**


Immunogen:	EFNA5 (Gln21-Asn203)
Isotype:	IgG
Specificity:	The antibody is a rabbit polyclonal antibody raised against EFNA5. It has been selected for its ability to recognize EFNA5 in immunohistochemical staining and western blotting.
Purification:	Antigen-specific affinity chromatography

## **Target Details**

Target:	Ephrin A5 (EFNA5)
Alternative Name:	Ephrin A5 (EFNA5 Products)
Background:	Alternative Names: AF1, EFL5, EPLG7, GLC1M, LERK7, RAGS, AL-1, EPH-related receptor


## **Target Details**

Immunohistochemistry in formalin fixed frozen section: 1:50-500 Immunohisto paraffin section: 1:10-100 Enzyme-linked Immunosorbent Assay: 1:100-1:5000 working dilutions must be determined by end user.  Comment: The thermal stability is described by the loss rate. The loss rate was determined by thermal degradation test, that is, incubate the protein at 37&degC for 48h, and no degradation and precipitation were observed. The loss rate is less than 5% within date under appropriate storage condition.  Restrictions: For Research Use only  Handling  Format: Liquid  Concentration: Lot specific  Buffer: PBS, pH 7.4, containing 0.02 % Sodium azide, 50 % glycerol.  Preservative: Sodium azide  Precaution of Use: WARNING: Reagents contain sodium azide. Sodium azide is very toxic if ingested Avoid contact with skin, eyes, or clothing. Wear eye or face protection when handle eye contact occurs, wash with copious amounts of water. If ingested or inhaled, on physician immediately. Sodium azide yields toxic hydrazoic acid under acidic contains contain sodium azide social under acidic contains contains and the physician immediately. Sodium azide yields toxic hydrazoic acid under acidic contains contains acide yields toxic hydrazoic acid under acidic contains acidic hydrazoic acid under acidic	rarget Details	
Application Details  - Western blotting: 1:50-400 Immunocytochemistry in formalin fixed cells: 1:50-5 - Immunohistochemistry in formalin fixed frozen section: 1:50-500 Immunohistochemistry in formalin fixed frozen section: 1:50-500 Immunohistochemistry in formalin fixed cells: 1:50-5 - Immunohistochemistry in formalin fixed frozen section: 1:50-500 Immunohistochemistry in formalin fixed cells: 1:50-5 - Immunohistochemistry in formalin fixed cells: 1:50-500 Immunohistry in formalin fixed cells: 1:50-50		tyrosine kinase ligand 7
Application Notes:  - Western blotting: 1:50-400 Immunocytochemistry in formalin fixed cells: 1:50-5 Immunohistochemistry in formalin fixed frozen section: 1:50-500 Immunohistochemistry in f	Pathways:	RTK Signaling
Immunohistochemistry in formalin fixed frozen section: 1:50-500 Immunohistochemistry in formalin fixed frozen section: 1:50-500 Immunohistochemistry in formalin fixed frozen section: 1:50-500 working dilutions must be determined by end user.  Comment: The thermal stability is described by the loss rate. The loss rate was determined by the thermal degradation test, that is, incubate the protein at 37&degC for 48h, and no degradation and precipitation were observed. The loss rate is less than 5% within date under appropriate storage condition.  Restrictions: For Research Use only  Handling  Format: Liquid  Concentration: Lot specific  Buffer: PBS, pH 7.4, containing 0.02 % Sodium azide, 50 % glycerol.  Preservative: Sodium azide  Precaution of Use: WARNING: Reagents contain sodium azide. Sodium azide is very toxic if ingested Avoid contact with skin, eyes, or clothing. Wear eye or face protection when handle eye contact occurs, wash with copious amounts of water. If ingested or inhaled, or physician immediately. Sodium azide yields toxic hydrazoic acid under acidic con azide-containing compounds in running water before discarding to avoid accumu potentially explosive deposits in lead or copper plumbing.  Handling Advice: Avoid repeated freeze-thaw cycles.  Storage: 4 °C	Application Details	
thermal degradation test, that is, incubate the protein at 37&degC for 48h, and no degradation and precipitation were observed. The loss rate is less than 5% within date under appropriate storage condition.  Restrictions: For Research Use only  Handling  Format: Liquid  Concentration: Lot specific  Buffer: PBS, pH 7.4, containing 0.02 % Sodium azide, 50 % glycerol.  Preservative: Sodium azide  Precaution of Use: WARNING: Reagents contain sodium azide. Sodium azide is very toxic if ingested Avoid contact with skin, eyes, or clothing. Wear eye or face protection when handl eye contact occurs, wash with copious amounts of water. If ingested or inhaled, ophysician immediately. Sodium azide yields toxic hydrazoic acid under acidic contazide-containing compounds in running water before discarding to avoid accumulative potentially explosive deposits in lead or copper plumbing.  Handling Advice: Avoid repeated freeze-thaw cycles.  Storage: 4 °C	Application Notes:	<ul> <li>Western blotting: 1:50-400 Immunocytochemistry in formalin fixed cells: 1:50-500         Immunohistochemistry in formalin fixed frozen section: 1:50-500 Immunohistochemistry in paraffin section: 1:10-100 Enzyme-linked Immunosorbent Assay: 1:100-1:5000 Optimal working dilutions must be determined by end user.     </li> </ul>
Format:  Liquid  Concentration:  Lot specific  Buffer:  PBS, pH 7.4, containing 0.02 % Sodium azide, 50 % glycerol.  Preservative:  Sodium azide  Precaution of Use:  WARNING: Reagents contain sodium azide. Sodium azide is very toxic if ingested Avoid contact with skin, eyes, or clothing. Wear eye or face protection when handle eye contact occurs, wash with copious amounts of water. If ingested or inhaled, or physician immediately. Sodium azide yields toxic hydrazoic acid under acidic contazide-containing compounds in running water before discarding to avoid accumulation potentially explosive deposits in lead or copper plumbing.  Handling Advice:  Avoid repeated freeze-thaw cycles.  Storage:  4 °C	Comment:	The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37&degC for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.
Format:  Liquid  Concentration:  Lot specific  Buffer:  PBS, pH 7.4, containing 0.02 % Sodium azide, 50 % glycerol.  Preservative:  Sodium azide  WARNING: Reagents contain sodium azide. Sodium azide is very toxic if ingested Avoid contact with skin, eyes, or clothing. Wear eye or face protection when handle eye contact occurs, wash with copious amounts of water. If ingested or inhaled, or physician immediately. Sodium azide yields toxic hydrazoic acid under acidic contazide-containing compounds in running water before discarding to avoid accumu potentially explosive deposits in lead or copper plumbing.  Handling Advice:  Avoid repeated freeze-thaw cycles.  Storage:  4 °C	Restrictions:	For Research Use only
Concentration:  Lot specific  Buffer:  PBS, pH 7.4, containing 0.02 % Sodium azide, 50 % glycerol.  Preservative:  Sodium azide  WARNING: Reagents contain sodium azide. Sodium azide is very toxic if ingested Avoid contact with skin, eyes, or clothing. Wear eye or face protection when handle eye contact occurs, wash with copious amounts of water. If ingested or inhaled, or physician immediately. Sodium azide yields toxic hydrazoic acid under acidic contazide-containing compounds in running water before discarding to avoid accumu potentially explosive deposits in lead or copper plumbing.  Handling Advice:  Avoid repeated freeze-thaw cycles.  Storage:  4 °C	Handling	
Buffer: PBS, pH 7.4, containing 0.02 % Sodium azide, 50 % glycerol.  Preservative: Sodium azide  WARNING: Reagents contain sodium azide. Sodium azide is very toxic if ingested Avoid contact with skin, eyes, or clothing. Wear eye or face protection when handle eye contact occurs, wash with copious amounts of water. If ingested or inhaled, or physician immediately. Sodium azide yields toxic hydrazoic acid under acidic contazide-containing compounds in running water before discarding to avoid accumulation potentially explosive deposits in lead or copper plumbing.  Handling Advice: Avoid repeated freeze-thaw cycles.  Storage: 4 °C	Format:	Liquid
Precaution of Use:  WARNING: Reagents contain sodium azide. Sodium azide is very toxic if ingested Avoid contact with skin, eyes, or clothing. Wear eye or face protection when handle eye contact occurs, wash with copious amounts of water. If ingested or inhaled, or physician immediately. Sodium azide yields toxic hydrazoic acid under acidic contazide-containing compounds in running water before discarding to avoid accumu potentially explosive deposits in lead or copper plumbing.  Handling Advice:  Avoid repeated freeze-thaw cycles.  Storage:  4 °C	Concentration:	Lot specific
Precaution of Use:  WARNING: Reagents contain sodium azide. Sodium azide is very toxic if ingested Avoid contact with skin, eyes, or clothing. Wear eye or face protection when handle eye contact occurs, wash with copious amounts of water. If ingested or inhaled, or physician immediately. Sodium azide yields toxic hydrazoic acid under acidic contazide-containing compounds in running water before discarding to avoid accumu potentially explosive deposits in lead or copper plumbing.  Handling Advice:  Avoid repeated freeze-thaw cycles.  Storage:  4 °C	Buffer:	PBS, pH 7.4, containing 0.02 % Sodium azide, 50 % glycerol.
Avoid contact with skin, eyes, or clothing. Wear eye or face protection when handle eye contact occurs, wash with copious amounts of water. If ingested or inhaled, or physician immediately. Sodium azide yields toxic hydrazoic acid under acidic contazide-containing compounds in running water before discarding to avoid accumu potentially explosive deposits in lead or copper plumbing.  Handling Advice:  Avoid repeated freeze-thaw cycles.  Storage:  4 °C	Preservative:	Sodium azide
Storage: 4 °C	Precaution of Use:	WARNING: Reagents contain sodium azide. Sodium azide is very toxic if ingested or inhaled.  Avoid contact with skin, eyes, or clothing. Wear eye or face protection when handling. If skin or eye contact occurs, wash with copious amounts of water. If ingested or inhaled, contact a physician immediately. Sodium azide yields toxic hydrazoic acid under acidic conditions. Dilute azide-containing compounds in running water before discarding to avoid accumulation of potentially explosive deposits in lead or copper plumbing.
	Handling Advice:	Avoid repeated freeze-thaw cycles.
Storage Comment: Store at 2-8 °C for one month. Aliquot and store at -80 °C for 12 months.	Storage:	4 °C
	Storage Comment:	Store at 2-8 °C for one month. Aliquot and store at -80 °C for 12 months.
Expiry Date: 12 months	Expiry Date:	12 months




## **Western Blotting**

Image 1.



## **Immunohistochemistry**

**Image 2.** Used in DAB staining on fromalin fixed paraffinembedded kidney tissue



## **Immunohistochemistry**

Image 3. DAB staining on IHC-P; Samples: Human Liver Tissue