

Datasheet for ABIN1882585

AnaTag™ HiLyte Fluor™ 750 Protein Labeling Kit



[Go to Product page](#)

Overview

Quantity:	1 kit
Conjugate:	HiLyte Fluor 750
Application:	Labeling (Lbl)

Product Details

Detection Method:	Fluorometric
-------------------	--------------

Target Details

Background:	HiLyte Fluor™ 750 acid, SE is the longest-wavelength amine-reactive HiLyte Fluor™ dye currently available. Its fluorescence emission maximum at 778 nm is well separated from commonly used far-red fluorophores.
-------------	---

Application Details

Comment:	Labeling kit
Restrictions:	For Research Use only

Handling

Handling Advice:	Keep Component A away from light and protected from moisture.
Storage:	4 °C