

Datasheet for ABIN1882599

AnaTag™ 5 - FITC Protein Labeling Kit



[Go to Product page](#)

Overview

Quantity:	1 kit
Application:	Labeling (Lbl)

Product Details

Detection Method:	Fluorometric
-------------------	--------------

Target Details

Background:	Fluorescein isothiocyanate (FITC) remains one of the most popular fluorescent labeling reagents, probably due to the low cost of the material. 5-FITC and 6-FITC have very similar absorption and fluorescence spectra.
-------------	---

Application Details

Comment:	Labeling kit
Restrictions:	For Research Use only

Handling

Handling Advice:	Keep Component A away from light and protected from moisture.
Storage:	4 °C