antibodies -online.com





Go to Product page

Datasheet for ABIN1925744

anti-PPM1E antibody (C-Term) (FITC)

| Overview | |
|----------------------|--------------------------------------------------------------------------------------|
| Quantity: | 200 μL |
| Target: | PPM1E |
| Binding Specificity: | C-Term |
| Reactivity: | Human |
| Host: | Rabbit |
| Clonality: | Polyclonal |
| Conjugate: | This PPM1E antibody is conjugated to FITC |
| Application: | ELISA, Immunohistochemistry (IHC) |
| Product Details | |
| Immunogen: | KLH conjugated synthetic peptide selected from the C-terminal region of human PPM1E. |
| | Type of Immunogen: Synthetic peptide |
| Isotype: | IgG |
| Purification: | Affinity purified |
| Target Details | |
| Target: | PPM1E |
| Alternative Name: | PPM1E (PPM1E Products) |
| Background: | Name/Gene ID: PPM1E |

Target Details

| rarget Details | |
|---------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------|
| | Synonyms: PPM1E, CaMKN, CaMKP-N, CaMKP-nucleus, Partner of PIX 1, Partner of PIXA, PP2CH, POPX1, KIAA1072, Partner of PIX-alpha, Protein phosphatase 1E |
| Gene ID: | 22843 |
| | |
| Application Details | |
| Application Notes: | Approved: ELISA, IHC |
| | Usage: The applications listed have been tested for the unconjugated form of this product. Other forms have not been tested. |
| Comment: | Target Species of Antibody: Human |
| Restrictions: | For Research Use only |
| Handling | |
| Format: | Liquid |
| Concentration: | Lot specific |
| Buffer: | PBS, pH 7.2, 0.09 % sodium azide. |
| Preservative: | Sodium azide |
| Precaution of Use: | This product contains Sodium azide: a POISONOUS AND HAZARDOUS SUBSTANCE which should be handled by trained staff only. |
| Handling Advice: | Aliquot to avoid repeated freezing and thawing. |
| Storage: | 4 °C,-20 °C |
| Storage Comment: | Short term: store at 4°C. Long term: aliquot and store -20°C for up to 6 months. Avoid freeze-thaw cycles. Protect from light. |
| Expiry Date: | 6 months |