.-online.com antibodies

Datasheet for ABIN1981732 TYK2 ELISA Kit

3 Images



Overview

Quantity:	96 tests
Target:	TYK2
Binding Specificity:	pTyr
Reactivity:	Human
Method Type:	Sandwich ELISA
Application:	ELISA

Product Details

Purpose:	Human Phosphotyrosine TYK2 ELISA Kit. This assay semi-quantitatively measures phosphotyrosine TYK2 in lysate samples.				
Sample Type:	Cell Lysate, Tissue Lysate				
Analytical Method:	Semi-Quantitative				
Detection Method:	Colorimetric				
Specificity:	The antibody pair provided in this kit recognizes Human Tyrosine-Phosphorylated-TYK2.				
Characteristics:	 Rapidly measure phosphorylated protein in lysates Screen numerous different cell lysates without performing a Western Blot analysis Minimal hands-on time, convenient, and non-radioactive material 				
Components:	 Pre-Coated 96-well Strip Microplate Wash Buffer Biotinylated Anti-Phosphotyrosine Antibody Stop Solution 				

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/6 | Product datasheet for ABIN1981732 | 09/11/2023 | Copyright antibodies-online. All rights reserved.

	 Assay Diluent(s) Positive Control Sample Lysis Buffer Streptavidin-Conjugated HRP TMB One-Step Substrate
Material not included:	 Distilled or deionized water 100 mL and 1 liter graduated cylinders Tubes to prepare sample dilutions Protease and Phosphatase inhibitors Precision pipettes to deliver 2 µL to 1 mL volumes Adjustable 1-25 mL pipettes for reagent preparation Benchtop rocker or shaker Microplate reader capable of measuring absorbance at 450 nm

Target Details

Target:	TYK2
Alternative Name:	TYK2 (TYK2 Products)
Gene ID:	7297
UniProt:	P29597
Pathways:	JAK-STAT Signaling, Hepatitis C

Application Details

Sample Volume:	100 μL
Plate:	Pre-coated
Protocol:	1. Prepare all reagents and samples as instructed in the manual.
	2. Add 100 μ L of sample or positive control to each well.
	3. Incubate 2.5 h at RT or O/N at 4 °C.
	4. Add 100 μ L of prepared primary antibody to each well.
	5. Incubate 1 h at RT.
	6. Add 100 μ L of prepared 1X HRP-Streptavidin to each well.
	7. Incubate 1 h at RT.
	8. Add 100 μ L of TMB One-Step Substrate Reagent to each well.
	9. Incubate 30 min at RT.
	10. Add 50 μL of Stop Solution to each well.
	11. Read at 450 nm immediately.

Application	Details
-------------	---------

Reagent Preparation:

2. Item E, Assay Diluent should be diluted 5-fold with deionized or distilled water before use. 3. Preparation of Positive Control: Briefly spin the Positive Control vial of Item K. Add 250 μ L 1x Assay Diluent (Item E, Assay Diluent should be diluted 5-fold with deionized or distilled water before use) into Item K vial to prepare a Positive Control Stock Solution. Dissolve the powder thoroughly by a gentle mix. Add 150 μ L prepared Positive Control Stock Solution from the vial of Item K, into a tube with 300 μ L 1x Assay Diluent to prepare P-1 (See i. Positive control of part IX.for a typical result). Pipette 300 μ L 1x Assay Diluent into each tube. Transfer 150 μ L prepared P-1 into a tube with 300 μ L 1x Assay Diluent to produce a dilution series (shown below). Mix each tube thoroughly before the next transfer. 1x Assay Diluent serves as the background. Phospho-TYK2 ELISA 6

1. Bring all reagents and samples to room temperature (18 - 25 °C) before use.

4. If the Wash Concentrate (20x) (Item B) contains visible crystals, warm to room temperature and mix gently until dissolved. Dilute 20 mL of Wash Buffer Concentrate into deionized or distilled water to yield 400 mL of 1x Wash Buffer.

5. Briefly spin the biotinylated antibody (Item C) before use. Add 100 μ L of 1x Assay Diluent into the vial to prepare a biotinylated anti-phosphotyrosine antibody concentrate. Pipette up and down to mix gently (the concentrate can be stored at 4 °C for 5 days or at -80 °C for one month). The biotinylated phosphotyrosine antibody should be diluted 80x with 1x Assay Diluent and used in step 4 of Part VII Assay Procedure.

6. Briefly spin the HRP-Streptavidin concentrate vial (Item G) and pipette up and down to mix gently before use. HRP-Streptavidin concentrate should be diluted 600 fold with 1x Assay Diluent. For example: Briefly spin the vial (Item G) and pipette up and down to mix gently . Add 20 μL of HRP-Streptavidin concentrate into a tube with 12 mL 1x Assay Diluent to prepare a 600-fold diluted HRP-Streptavidin solution (don't store the diluted solution for next P-1 P-2 P-3 P-4 P-5 0 150 μL 150μl Positive Control Stock Solution 300 μL 1x Assay Diluent 150μl 150 μL 150 μL Phospho-TYK2 ELISA 7 day use). Mix well.

7. Cell Lysate Buffer should be diluted 2-fold with deionized or distilled water before use (recommend to add protease and phosphatase inhibitors). VII.

Sample Preparation: Cell lysates - Rinse cells with PBS, making sure to remove any remaining PBS before adding the lysis buffer. Solubilize cells at 4 x 107 cells/mL in 1x Lysis Buffer (we recommend adding protease and phosphatase inhibitors to lysis buffer prior to sample preparation). Pipette up and down to resuspend and incubate the lysates with shaking at 2 - 8° C for 30 minutes. Microcentrifuge at 13,000 rpm for 10 minutes at 2 - 8° C, and transfer the supernates into a clean test tube. Lysates should be used immediately or aliquoted and stored at -70 °C. Avoid repeated freeze-thaw cycles. Thawed lysates should be kept on ice prior to use.

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 3/6 | Product datasheet for ABIN1981732 | 09/11/2023 | Copyright antibodies-online. All rights reserved.

Application Details

	For the initial experiment, we recommend to do a serial dilution testing such as 5-fold and 100-
	fold dilution for your cell lysates with Assay Diluent (Item E) before use.
	Note: The fold dilution of sample used depends on the abundance of phosphorylated proteins
	and should be determined empiricallys. Phospho-TYK2 ELISA 5 More of the sample can be
	used if signals are too weak. If signals are too strong, the sample can be diluted further.
	Cell lysate buffer should be diluted 2-fold with deionized or distilled water before use
	(recommend to add protease and phosphatase inhibitors).
Assay Procedure:	1. Bring all reagents to room temperature (18 - 25 °C) before use. It is recommended that all
	samples or Positive Control should be run at least in duplicate.
	2. Add 100 μ L of each sample or positive control into appropriate wells. Cover well with plate
	holder and incubate for 2.5 hours at room temperature or over night at 4 °C with shaking.
	3. Discard the solution and wash 4 times with 1x Wash Solution. Wash by filling each well with
	Wash Buffer (300 μ L) using a multi-channel pipette or autowasher. Complete removal of liquid
	at each step is essential to good performance. After the last wash, remove any remaining Wash
	Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
	4. Add 100 µL of prepared 1X biotinylated anti-phosphotyrosine antibody (Reagent Preparation
	step 5) to each well. Incubate for 1 hour at room temperature with shaking.
	5. Discard the solution. Repeat the wash as in step3. Phospho-TYK2 ELISA 8
	6. Add 100 µL of prepared 1X HRP-Streptavidin solution (see Reagent Preparation step 6) to
	each well. Incubate for 45 minutes at room temperature with shaking.
	7. Discard the solution. Repeat the wash as in step3.
	8. Add 100 µL of TMB One-Step Substrate Reagent (Item H) to each well. Incubate for 30
	minutes at room temperature in the dark with shaking.
	9. Add 50 μL of Stop Solution (Item I) to each well. Read at 450 nm immediately.
Calculation of Results:	ELISA data analysis: Average the duplicate readings for each sample or positive control then
	subtract the average blank optical density.
	i. Positive Control Jurkat cells were treated with Pervanadate at 37 °C for 10 min. Solubilize
	cells at 4 x 107 cells/mL in lysis buffer. Serial dilutions of lysates were analyzed in this ELISA.
	Please see step 3 of Part VI Reagent Preparation for detail. 0 0.5 1 1.5 2 2.5 3 0D =4 50 n m P-1
	P-2 P-3 P-4 P-5 Positive control dilution series Assay Diluent Phospho-TYK2 ELISA 10
	ii. Pervanadate Stimulation of Jurkat Cell Line Jurkat cells were treated or untreated with
	Pervanadate for 10 min at 37 °C. Cell lysates were analyzed using this phosphoELISA: 0 1 2 3 $$
	Untreated Pervanadate OD = 4 50 n m Phospho-TYK2 ELISA 11 X
Restrictions [.]	For Research Use only

Restrictions:


For Research Use only

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 4/6 | Product datasheet for ABIN1981732 | 09/11/2023 | Copyright antibodies-online. All rights reserved.

Handling

Handling Advice:	Avoid repeated freeze- thaw cycles.
Storage:	-20 °C
Storage Comment:	Upon receipt, the kit should be stored at -20 °C. Please use within 6 months from the date of
	shipment. After initial use, Wash Buffer Concentrate (Item B), Assay Diluent (Item E), TMB One-
	Step Substrate Reagent (Item H), HRP-Streptavidin (Item G), Stop Solution (Item I) and Cell
	Lysate Buffer (Item J) should be stored at 4 $^\circ$ C to avoid repeated freeze-thaw cycles. Return
	unused wells to the pouch containing desiccant pack, reseal along entire edge and store at -20
	°C. Reconstituted Positive Control (Item K) should be stored at -70 °C.
Expiry Date:	6 months

Images



Positive control dilution series

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 5/6 | Product datasheet for ABIN1981732 | 09/11/2023 | Copyright antibodies-online. All rights reserved.

150µl Positive Control Stock Solution + 300 µl	1x	150µl	150 µl	150 µl	150 µl	
Assay Diluent	\mathbb{N}	Da	Da	D a	D	
	6	6	$\overline{\bigcirc}$	$\overline{\bigcirc}$	6	\bigcirc
	P-1	P-2	P-3	P-4	P-5	0

Image 3. This picture shows the preparation of the positive control.

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 6/6 | Product datasheet for ABIN1981732 | 09/11/2023 | Copyright antibodies-online. All rights reserved.