# .-online.com antibodies

# Datasheet for ABIN2472083 Selectin E/CD62e ELISA Kit

Image



### Overview

Quantity:	96 tests
Target:	Selectin E/CD62e (SELE)
Binding Specificity:	Soluble
Reactivity:	Human
Method Type:	Sandwich ELISA
Detection Range:	93.50-6000 pg/mL
Minimum Detection Limit:	93.50 pg/mL
Application:	ELISA

# Product Details

Purpose:	This assay employs the quantitative sandwich enzyme immunoassay technique for quantitative detection.
Sample Type:	Serum, Plasma, Cell Culture Supernatant
Analytical Method:	Quantitative
Detection Method:	Colorimetric
Specificity:	SE-Selectin, Human
Sensitivity:	1.94 pg/mL
Components:	Plate, Standard, Diluent
Material not included:	Microplate reader capable of measuring absorbance at 450 nm, with correction wavelength set at 570 nm or 630 nm.

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/8 | Product datasheet for ABIN2472083 | 09/11/2023 | Copyright antibodies-online. All rights reserved. Pipettes and pipette tips.

50  $\mu L$  to 300  $\mu L$  adjustable multichannel micropipette with disposable tips.

Multichannel micropipette reservoir.

Beakers, flasks, cylinders necessary for preparation of reagents.

Deionized or distilled water.

Polypropylene test tubes for dilution.

# Target Details

Target:	Selectin E/CD62e (SELE)
Alternative Name:	E-Selectin (SELE Products)
Background:	E-selectin, also known as CD62 antigen-like family member E (CD62E), endothelial-leukocyte adhesion molecule 1 (ELAM-1), or leukocyte-endothelial cell adhesion molecule 2 (LECAM2), is a 115 kDa glycoprotein expressed only after activation by inflammatory cytokines (IL-1 $\beta$ or TNF- $\alpha$ ) or endotoxins on endothelial cells.
Gene ID:	6401
NCBI Accession:	NM_000450
UniProt:	P16581
Pathways:	Thromboxane A2 Receptor Signaling

## Application Details

Sample Volume:	100 µL
Assay Time:	3 - 4 h
Plate:	Pre-coated
Protocol:	1. Prepare all reagents and standards as directed.
	2. Add 50 µL Assay Buffer to each well.
	3. Add 50 $\mu$ L Standard or sample per well within 15 minutes.
	4. Add 50 µL Detect Antibody to each well.
	5. Incubate for 2 hours at RT.
	6. Aspirate and wash 6 times.
	7. Add 100 µL Streptavidin-HRP to each well.
	8. Incubate for 45 min at RT.
	9. Aspirate and wash 6 times.

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 2/8 | Product datasheet for ABIN2472083 | 09/11/2023 | Copyright antibodies-online. All rights reserved. 10. Add 100 µL Substrate Solution to each well. Incubate for 10 - 30 minutes at RT. Protect from light.
11. Add 100 µL Stop Solution to each well.

12. Read at 450 nm within 30 minutes. Correction 570 or 630 nm

Reagent Preparation: If crystals form in the Buffer Concentrates, warm and gently stir them until completely dissolved. Washing Buffer (1x) Pour entire contents (50 mL) of the Washing Buffer (20x) into a clean 1000 mL graduated cylinder. Bring to final volume of 1000 mL with pure or deionized water. Mix gently to avoid foaming. Transfer to a clean wash bottle and store at 2 to 25 °C. Washing Buffer (1x) is stable for 30 days. Assay Buffer (1x) Pour the entire contents (5 mL) of the Assay Buffer (10x) into a clean 100 mL graduated cylinder. Bring to final volume of 50 mL with distilled water. Mix gently to avoid foaming. Store at 2 to 8 °C. Assay Buffer (1x) is stable for 30 days. Detect Antibody Mix well prior to making dilutions. Make a 1:100 dilution of the concentrated Detect Antibody solution with Assay Buffer (1x) in a clean plastic tube as needed. The diluted Detect Antibody should be used within 30 minutes after dilution. Streptavidin-HRP Mix well prior to making dilutions. Make a 1:100 dilution of the concentrated Streptavidin-HRP solution with Assay Buffer (1x) in a clean plastic tube as needed. The diluted Streptavidin-HRP should be used within 30 minutes after dilution. Sample Dilution If your samples have high sE-Selectin content, dilute serum/plasma samples with Assay Buffer (1x). For cell culture supernates, dilute with cell culture medium. Human sE-Selectin Standard Reconstitute Human sE-Selectin Standard by addition of distilled water. Reconstitution volume is stated on the label of the standard vial. Swirl or mix gently to insure complete and homogeneous solubilization (concentration of reconstituted standard = 24000 pg/mL). Allow the standard to reconstitute for 10 - 30 minutes. Mix well prior to making dilutions. Use polypropylene tubes. Human sE-Selectin ELISA Kit 7 For serum/plasma samples, mixing concentrated human sE-Selectin standard (150 µL) with 150 µL of Assay Buffer (1x) creates the high standard (12000 pg/ mL). Pipette 150 µL of Assay Buffer (1x) into each tube. Use the high standard to produce a 1:1 dilution series . Mix each tube thoroughly before the next transfer. Assay Buffer (1x) serves as the zero standard (0 pg/mL). For cell culture supernates, mixing concentrated human sE-Selectin standard (150 µL) with 150 µL of cell culture medium creates the high standard (12000 pg/mL). Pipette 150 µL of cell culture medium into each tube. Use the high standard to produce a 1:1 dilution series. Mix each tube thoroughly before the next transfer. Cell culture medium serves as the zero standard (0 pg/ mL). Sample Collection: Cell

Culture Supernates:

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 3/8 | Product datasheet for ABIN2472083 | 09/11/2023 | Copyright antibodies-online. All rights reserved. Remove particulates by centrifugation and assay freshly prepared samples immediately or aliquot and store samples at  $\leq$  20 °C for later use. Avoid repeated freeze thaw cycles. Serum:

Use a serum separator tube (SST) and allow samples to clot for 30 minutes before centrifugation for 10 minutes at 1000 g. Remove serum and assay freshly prepared samples immediately or aliquot and store samples at  $\leq$  20 °C for later use. Avoid repeated freeze thaw cycles. Plasma:

	Collect plasma using EDTA, citrate or heparin as anticoagulant. Centrifuge at 1000 g within 30
	minutes of collection. Assay freshly prepared samples immediately or aliquot and store
	samples at $\leq$ 20 °C for later use. Avoid repeated freeze thaw cycles. Other biological samples
	might be suitable for use in the assay. Cell culture supernates, serum and plasma were tested
	with this assay. Dilution with Assay Buffer may be needed. Note: Samples containing a visible
	precipitate must be clarified prior to use in the assay. Do not use grossly hemolyzed or lipemic
	specimens. If samples are to be run within 24 hours, they may be stored at 2 to 8 °C. For longer
	storage, aliquot samples and store frozen at 20 °C to avoid loss of bioactive human sE:
	Selectin. Avoid repeated freeze thaw cycles. Human sE:
	Selectin ELISA Kit 6 SAMPLE PREPARATION Serum and plasma samples require a 5:
	fold dilution. A suggested 5:
	fold dilution is 50 $\mu$ L sample + 200 $\mu$ L Assay Buffer (1x)
Assay Procedure:	Bring all reagents and samples to room temperature before use.
	1. Prepare all reagents including microplate, samples, standards and working solution as
	described in the previous sections.
	2. Remove excess microplate strips and return them to the foil pouch containing the desiccant
	pack, and reseal for further use.
	3. Add 300 $\mu L$ Washing Buffer (1x) per well, and allow it for about 30 seconds before aspiration.
	Soaking is highly recommended to obtain a good test performance. Empty wells and tap
	microwell strips on absorbent pad or paper towel to remove excess Washing Buffer (1x). Use
	the microwell strips immediately after washing. Do not allow wells to dry.
	4. Add 50 µL of Assay Buffer (1x) to each well.
	5. Add 50 $\mu L$ of Standard or sample to each well. Ensure reagent addition is uninterrupted and
	completed within 15 minutes.
	6. Add 50 µL of Detect Antibody to each well.
	7. Seal the plate with an adhesive film. Incubate at room temperature (18 to 25 $^\circ$ C) for 2 hours
	on a microplate shaker set at 100 rpm.

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 4/8 | Product datasheet for ABIN2472083 | 09/11/2023 | Copyright antibodies-online. All rights reserved.

	8. Aspirate each well and wash by filling each well with 300 $\mu$ L Washing Buffer (1x), repeat five
	times for a total six washes. Complete removal of liquid at each step is essential to the best
	performance. After the last wash, remove any remaining Washing Buffer (1x) by aspirating or
	decanting. Invert the plate and tap it against clean paper towels.
	9. Add 100 µL of Streptavidin-HRP to each well.
	10. Seal the plate with a fresh adhesive film. Incubate at room temperature (18 to 25 °C) for
	45 minutes on a microplate shaker set at 100 rpm. 11. Repeat aspiration/wash as in step 8.
	12. Add 100 $\mu$ L of Substrate Solution to each well. Incubate for 10 - 30 minutes at room
	temperature. Protect from light.
	13. Add 100 $\mu$ L of Stop Solution to each well. The color will turn yellow. If the color in the well is
	green or if the color change does not appear uniform, gently tap the plate to ensure thorough
	mixing.
	14. Measure the optical density value within 30 minutes by microplate reader set to 450 nm. If
	wavelength correction is available, set to 570 nm or 630 nm. If wavelength correction is not
	available, subtract readings at 570 nm or 630 nm from the readings at 450 nm. This subtractior
	will correct for optical imperfections in the plate. Reading directly at 450 nm without correction
	may generate higher concentration than true value.
Calculation of Results:	Average the duplicate optical density readings for each standards and sample, then subtract
	the average optical density value of the zero standard. Standard Concentration as horizontal
	axis, optical density (OD) Value as the vertical axis, regressing the data and create a standard
	curve using computer software. The data may be linearized by plotting the log of the sE-
	Selectin concentrations versus the log of the OD and the best fit line can be determined by
	regression analysis. This procedure will produce an adequate but less precise fit of the data.
	Note: The finally concentration of top standard is 6000 pg/mL. If instruction in this protocol
	have been followed samples have been diluted by 1:1 ratio (50 μL sample + 50 μL Assay Buffer)
	the concentration read from the standard curve must be multiplied by the dilution factor (x2). If
	serum/plasma samples have been diluted following the instruction, the final dilution factor is
	10.
Assay Precision:	Intra-assay precision (precision within an assay): Three serum-based and buffer-based samples
	of known concentration were tested twenty times on one plate to assess intra-assay precision.
	Inter-assay precision (precision between assays): Three serum-based and buffer-based
	samples of known concentration were tested in six separate assays to assess inter-assay
	precision.
	Spike recovery:

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 5/8 | Product datasheet for ABIN2472083 | 09/11/2023 | Copyright antibodies-online. All rights reserved.

Spike recovery:

The spike recovery was evaluated by spiking 3 levels of human sE-Selectin into five health human serum samples. The un-spiked serum was used as blank in these experiments. The recovery ranged from 86 % to 111 % with an overall mean recovery of 99 %.
For Research Use only
Buffer contains: 0.02 % sodium azide
Sodium azide
This product contains Sodium azide: a POISONOUS AND HAZARDOUS SUBSTANCE which should be handled by trained staff only.
<ul> <li>Intended for research use only and are not for use in diagnostic or therapeutic procedures.</li> <li>Treat all chemicals with caution because they can be potentially hazardous.</li> <li>It is recommended that this product is handled only by persons who have been trained in</li> <li>laboratory techniques and in accordance with the principles of good laboratory practice. Wear</li> <li>personal protection equipment such as laboratory coat, safety glasses and gloves.</li> <li>Avoid direct contact with skin or eyes. Wash immediately with water in the case of contact with skin or eyes. Avoid contact of skin or mucous membranes with kit reagents or specimens. See material safety data sheet(s) for specific advice.</li> <li>Pure water or deionized water must be used for reagent preparation.</li> <li>The Stop Solution provided with this kit is an acid solution. Wear personal protection equipment with caution.</li> </ul>

Rubber or disposable latex gloves should be worn while handling kit reagents or specimens.

Avoid contact of substrate solution with oxidizing agents and metal.

Avoid splashing or generation of aerosols.

Use disposable pipette tips and/or pipettes to avoid microbial or cross-contamination of reagents or specimens which may invalidate the test.

Use clean, dedicated reagent trays for dispensing the conjugate and substrate reagent.

Exposure to acid inactivates the HRP and antibody conjugate.

Substrate solution must be warmed to room temperature prior to use.


Decontaminate and dispose specimens and all potentially contaminated materials as they

could contain infectious agents. The preferred method of decontamination is autoclaving for a minimum of 1 hour at 121.5 °C.

Liquid wastes not containing acid and neutralized waste may be mixed with sodium

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 6/8 | Product datasheet for ABIN2472083 | 09/11/2023 | Copyright antibodies-online. All rights reserved.

#### Images



### ELISA

Image 1.

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 8/8 | Product datasheet for ABIN2472083 | 09/11/2023 | Copyright antibodies-online. All rights reserved.