

Datasheet for ABIN2533176

EXT1 ELISA Kit



Overview

Overview	
Quantity:	96 tests
Target:	EXT1
Reactivity:	Human
Method Type:	Sandwich ELISA
Detection Range:	0.781 ng/mL - 50 ng/mL
Minimum Detection Limit:	0.781 ng/mL
Application:	ELISA
Product Details	
Purpose:	Human Exostoses 1 ELISA Kit is a sandwich ELISA kit for use with Tissue homogenates and
	other biological fluids. This assay has high sensitivity and excellent specificity for detection of
	Exostoses 1 (EXT1)
	No significant cross-reactivity or interference between Exostoses 1 (EXT1) and analogues was
	observed.
Sample Type:	Tissue Homogenate
Analytical Method:	Quantitative
Detection Method:	Colorimetric
Specificity:	This assay has high sensitivity and excellent specificity for detection of Exostoses 1 (EXT1)
Sensitivity:	< 0.31 ng/mL

Target Details

Expiry Date:

6 months

Target:	EXT1
Alternative Name:	Exostoses 1 (EXT1) (EXT1 Products)
Pathways:	Glycosaminoglycan Metabolic Process
Application Details	
Application Notes:	Stability: The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5 % within the expiration date under appropriate storage conditions. To minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout. Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user. Standard Form: Lyophilized
Comment:	The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5% within the expiration date under appropriate storage conditions minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout.
Plate:	Pre-coated
Restrictions:	For Research Use only
Handling	
Storage:	4 °C/-20 °C
Storage Comment:	Upon receipt, store the kit according to the storage instruction in the kit's manual.