antibodies -online.com





Datasheet for ABIN2534217

Haptoglobin beta Chain (HP) ELISA Kit



Overview

Overview	
Quantity:	96 tests
Target:	Haptoglobin beta Chain (HP)
Reactivity:	Mouse
Method Type:	Sandwich ELISA
Detection Range:	31.25 ng/mL - 2000 ng/mL
Minimum Detection Limit:	31.25 ng/mL
Application:	ELISA
Product Details	
Purpose:	Mouse Haptoglobin ELISA Kit is a sandwich ELISA kit for use with Serum, plasma and other
	biological fluids. This assay has high sensitivity and excellent specificity for detection of
	Haptoglobin (Hpt)
	No significant cross-reactivity or interference between Haptoglobin (Hpt) and analogues was
	observed.
Sample Type:	Plasma, Serum
Analytical Method:	Quantitative
Detection Method:	Colorimetric
Specificity:	This assay has high sensitivity and excellent specificity for detection of Haptoglobin (Hpt)
Sensitivity:	< 13.1 ng/mL

Target Details

Target:	Haptoglobin beta Chain (HP)
Alternative Name:	Haptoglobin beta (HP Products)
Application Details	
Application Notes:	Stability: The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5 % within the expiration date under appropriate storage conditions. To minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout. Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user. Standard Form: Lyophilized
Comment:	The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5% within the expiration date under appropriate storage conditions minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout.
Plate:	Pre-coated
Restrictions:	For Research Use only
Handling	
Storage:	4 °C/-20 °C
Storage Comment:	Upon receipt, store the kit according to the storage instruction in the kit's manual.
Expiry Date:	6 months