

Datasheet for ABIN2534842

ITK ELISA Kit



Overview

Quantity:	96 tests
Target:	ITK
Reactivity:	Mouse
Method Type:	Sandwich ELISA
Detection Range:	0.312 ng/mL - 20 ng/mL
Minimum Detection Limit:	0.312 ng/mL
Application:	ELISA
Product Details	
Purpose:	Mouse ITK ELISA Kit is a sandwich ELISA kit for use with Tissue homogenates, cell lysates and
	other biological fluids. This assay has high sensitivity and excellent specificity for detection of
	IL2 Inducible T-Cell Kinase (ITK)
	No significant cross-reactivity or interference between IL2 Inducible T-Cell Kinase (ITK) and
	analogues was observed.
Sample Type:	Cell Lysate, Tissue Homogenate
Analytical Method:	Quantitative
Detection Method:	Colorimetric
Specificity:	This assay has high sensitivity and excellent specificity for detection of IL2 Inducible T-Cell
	Kinase (ITK)
Sensitivity:	< 0.119 ng/mL

Target Details

l arget Details	
Target:	ITK
Alternative Name:	IL2 Inducible T-Cell Kinase (ITK) (ITK Products)
Pathways:	TCR Signaling, Fc-epsilon Receptor Signaling Pathway
Application Details	
Application Notes:	Stability: The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5 % within the expiration date under appropriate storage conditions. To minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout. Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user. Standard Form: Lyophilized
Comment:	The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5% within the expiration date under appropriate storage conditions minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout.
Plate:	Pre-coated
Restrictions:	For Research Use only
Handling	
Storage:	4 °C/-20 °C
Storage Comment:	Upon receipt, store the kit according to the storage instruction in the kit's manual.
Expiry Date:	6 months