

## Datasheet for ABIN2535733

## **JAM2 ELISA Kit**



$\sim$				
( )	ve	r\/		Λ/
$\cup$	$\vee \subset$	1 V I	$\Box$	٧V

Quantity:	96 tests	
Target:	JAM2	
Reactivity:	Human	
Method Type:	Sandwich ELISA	
Detection Range:	31.25 pg/mL - 2000 pg/mL	
Minimum Detection Limit:	31.25 pg/mL	
Application:	ELISA	
Product Details		
Purpose:	Human JAM2 ELISA Kit is a sandwich ELISA kit for use with Serum, plasma, tissue	
	homogenates and other biological fluids. This assay has high sensitivity and excellent	
	specificity for detection of Junctional Adhesion Molecule 2 (JAM2)	
	No significant cross-reactivity or interference between Junctional Adhesion Molecule 2 (JAM2)	
	and analogues was observed.	
Sample Type:	Plasma, Serum, Tissue Homogenate	
Analytical Method:	Quantitative	
Detection Method:	Colorimetric	
Specificity:	This assay has high sensitivity and excellent specificity for detection of Junctional Adhesion	
	Molecule 2 (JAM2)	
Sensitivity:	< 12.7 pg/mL	

## Target Details

Target:	JAM2	
Abstract:	JAM2 Products	
Application Details		
Application Notes:	Stability: The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5 % within the expiration date under appropriate storage conditions. To minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout. Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user.  Standard Form: Lyophilized	
Comment:	The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5% within the expiration date under appropriate storage conditions minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout.	
Plate:	Pre-coated	
Restrictions:	For Research Use only	
Handling		
Storage:	4 °C/-20 °C	
Storage Comment:	Upon receipt, store the kit according to the storage instruction in the kit's manual.	
Expiry Date:	6 months	