

Datasheet for ABIN2535971

LMAN2 ELISA Kit



Overview

Quantity:	96 tests
Target:	LMAN2
Reactivity:	Human
Method Type:	Sandwich ELISA
Detection Range:	0.156 ng/mL - 10 ng/mL
Minimum Detection Limit:	0.156 ng/mL
Application:	ELISA
Product Details	
Purpose:	Human LMAN2 ELISA Kit is a sandwich ELISA kit for use with Serum, plasma, tissue
	homogenates and other biological fluids. This assay has high sensitivity and excellent
	specificity for detection of Lectin, Mannose Binding 2 (LMAN2)
	No significant cross-reactivity or interference between Lectin, Mannose Binding 2 (LMAN2) and
	analogues was observed.
Sample Type:	Plasma, Serum, Tissue Homogenate
Analytical Method:	Quantitative
Detection Method:	Colorimetric
Specificity:	This assay has high sensitivity and excellent specificity for detection of Lectin, Mannose
	Binding 2 (LMAN2)
Sensitivity:	< 0.055 ng/mL

Target Details

Expiry Date:

6 months

Target Details	
Target:	LMAN2
Alternative Name:	Lectin, Mannose Binding 2 (LMAN2) (LMAN2 Products)
Pathways:	SARS-CoV-2 Protein Interactome
Application Details	
Application Notes:	Stability: The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5 % within the expiration date under appropriate storage conditions. To minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout. Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user. Standard Form: Lyophilized
Comment:	The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5% within the expiration date under appropriate storage conditions minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout.
Plate:	Pre-coated
Restrictions:	For Research Use only
Handling	
Storage:	4 °C/-20 °C
Storage Comment:	Upon receipt, store the kit according to the storage instruction in the kit's manual.