

Datasheet for ABIN2536683

MFGE8 ELISA Kit



C	۱۱ /	\cap	~\ /	ic	11/	1
	V	CI	V	IF	٧,	۷

Quantity:	96 tests
Target:	MFGE8
Reactivity:	Mouse
Method Type:	Sandwich ELISA
Detection Range:	15.62 pg/mL - 1000 pg/mL
Minimum Detection Limit:	15.62 pg/mL
Application:	ELISA
Product Details	
Purpose:	Mouse MFGE8 ELISA Kit is a sandwich ELISA kit for use with Serum, plasma and other
	biological fluids. This assay has high sensitivity and excellent specificity for detection of Milk
	Fat Globule EGF Factor 8 (MFGE8)
	No significant cross-reactivity or interference between Milk Fat Globule EGF Factor 8 (MFGE8)
	and analogues was observed.
Sample Type:	Plasma, Serum
Analytical Method:	Quantitative
Detection Method:	Colorimetric
Specificity:	This assay has high sensitivity and excellent specificity for detection of Milk Fat Globule EGF
	Factor 8 (MFGE8)
Sensitivity:	< 6.0 pg/mL

Target Details

. a. got z otalio	
Target:	MFGE8
Alternative Name:	Milk Fat Globule EGF Factor 8 (MFGE8) (MFGE8 Products)
Pathways:	SARS-CoV-2 Protein Interactome
Application Details	
Application Notes:	Stability: The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5 % within the expiration date under appropriate storage conditions. To minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout. Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user. Standard Form: Lyophilized
Comment:	The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5% within the expiration date under appropriate storage conditions minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout.
Plate:	Pre-coated
Restrictions:	For Research Use only
Handling	
Storage:	4 °C/-20 °C
Storage Comment:	Upon receipt, store the kit according to the storage instruction in the kit's manual.
Expiry Date:	6 months