

Datasheet for ABIN2537733

PMP22 ELISA Kit



C	۱۱ /	\cap	~\ /	ic	11/	1
	V	CI	V	IF	٧,	٧

Quantity:	96 tests
Target:	PMP22
Reactivity:	Mouse
Method Type:	Sandwich ELISA
Detection Range:	0.156 ng/mL - 10 ng/mL
Minimum Detection Limit:	0.156 ng/mL
Application:	ELISA
Product Details	
Purpose:	Mouse PMP22 ELISA Kit is a sandwich ELISA kit for use with Tissue homogenates, cell lysates and other biological fluids. This assay has high sensitivity and excellent specificity for detection
	of Peripheral Myelin Protein 22 (PMP22)
	No significant cross-reactivity or interference between Peripheral Myelin Protein 22 (PMP22)
	and analogues was observed.
Sample Type:	Cell Lysate, Tissue Homogenate
Analytical Method:	Quantitative
Detection Method:	Colorimetric
Specificity:	This assay has high sensitivity and excellent specificity for detection of Peripheral Myelin Protein 22 (PMP22)
Sensitivity:	< 0.058 ng/mL

Target Details

Target:	PMP22	
Abstract:	PMP22 Products	
Application Details		
Application Notes:	Stability: The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5 % within the expiration date under appropriate storage conditions. To minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout. Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user. Standard Form: Lyophilized	
Comment:	The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5% within the expiration date under appropriate storage conditions minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout.	
Plate:	Pre-coated	
Restrictions:	For Research Use only	
Handling		
Storage:	4 °C/-20 °C	
Storage Comment:	Upon receipt, store the kit according to the storage instruction in the kit's manual.	
Expiry Date:	6 months	