

## Datasheet for ABIN2537768 **PES1 ELISA Kit**



[Go to Product page](#)

### Overview

Quantity:	96 tests
Target:	PES1
Reactivity:	Human
Method Type:	Sandwich ELISA
Detection Range:	0.156 ng/mL - 10 ng/mL
Minimum Detection Limit:	0.156 ng/mL
Application:	ELISA

### Product Details

Purpose:	Human Pescadillo Homolog 1 ELISA Kit is a sandwich ELISA kit for use with Tissue homogenates, cell lysates and other biological fluids. This assay has high sensitivity and excellent specificity for detection of Pescadillo Homolog 1 (PES1) No significant cross-reactivity or interference between Pescadillo Homolog 1 (PES1) and analogues was observed.
Sample Type:	Cell Lysate, Tissue Homogenate
Analytical Method:	Quantitative
Detection Method:	Colorimetric
Specificity:	This assay has high sensitivity and excellent specificity for detection of Pescadillo Homolog 1 (PES1)
Sensitivity:	< 0.063 ng/mL

## Target Details

---

Target: PES1

---

Abstract: [PES1 Products](#)

## Application Details

---

Application Notes: Stability: The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5 % within the expiration date under appropriate storage conditions. To minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout. Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user.  
Standard Form: Lyophilized

---

Comment: The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5% within the expiration date under appropriate storage conditions. - minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout.

---

Plate: Pre-coated

---

Restrictions: For Research Use only

---

## Handling

---

Storage: 4 °C/-20 °C

---

Storage Comment: Upon receipt, store the kit according to the storage instruction in the kit's manual.

---

Expiry Date: 6 months