

Datasheet for ABIN2538060 **WBP2NL ELISA Kit**



[Go to Product page](#)

Overview

Quantity: 96 tests

Target: WBP2NL

Reactivity: Cow

Method Type: Sandwich ELISA

Detection Range: 0.156 ng/mL - 10 ng/mL

Minimum Detection Limit: 0.156 ng/mL

Application: ELISA

Product Details

Purpose: Cow PAWP ELISA Kit is a sandwich ELISA kit for use with Tissue homogenates and other biological fluids. This assay has high sensitivity and excellent specificity for detection of Postacrosomal Sheath WW Domain Binding Protein (PAWP)
No significant cross-reactivity or interference between Postacrosomal Sheath WW Domain Binding Protein (PAWP) and analogues was observed.

Sample Type: Tissue Homogenate

Analytical Method: Quantitative

Detection Method: Colorimetric

Specificity: This assay has high sensitivity and excellent specificity for detection of Postacrosomal Sheath WW Domain Binding Protein (PAWP)

Sensitivity: < 0.059 ng/mL

Target Details

Target:	WBP2NL
Alternative Name:	Postacrosomal Sheath WW Domain Binding Protein (PAWP) (WBP2NL Products)

Application Details

Application Notes:	<p>Stability: The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5 % within the expiration date under appropriate storage conditions. To minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout.</p> <p>Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user.</p> <p>Standard Form: Lyophilized</p>
Comment:	<p>The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5% within the expiration date under appropriate storage conditions. - minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout.</p>
Plate:	Pre-coated
Restrictions:	For Research Use only

Handling

Storage:	4 °C/-20 °C
Storage Comment:	Upon receipt, store the kit according to the storage instruction in the kit's manual.
Expiry Date:	6 months