

Datasheet for ABIN2539819

TALDO1 ELISA Kit



Go to Product page

_	verview				
	1//	r	1//	\triangle	۸/
	V		VI		/ V

Quantity:	96 tests	
Target:	TALD01	
Reactivity:	Human	
Method Type:	Sandwich ELISA	
Detection Range:	0.156 ng/mL - 10 ng/mL	
Minimum Detection Limit:	0.156 ng/mL	
Application:	ELISA	
Product Details		
Purpose:	Human Transaldolase ELISA Kit is a sandwich ELISA kit for use with Serum, plasma, tissue	
	homogenates and other biological fluids. This assay has high sensitivity and excellent	
	specificity for detection of Transaldolase (TALDO1)	
	No significant cross-reactivity or interference between Transaldolase (TALDO1) and analogues	
	was observed.	
Sample Type:	Plasma, Serum, Tissue Homogenate	
Analytical Method:	Quantitative	
Detection Method:	Colorimetric	
Specificity:	This assay has high sensitivity and excellent specificity for detection of Transaldolase	
	(TALDO1)	
Sensitivity:	< 0.058 ng/mL	

Target Details

Target:	TALDO1	
Alternative Name:	Transaldolase (TALDO1) (TALDO1 Products)	
Application Details		
Application Notes:	Stability: The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5 % within the expiration date under appropriate storage conditions. To minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout. Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user. Standard Form: Lyophilized	
Comment:	The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5% within the expiration date under appropriate storage conditions minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout.	
Plate:	Pre-coated	
Restrictions:	For Research Use only	
Handling		
Storage:	4 °C/-20 °C	
Storage Comment:	Upon receipt, store the kit according to the storage instruction in the kit's manual.	
Expiry Date:	6 months	