

Datasheet for ABIN2713727

CD47 Protein (CD47) (Transcript Variant 2) (Fc Tag)



Overview

Quantity:	10 μg
Target:	CD47
Protein Characteristics:	Transcript Variant 2
Origin:	Human
Source:	HEK-293 Cells
Protein Type:	Recombinant
Purification tag / Conjugate:	This CD47 protein is labelled with Fc Tag.
Application:	Antibody Production (AbP), Standard (STD)
Product Details	
Characteristics:	 Recombinant human CD47 (transcript variant 2) protein expressed in HEK293 cells. Produced with end-sequenced ORF clone
Purity:	> 95 % as determined by SDS-PAGE and Coomassie blue staining
Endotoxin Level:	Endotoxin level is <0.1 ng/μg of protein (<1EU/μg).
Target Details	
Target:	CD47
Alternative Name:	CD47 (CD47 Products)
Background:	This gene encodes a membrane protein, which is involved in the increase in intracellular calcium concentration that occurs upon cell adhesion to extracellular matrix. The encoded

Target Details

- Target Details	
	protein is also a receptor for the C-terminal cell binding domain of thrombospondin, and it may
	play a role in membrane transport and signal transduction. This gene has broad tissue
	distribution, and is reduced in expression on Rh erythrocytes. Alternatively spliced transcript
	variants have been found for this gene.
Molecular Weight:	40.8kD
NCBI Accession:	NP_942088
Application Details	
Application Notes:	Recombinant human proteins can be used for:
	Native antigens for optimized antibody production
	Positive controls in ELISA and other antibody assays
Comment:	The tag is located at the C-terminal.
Restrictions:	For Research Use only
Handling	
Buffer:	Lyophilized from a 0.2 µM filtered solution of PBS, pH 7.4
Handling Advice:	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not
	recommended to reconstitute to a concentration less than 100 µg/mL. Dissolve the lyophilized
	protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Storage:	-80 °C
Storage Comment:	Store at -80°C. Thaw on ice, aliquot to individual single-use tubes, and then re-freeze

immediately. Only 2-3 freeze thaw cycles are recommended.