# antibodies .- online.com





## MARCH10 Protein (His tag)





( )	11/	IN	/ie	A .
	/ // <del> </del>	۱ ات	/   (−	' \/\/

Quantity:	50 µg
Target:	MARCH10
Origin:	Human
Source:	Escherichia coli (E. coli)
Protein Type:	Recombinant
Purification tag / Conjugate:	This MARCH10 protein is labelled with His tag.
Application:	Antibody Production (AbP), Standard (STD)
Product Details	
Characteristics:	<ul> <li>Recombinant human MARCH10 / RNF190 (full length, N-term HIS tag, transcript variant 1) protein expressed in E.coli.</li> </ul>
Characteristics:	<ul> <li>Recombinant human MARCH10 / RNF190 (full length, N-term HIS tag, transcript variant 1) protein expressed in E.coli.</li> <li>Produced with end-sequenced ORF clone</li> </ul>
Characteristics:  Purity:	protein expressed in E.coli.
	protein expressed in E.coli.  • Produced with end-sequenced ORF clone
Purity:	protein expressed in E.coli.  • Produced with end-sequenced ORF clone
Purity: Target Details	<ul> <li>protein expressed in E.coli.</li> <li>Produced with end-sequenced ORF clone</li> <li>&gt; 80 % as determined by SDS-PAGE and Coomassie blue staining</li> </ul>
Purity: Target Details Target:	protein expressed in E.coli.  • Produced with end-sequenced ORF clone  > 80 % as determined by SDS-PAGE and Coomassie blue staining  MARCH10

proteins, thereby signaling their vesicular transport between membrane compartments

(Morokuma et al., 2007 [PubMed 17604280]).[supplied by OMIM, Apr 2010]

### **Target Details**

Molecular Weight:	90.3 kDa
NCBI Accession:	NP_689811


## **Application Details**

Application Notes:	Recombinant human proteins can be used for:
	Native antigens for optimized antibody production
	Positive controls in ELISA and other antibody assays
Comment:	The tag is located at the N-terminal.
Restrictions:	For Research Use only

### Handling

Concentration:	50 μg/mL
Buffer:	25 mM Tris, pH 8.0, 150 mM NaCl, 10 % glycerol, 1 % Sarkosyl.
Storage:	-80 °C
Storage Comment:	Store at -80°C. Thaw on ice, aliquot to individual single-use tubes, and then re-freeze immediately. Only 2-3 freeze thaw cycles are recommended.

#### **Images**



#### **Western Blotting**

Image 1. Validation with Western Blot