antibodies

Datasheet for ABIN2732605 SPIRE2 Protein (Myc-DYKDDDDK Tag)

Image




Overview

| Quantity: | 20 µg |
|---|--|
| Target: | SPIRE2 |
| Origin: | Human |
| Source: | HEK-293 Cells |
| Protein Type: | Recombinant |
| Purification tag / Conjugate: | This SPIRE2 protein is labelled with Myc-DYKDDDDK Tag. |
| Application: | Antibody Production (AbP), Standard (STD) |
| | |
| Product Details | |
| | |
| Characteristics: | Recombinant human SPIRE2 protein expressed in HEK293 cells. |
| Characteristics: | |
| Characteristics: | |
| Characteristics: Purity: | |
| | Produced with end-sequenced ORF clone |
| | Produced with end-sequenced ORF clone |
| Purity: | Produced with end-sequenced ORF clone |
| Purity: Target Details | Produced with end-sequenced ORF clone > 80 % as determined by SDS-PAGE and Coomassie blue staining |
| Purity: Target Details Target: | Produced with end-sequenced ORF clone > 80 % as determined by SDS-PAGE and Coomassie blue staining SPIRE2 |
| Purity: Target Details Target: Alternative Name: | Produced with end-sequenced ORF clone > 80 % as determined by SDS-PAGE and Coomassie blue staining SPIRE2 Spire2 (SPIRE2 Products) |

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/2 | Product datasheet for ABIN2732605 | 09/11/2023 | Copyright antibodies-online. All rights reserved.

| Application Details | |
|---------------------|--|
| Application Notes: | Recombinant human proteins can be used for: Native antigens for optimized antibody production Positive controls in ELISA and other antibody assays |
| Comment: | The tag is located at the C-terminal. |
| Restrictions: | For Research Use only |
| Handling | |
| Concentration: | 50 μg/mL |
| Buffer: | 25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10 % glycerol. |
| Storage: | -80 °C |
| Storage Comment: | Store at -80°C. Thaw on ice, aliquot to individual single-use tubes, and then re-freeze immediately. Only 2-3 freeze thaw cycles are recommended. |

Images



Western Blotting

Image 1. Validation with Western Blot