

Datasheet for ABIN2811835

anti-Complement C4 antibody (AA 21-120) (AbBy Fluor® 594)



Go to Product page

This Complement C4 antibody is conjugated to AbBy Fluor® 594	
mbedded	
- Indicate	

	C4B3, C4B5, C4BD, C4B12, C4B_2, CPAMD3, Complement C4-A, Acidic complement C4, C3 a		
	PZP-like alpha-2-macroglobulin domain-containing protein 2, C4A		
	Background: Non-enzymatic component of C3 and C5 convertases and thus essential for the		
	propagation of the classical complement pathway. Covalently binds to immunoglobulins and		
	immune complexes and enhances the solubilization of immune aggregates and the clearance		
	of IC through CR1 on erythrocytes. C4A isotype is responsible for effective binding to form		
	amide bonds with immune aggregates or protein antigens, while C4B isotype catalyzes the		
	transacylation of the thioester carbonyl group to form ester bonds with carbohydrate antigens.		
	Derived from proteolytic degradation of complement C4, C4a anaphylatoxin is a mediator of		
	local inflammatory process. It induces the contraction of smooth muscle, increases vascular		
	permeability and causes histamine release from mast cells and basophilic leukocytes.		
Gene ID:	720, 721		
UniProt:	P0C0L4		
Pathways:	Complement System		
Application Details			
Application Notes:	IF(IHC-P) 1:50-200		
	IF(IHC-F) 1:50-200		
	IF(ICC) 1:50-200		
Restrictions:	For Research Use only		
Handling			
Format:	Liquid		
Concentration:	1 μg/μL		
Buffer:	Aqueous buffered solution containing 0.01M TBS (pH 7.4) with 1 % BSA, 0.03 % Proclin300 and		
	50 % Glycerol.		
Preservative:	ProClin		
Precaution of Use:	This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be		
	handled by trained staff only.		
Storage:	-20 °C		
Storage Comment:	Store at -20°C. Aliquot into multiple vials to avoid repeated freeze-thaw cycles.		

	1	l:
-	ココロロ	urnen
	Hand	III IU

Expiry Date:

12 months