antibodies .- online.com







anti-CYP46A1 antibody (AA 51-150) (AbBy Fluor® 594)



()	1/0	r\ /1	014	
()	ve	I V I	-v	V

Quantity:	100 μL	
Target:	CYP46A1	
Binding Specificity:	AA 51-150	
Reactivity:	Mouse	
Host:	Rabbit	
Clonality:	Polyclonal	
Conjugate:	This CYP46A1 antibody is conjugated to AbBy Fluor® 594	
Application:	Western Blotting (WB), Immunofluorescence (Cultured Cells) (IF (cc)), Immunofluorescence (Paraffin-embedded Sections) (IF (p))	

Product Details

Immunogen:	KLH conjugated synthetic peptide derived from human CYP46
Isotype:	IgG
Cross-Reactivity:	Mouse
Predicted Reactivity:	Human,Rat,Dog,Cow,Sheep,Horse,Rabbit
Purification:	Purified by Protein A.

Target Details

Target:	CYP46A1
Alternative Name:	CYP46 (CYP46A1 Products)

Target Details

Background:	Synonyms: Cholesterol 24-hydroxylase, CP46, CYP46, CP46A_HUMAN, cytochrome P450 family	
	46.	
	Background: The Cyp46 gene encodes a member of the cytochrome P450 superfamily of	
	enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions	
	involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This	
	endoplasmic reticulum protein is expressed exclusively in the brain, where it converts	
	cholesterol to 24S-hydroxycholesterol by adding a hydroxyl group to cholesterol, producing a	
	product that is more soluble than cholesterol and able to be exported from the brain.Cyp46 is	
	also known as 24S-cholesterol hydroxylase.	
Gene ID:	10858	
Application Details		
Application Notes:	IF(IHC-P) 1:50-200	
	IF(IHC-F) 1:50-200	
	IF(ICC) 1:50-200	
Restrictions:	For Research Use only	
Handling		
Format:	Liquid	
Concentration:	1 μg/μL	
Buffer:	Aqueous buffered solution containing 0.01M TBS (pH 7.4) with 1 % BSA, 0.03 % Proclin300 and	
	50 % Glycerol.	
Preservative:	ProClin	
Precaution of Use:	This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be	
	handled by trained staff only.	
Storage:	-20 °C	
Storage Comment:	Store at -20°C. Aliquot into multiple vials to avoid repeated freeze-thaw cycles.	
Expiry Date:	12 months	