

## Datasheet for ABIN2945536

### **BIRC2 ELISA Kit**



[Go to Product page](#)

#### Overview

|                          |                       |
|--------------------------|-----------------------|
| Quantity:                | 96 tests              |
| Target:                  | BIRC2                 |
| Reactivity:              | Human                 |
| Method Type:             | Sandwich ELISA        |
| Detection Range:         | 78 pg/mL - 5000 pg/mL |
| Minimum Detection Limit: | 78 pg/mL              |
| Application:             | ELISA                 |

#### Product Details

|                    |                                                                                                                                                |
|--------------------|------------------------------------------------------------------------------------------------------------------------------------------------|
| Purpose:           | Human Baculoviral IAP repeat-containing protein 2 ELISA Kit is an ELISA kit against Human Baculoviral IAP repeat-containing protein 2 (BIRC2). |
| Sample Type:       | Plasma, Serum                                                                                                                                  |
| Analytical Method: | Quantitative                                                                                                                                   |
| Detection Method:  | Colorimetric                                                                                                                                   |

#### Target Details

|                   |                                                                                                                                                                                |
|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Target:           | BIRC2                                                                                                                                                                          |
| Alternative Name: | Baculoviral IAP Repeat-Containing Protein 2 ( <a href="#">BIRC2 Products</a> )                                                                                                 |
| Pathways:         | <a href="#">Apoptosis</a> , <a href="#">Caspase Cascade in Apoptosis</a> , <a href="#">Activation of Innate immune Response</a> , <a href="#">Toll-Like Receptors Cascades</a> |

## Application Details

---

|                    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|--------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Application Notes: | Stability: The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5 % within the expiration date under appropriate storage conditions. To minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout. Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user. |
|--------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

---

|        |            |
|--------|------------|
| Plate: | Pre-coated |
|--------|------------|

---

|               |                       |
|---------------|-----------------------|
| Restrictions: | For Research Use only |
|---------------|-----------------------|

---

## Handling

---

|          |             |
|----------|-------------|
| Storage: | 4 °C/-20 °C |
|----------|-------------|

---

|                  |                                                                                       |
|------------------|---------------------------------------------------------------------------------------|
| Storage Comment: | Upon receipt, store the kit according to the storage instruction in the kit's manual. |
|------------------|---------------------------------------------------------------------------------------|

---

|              |          |
|--------------|----------|
| Expiry Date: | 6 months |
|--------------|----------|

---