

# Datasheet for ABIN2948397

# **ICOSLG ELISA Kit**



Go to Product page

|  | er |  |  |
|--|----|--|--|
|  |    |  |  |
|  |    |  |  |

| Quantity:                | 96 tests               |
|--------------------------|------------------------|
| Target:                  | ICOSLG                 |
| Reactivity:              | Mouse                  |
| Method Type:             | Sandwich ELISA         |
| Detection Range:         | 0.156 ng/mL - 10 ng/mL |
| Minimum Detection Limit: | 0.156 ng/mL            |
| Application:             | ELISA                  |

## **Product Details**

| Purpose:           | Mouse ICOS ligand ELISA Kit is an ELISA kit against Mouse ICOS ligand (Icoslg). |
|--------------------|---|
| Sample Type:       | Plasma, Serum   |
| Analytical Method: | Quantitative  |
| Detection Method:  | Colorimetric  |

# **Target Details**

| Target:  | ICOSLG                        |
|--|-------------------------------|
| Alternative Name:  | ICOS Ligand (ICOSLG Products) |
| Pathways: Production of Molecular Mediator of Immune Response Activated T Cell Proliferation |                               |

# **Application Details**

| Application Notes: | Stability: The stability of the kit is determined by the rate of activity loss. The loss rate is less |
|--------------------|---|
|                    | than 5 % within the expiration date under appropriate storage conditions. To minimize                 |
|                    | performance fluctuations, operation procedures and lab conditions should be strictly controlled.      |
|                    | It is also strongly suggested that the whole assay is performed by the same user throughout.          |
|                    | Recommended dilutions: Optimal dilutions/concentrations should be determined by the end               |
|                    | user.   |
| Plate:             | Pre-coated  |
| Restrictions:      | For Research Use only   |
| Handling           |   |
| Storage:           | 4 °C/-20 °C   |
| Storage Comment:   | Upon receipt, store the kit according to the storage instruction in the kit's manual.                 |
| Expiry Date:       | 6 months  |