

## Datasheet for ABIN2952238 **TYK2 ELISA Kit**



[Go to Product page](#)

### Overview

|                          |                        |
|--------------------------|------------------------|
| Quantity:                | 96 tests               |
| Target:                  | TYK2                   |
| Reactivity:              | Human                  |
| Method Type:             | Sandwich ELISA         |
| Detection Range:         | 0.313 ng/mL - 20 ng/mL |
| Minimum Detection Limit: | 0.313 ng/mL            |
| Application:             | ELISA                  |

### Product Details

|                    |  |
|--------------------|--|
| Purpose:           | Human Tyrosine Kinase 2 ELISA Kit is an ELISA Kit against Tyrosine Kinase 2. |
| Sample Type:       | Plasma, Serum  |
| Analytical Method: | Quantitative   |
| Detection Method:  | Colorimetric   |
| Sensitivity:       | 0.19 ng/mL   |

### Target Details

|                   |  |
|-------------------|--|
| Target:           | TYK2   |
| Alternative Name: | Tyrosine Kinase 2 ( <a href="#">TYK2 Products</a> )              |
| Pathways:         | <a href="#">JAK-STAT Signaling</a> , <a href="#">Hepatitis C</a> |

## Application Details

---

**Application Notes:**                      **Stability:** The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5 % within the expiration date under appropriate storage conditions. To minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout. Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user.

Standard Form: Lyophilized

---

**Plate:**                                      Pre-coated

---

**Restrictions:**                          For Research Use only

---

## Handling

---

**Storage:**                                  4 °C/-20 °C

---

**Storage Comment:**                      Upon receipt, store the kit according to the storage instruction in the kit's manual.

---

**Expiry Date:**                              6 months

---