# antibodies -online.com





## **HGSNAT Protein (AA 1-190) (His tag)**



#### Overview

Quantity:	1 mg
Target:	HGSNAT
Protein Characteristics:	AA 1-190
Origin:	Human
Source:	Insect Cells
Protein Type:	Recombinant
Purification tag / Conjugate:	This HGSNAT protein is labelled with His tag.
Application:	ELISA, Western Blotting (WB), SDS-PAGE (SDS), Crystallization (Crys)
Product Details	

Application:	ELISA, Western Blotting (WB), SDS-PAGE (SDS), Crystallization (Crys)
Product Details	
Sequence:	MTGARASAAE QRRAGRSGQA RAAERAAGMS GAGRALAALL LAASVLSAAL LAPGGSSGRD
	AQAAPPRDLD KKRHAELKMD QALLLIHNEL LWTNLTVYWK SECCYHCLFQ VLVNVPQSPK
	AGKPSAAAAS VSTQHGSILQ LNDTLEEKEV CRLEYRFGEF GNYSLLVKNI HNGVSEIACD
	LAVNEDPVDS
	Sequence without tag. Tag location is at the discretion of the manufacturer. If you have a
	special request, please contact us.
Characteristics:	<ul> <li>Made in Germany - from design to production - by highly experienced protein experts.</li> <li>Human HGSNAT Protein (raised in Insect Cells) purified by multi-step, protein-specific process to ensure crystallization grade.</li> <li>State-of-the-art algorithm used for plasmid design (Gene synthesis).</li> </ul>
	This protein is a made to order protein and will be made for the first time for your order. Our

experts in the lab will ensure that you receive a correctly folded protein.

The big advantage of ordering our made-to-order proteins in comparison to ordering custom made proteins from other companies is that there is no financial obligation in case the protein cannot be expressed or purified.

In the unlikely event that the protein cannot be expressed or purified we do not charge anything (other companies might charge you for any performed steps in the expression process for custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression experiments or purification optimization).

When you order this made-to-order protein you will only pay upon receival of the correctly folded protein. With no financial risk on your end you can rest assured that our experienced protein experts will do everything to make sure that you receive the protein you ordered.

The concentration of our recombinant proteins is measured using the absorbance at 280nm.

The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.

The concentration of the protein is calculated using its specific absorption coefficient. We use the Expasy's protparam tool to determine the absorption coefficient of each protein.

#### Purification:

Two step purification of proteins expressed in baculovirus infected SF9 insect cells:

- 1. In a first purification step, the protein is purified from the cleared cell lysate using three different His-tag capture materials: high yield, EDTA resistant, or DTT resistant. Eluate fractions are analyzed by SDS-PAGE.
- 2. Protein containing fractions of the best purification are subjected to second purification step through size exclusion chromatography. Eluate fractions are analyzed by SDS-PAGE and Western blot.

Purity:

>95 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.

Sterility:

0.22 µm filtered

Endotoxin Level:

Protein is endotoxin free.

Grade:

Crystallography grade

#### **Target Details**

Target:	HGSNAT
Alternative Name:	HGSNAT (HGSNAT Products)
Background:	Lysosomal acetyltransferase that acetylates the non-reducing terminal alpha-glucosamine
	residue of intralysosomal heparin or heparan sulfate, converting it into a substrate for luminal
	alpha-N-acetyl glucosaminidase. {EC0:0000269 PubMed:16960811,

### **Target Details**

Expiry Date:

69 PubMed:17033958, ECO:0000269 PubMed:19823584,
69 PubMed:20650889}.
cluding tag.
noglycan Metabolic Process
to the applications listed above we expect the protein to work for functional studies
the protein has not been tested for functional studies yet we cannot offer a
though.
which it is highly likely that the recombinant protein with the default tag will be
ur protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to
olubility. We will discuss all possible options with you in detail to assure that you
r protein of interest.
ch Use only
aCL, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer.
ated freeze-thaw cycles.
9°C.

Unlimited (if stored properly)