# antibodies -online.com





# Pejvakin Protein (AA 1-352) (Strep Tag)



**Image** 



Go to Product page

#### Overview

Quantity:	1 mg
Target:	Pejvakin (DFNB59)
Protein Characteristics:	AA 1-352
Origin:	Human
Source:	Tobacco (Nicotiana tabacum)
Protein Type:	Recombinant
Purification tag / Conjugate:	This Pejvakin protein is labelled with Strep Tag.
Application:	ELISA, SDS-PAGE (SDS), Western Blotting (WB)

#### **Product Details**

Sequence:

MFAAATKSFV KQVGDGGRLV PVPSLSEADK YQPLSLVVKK KRCFLFPRYK FTSTPFTLKD ILLGDREISA GISSYQLLNY EDESDVSLYG RRGNHIVNDV GINVAGSDSI AVKASFGIVT KHEVEVSTLL KEITTRKINF DHSLIRQSRS SRKAVLCVVM ESIRTTRQCS LSVHAGIRGE AMRFHFMDEQ NPKGRDKAIV FPAHTTIAFS VFELFIYLDG AFDLCVTSVS KGGFEREETA TFALLYRLRN ILFERNRRVM DVISRSQLYL DDLFSDYYDK PLSMTDISLK EGTHIRVNLL NHNIPKGPCI LCGMGNFKRE TVYGCFQCSV DGQKYVRLHA VPCFDIWHKR MK

Sequence without tag. The proposed Strep-Tag is based on experience s with the expression system, a different complexity of the protein could make another tag necessary. In case you have a special request, please contact us.

Characteristics:

Key Benefits:

- Made in Germany from design to production by highly experienced protein experts.
- · Protein expressed with ALiCE® and purified by multi-step, protein-specific process to ensure

correct folding and modification.

- These proteins are normally active (enzymatically functional) as our customers have reported (not tested by us and not guaranteed).
- State-of-the-art algorithm used for plasmid design (Gene synthesis).

This protein is a **made-to-order protein** and will be made for the first time for your order. Our experts in the lab will ensure that you receive a correctly folded protein.

The big advantage of ordering our **made-to-order proteins** in comparison to ordering custom made proteins from other companies is that there is no financial obligation in case the protein cannot be expressed or purified.

#### Expression System:

- ALiCE®, our Almost Living Cell-Free Expression System is based on a lysate obtained from Nicotiana tabacum c.v.. This contains all the protein expression machinery needed to produce even the most difficult-to-express proteins, including those that require posttranslational modifications.
- During lysate production, the cell wall and other cellular components that are not required for
  protein production are removed, leaving only the protein production machinery and the
  mitochondria to drive the reaction. During our lysate completion steps, the additional
  components needed for protein production (amino acids, cofactors, etc.) are added to
  produce something that functions like a cell, but without the constraints of a living system all that's needed is the DNA that codes for the desired protein!

#### Concentration:

- The concentration of our recombinant proteins is measured using the absorbance at 280nm.
- The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.
- We use the Expasy's ProtParam tool to determine the absorption coefficient of each protein.

#### Purification:

Two step purification of proteins expressed in Almost Living Cell-Free Expression System (ALiCE®):

- 1. In a first purification step, the protein is purified from the cleared cell lysate using StrepTag capture material. Eluate fractions are analyzed by SDS-PAGE.
- Protein containing fractions of the best purification are subjected to second purification step through size exclusion chromatography. Eluate fractions are analyzed by SDS-PAGE and Western blot.

Purity:

>80 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.

# **Product Details** Low Endotoxin less than 1 EU/mg (< 0.1 ng/mg) Endotoxin Level: Grade: Crystallography grade Target Details Target: Pejvakin (DFNB59) PJVK (DFNB59 Products) Alternative Name: Background: Pejvakin (Autosomal recessive deafness type 59 protein), FUNCTION: Peroxisome-associated protein required to protect auditory hair cells against noise-induced damage. Acts by regulating noise-induced peroxisome proliferation in auditory hair cells and neurons, and promoting autophagic degradation of damaged peroxisomes (pexophagy). Noise overexposure increases reactive oxygen species (ROS) levels, causing oxidative damage to auditory hair cells and resulting in hearing loss. PJVK acts as a ROS sensor that recruits the autophagy machinery to trigger pexophagy of peroxisomes damaged by oxidative stress. In addition to pexophagy, also required to promote peroxisome proliferation in response to sound overstimulation. {ECO:0000250|UniProtKB:Q0ZLH2}. Molecular Weight: 39.9 kDa UniProt: Q0ZLH3 Pathways: Sensory Perception of Sound **Application Details** In addition to the applications listed above we expect the protein to work for functional studies **Application Notes:** as well. As the protein has not been tested for functional studies yet we cannot offer a guarantee though. Comment: ALiCE®, our Almost Living Cell-Free Expression System is based on a lysate obtained from Nicotiana tabacum c.v.. This contains all the protein expression machinery needed to produce even the most difficult-to-express proteins, including those that require post-translational modifications. During lysate production, the cell wall and other cellular components that are not required for protein production are removed, leaving only the protein production machinery and the mitochondria to drive the reaction. During our lysate completion steps, the additional


components needed for protein production (amino acids, cofactors, etc.) are added to produce

something that functions like a cell, but without the constraints of a living system - all that's

## **Application Details**

	needed is the DNA that codes for the desired protein!
Restrictions:	For Research Use only
Handling	
Format:	Liquid
Buffer:	The buffer composition is at the discretion of the manufacturer. If you have a special request, please contact us.
Handling Advice:	Avoid repeated freeze-thaw cycles.
Storage:	-80 °C
Storage Comment:	Store at -80°C.
Expiry Date:	Unlimited (if stored properly)

### Images



**Image 1.** "Crystallography Grade" protein due to multi-step, protein-specific purification process