# antibodies .- online.com





# LACTBL1 Protein (AA 1-500) (Strep Tag)



**Image** 



Go to Product page

#### Overview

Quantity:	1 mg
Target:	LACTBL1
Protein Characteristics:	AA 1-500
Origin:	Human
Source:	Tobacco (Nicotiana tabacum)
Protein Type:	Recombinant
Purification tag / Conjugate:	This LACTBL1 protein is labelled with Strep Tag.
Application:	ELISA, SDS-PAGE (SDS), Western Blotting (WB)

#### **Product Details**

Sequence:

MCPRHPEPVP LAHPLPVLKE ALEKVDQILR QAMSAPGVAA MSAVVIHNDT VLWTGNFGKK NGSDPASGAP NEYTMYRISS ISKIFPVLML YRLWEEGIVA SLDDPLERYA STFTINNPLG LASAEQQGLI LRRMASQLSG LPRRLRSTSL LWKGSTQEAL NLLKDDVLVV DPGTRCHYST LAFSLLAHVL AAHTAQGDYQ RWVSENVLEP LGMADTGFDL TPDVRARLAA GFYGSGRPAP LYDLGWYRPS GQMYSTAADL AKLAVALLGG GPRRLLRPDA AKTLLAPLLA CPGAYFANET GTPWEFHAQR GYRVVRKDGD LDGYAATFSL VPPLRLGLVL LLAGPRPPGP DLVARAYDEL LPALERALRE AEPGPAPPPT AHPFAGYFTF ANLTFYEVRA GPAGELRLRQ FGPRVEALVP PAFRTLALRH LHGRVFQLHV AHEFPCALPL GDAWLSLEAQ HGQLVNFYPL DHHGLSPGFD VPGLNTYRVL RLRGKPVFKT

Sequence without tag. The proposed Strep-Tag is based on experience s with the expression system, a different complexity of the protein could make another tag necessary. In case you have a special request, please contact us.

Characteristics:

#### Key Benefits:

- Made in Germany from design to production by highly experienced protein experts.
- Protein expressed with ALiCE® and purified by multi-step, protein-specific process to ensure correct folding and modification.
- These proteins are normally active (enzymatically functional) as our customers have reported (not tested by us and not guaranteed).
- State-of-the-art algorithm used for plasmid design (Gene synthesis).

This protein is a **made-to-order protein** and will be made for the first time for your order. Our experts in the lab will ensure that you receive a correctly folded protein.

The big advantage of ordering our **made-to-order proteins** in comparison to ordering custom made proteins from other companies is that there is no financial obligation in case the protein cannot be expressed or purified.

#### Expression System:

- ALiCE®, our Almost Living Cell-Free Expression System is based on a lysate obtained from Nicotiana tabacum c.v.. This contains all the protein expression machinery needed to produce even the most difficult-to-express proteins, including those that require posttranslational modifications.
- During lysate production, the cell wall and other cellular components that are not required for
  protein production are removed, leaving only the protein production machinery and the
  mitochondria to drive the reaction. During our lysate completion steps, the additional
  components needed for protein production (amino acids, cofactors, etc.) are added to
  produce something that functions like a cell, but without the constraints of a living system all that's needed is the DNA that codes for the desired protein!

#### Concentration:

- The concentration of our recombinant proteins is measured using the absorbance at 280nm.
- The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.
- We use the Expasy's ProtParam tool to determine the absorption coefficient of each protein.

Purification:

Two step purification of proteins expressed in Almost Living Cell-Free Expression System (ALiCE®):

- 1. In a first purification step, the protein is purified from the cleared cell lysate using StrepTag capture material. Eluate fractions are analyzed by SDS-PAGE.
- 2. Protein containing fractions of the best purification are subjected to second purification step through size exclusion chromatography. Eluate fractions are analyzed by SDS-PAGE and


### **Product Details**

Purity: >80 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Bloendotoxin Level:  Low Endotoxin less than 1 EU/mg (< 0.1 ng/mg)  Grade:  Crystallography grade  Farget Details  Farget:  LACTBL1  ACTBL1 (LACTBL1 Products)  Background:  Putative beta-lactamase-like 1  Molecular Weight:  54.7 kDa  Application Details  Application Notes:  In addition to the applications listed above we expect the protein to work for functional as well. As the protein has not been tested for functional studies yet we cannot offer a guarantee though.  Comment:  ALICE®, our Almost Living Cell-Free Expression System is based on a lysate obtained Nicotiana tabacum c.v This contains all the protein expression machinery needed to even the most difficult-to-express proteins, including those that require post-translation modifications.  During lysate production, the cell wall and other cellular components that are not required protein production are removed, leaving only the protein production machinery and the mitochondria to drive the reaction. During our lysate completion steps, the additional components needed for protein production (amino acids, cofactors, etc.) are added to	
Farget Details  Farget: LACTBL1  Alternative Name: LACTBL1 (LACTBL1 Products)  Background: Putative beta-lactamase-like 1  Molecular Weight: 54.7 kDa  Application Details  Application Notes: In addition to the applications listed above we expect the protein to work for functiona as well. As the protein has not been tested for functional studies yet we cannot offer a guarantee though.  Comment: ALICE®, our Almost Living Cell-Free Expression System is based on a lysate obtained Nicotiana tabacum c.v This contains all the protein expression machinery needed to even the most difficult-to-express proteins, including those that require post-translation modifications.  During lysate production, the cell wall and other cellular components that are not require protein production are removed, leaving only the protein production machinery and the mitochondria to drive the reaction. During our lysate completion steps, the additional	
Farget: LACTBL1  Alternative Name: LACTBL1 (LACTBL1 Products)  Background: Putative beta-lactamase-like 1  Molecular Weight: 54.7 kDa  Application Details  Application Notes: In addition to the applications listed above we expect the protein to work for functional as well. As the protein has not been tested for functional studies yet we cannot offer a guarantee though.  Comment: ALICE®, our Almost Living Cell-Free Expression System is based on a lysate obtained Nicotiana tabacum c.v This contains all the protein expression machinery needed to even the most difficult-to-express proteins, including those that require post-translation modifications.  During lysate production, the cell wall and other cellular components that are not require protein production machinery and the mitochondria to drive the reaction. During our lysate completion steps, the additional	
Alternative Name:  LACTBL1 (LACTBL1 Products)  Background:  Putative beta-lactamase-like 1  Molecular Weight:  54.7 kDa  ABMY62  Application Details  Application Notes:  In addition to the applications listed above we expect the protein to work for functional as well. As the protein has not been tested for functional studies yet we cannot offer a guarantee though.  Comment:  ALICE®, our Almost Living Cell-Free Expression System is based on a lysate obtained Nicotiana tabacum c.v This contains all the protein expression machinery needed to even the most difficult-to-express proteins, including those that require post-translation modifications.  During lysate production, the cell wall and other cellular components that are not require protein production are removed, leaving only the protein production machinery and the mitochondria to drive the reaction. During our lysate completion steps, the additional	
Alternative Name:  LACTBL1 (LACTBL1 Products)  Putative beta-lactamase-like 1  Sackground:  Putative beta-lactamase-like 1  S4.7 kDa  Application Details  Application Details  In addition to the applications listed above we expect the protein to work for functions as well. As the protein has not been tested for functional studies yet we cannot offer a guarantee though.  Comment:  ALICE®, our Almost Living Cell-Free Expression System is based on a lysate obtained Nicotiana tabacum c.v This contains all the protein expression machinery needed to even the most difficult-to-express proteins, including those that require post-translatic modifications.  During lysate production, the cell wall and other cellular components that are not require protein production are removed, leaving only the protein production machinery and the mitochondria to drive the reaction. During our lysate completion steps, the additional	
Background:  Putative beta-lactamase-like 1  Molecular Weight:  54.7 kDa  JniProt:  A8MY62  Application Details  In addition to the applications listed above we expect the protein to work for functions as well. As the protein has not been tested for functional studies yet we cannot offer a guarantee though.  Comment:  ALICE®, our Almost Living Cell-Free Expression System is based on a lysate obtained Nicotiana tabacum c.v This contains all the protein expression machinery needed to even the most difficult-to-express proteins, including those that require post-translation modifications.  During lysate production, the cell wall and other cellular components that are not require protein production are removed, leaving only the protein production machinery and the mitochondria to drive the reaction. During our lysate completion steps, the additional	
Molecular Weight: 54.7 kDa  Application Details  Application Notes: In addition to the applications listed above we expect the protein to work for functions as well. As the protein has not been tested for functional studies yet we cannot offer a guarantee though.  Comment: ALiCE®, our Almost Living Cell-Free Expression System is based on a lysate obtained Nicotiana tabacum c.v This contains all the protein expression machinery needed to even the most difficult-to-express proteins, including those that require post-translation modifications.  During lysate production, the cell wall and other cellular components that are not require protein production are removed, leaving only the protein production machinery and the mitochondria to drive the reaction. During our lysate completion steps, the additional	
Application Details  Application Notes:  In addition to the applications listed above we expect the protein to work for functional as well. As the protein has not been tested for functional studies yet we cannot offer a guarantee though.  Comment:  ALICE®, our Almost Living Cell-Free Expression System is based on a lysate obtained Nicotiana tabacum c.v This contains all the protein expression machinery needed to even the most difficult-to-express proteins, including those that require post-translation modifications.  During lysate production, the cell wall and other cellular components that are not require protein production are removed, leaving only the protein production machinery and the mitochondria to drive the reaction. During our lysate completion steps, the additional	
Application Details  Application Notes:  In addition to the applications listed above we expect the protein to work for functions as well. As the protein has not been tested for functional studies yet we cannot offer a guarantee though.  Comment:  ALICE®, our Almost Living Cell-Free Expression System is based on a lysate obtained Nicotiana tabacum c.v This contains all the protein expression machinery needed to even the most difficult-to-express proteins, including those that require post-translation modifications.  During lysate production, the cell wall and other cellular components that are not require protein production are removed, leaving only the protein production machinery and the mitochondria to drive the reaction. During our lysate completion steps, the additional	
Application Notes:  In addition to the applications listed above we expect the protein to work for functional as well. As the protein has not been tested for functional studies yet we cannot offer a guarantee though.  ALICE®, our Almost Living Cell-Free Expression System is based on a lysate obtained Nicotiana tabacum c.v This contains all the protein expression machinery needed to even the most difficult-to-express proteins, including those that require post-translation modifications.  During lysate production, the cell wall and other cellular components that are not require protein production are removed, leaving only the protein production machinery and the mitochondria to drive the reaction. During our lysate completion steps, the additional	
as well. As the protein has not been tested for functional studies yet we cannot offer a guarantee though.  ALICE®, our Almost Living Cell-Free Expression System is based on a lysate obtained Nicotiana tabacum c.v This contains all the protein expression machinery needed to even the most difficult-to-express proteins, including those that require post-translatic modifications.  During lysate production, the cell wall and other cellular components that are not require protein production are removed, leaving only the protein production machinery and the mitochondria to drive the reaction. During our lysate completion steps, the additional	
Nicotiana tabacum c.v This contains all the protein expression machinery needed to even the most difficult-to-express proteins, including those that require post-translation modifications.  During lysate production, the cell wall and other cellular components that are not require protein production are removed, leaving only the protein production machinery and the mitochondria to drive the reaction. During our lysate completion steps, the additional	
something that functions like a cell, but without the constraints of a living system - all needed is the DNA that codes for the desired protein!	oroduce nal red for e
Restrictions: For Research Use only	
Handling	
Format: Liquid	
Buffer: The buffer composition is at the discretion of the manufacturer. If you have a special	equest,

# Handling

	please contact us.
Handling Advice:	Avoid repeated freeze-thaw cycles.
Storage:	-80 °C
Storage Comment:	Store at -80°C.
Expiry Date:	Unlimited (if stored properly)

## **Images**



**Image 1.** "Crystallography Grade" protein due to multi-step, protein-specific purification process