
# antibodies .- online.com





# AGL Protein (AA 1-1532) (Strep Tag)



**Image** 



Go to Product page

#### Overview

Quantity:	1 mg
Target:	AGL
Protein Characteristics:	AA 1-1532
Origin:	Human
Source:	Tobacco (Nicotiana tabacum)
Protein Type:	Recombinant
Purification tag / Conjugate:	This AGL protein is labelled with Strep Tag.
Application:	SDS-PAGE (SDS), Western Blotting (WB), ELISA

## **Product Details**

Sequence:

MGHSKQIRIL LLNEMEKLEK TLFRLEQGYE LQFRLGPTLQ GKAVTVYTNY PFPGETFNRE
KFRSLDWENP TEREDDSDKY CKLNLQQSGS FQYYFLQGNE KSGGGYIVVD PILRVGADNH
VLPLDCVTLQ TFLAKCLGPF DEWESRLRVA KESGYNMIHF TPLQTLGLSR SCYSLANQLE
LNPDFSRPNR KYTWNDVGQL VEKLKKEWNV ICITDVVYNH TAANSKWIQE HPECAYNLVN
SPHLKPAWVL DRALWRFSCD VAEGKYKEKG IPALIENDHH MNSIRKIIWE DIFPKLKLWE
FFQVDVNKAV EQFRRLLTQE NRRVTKSDPN QHLTIIQDPE YRRFGCTVDM NIALTTFIPH
DKGPAAIEEC CNWFHKRMEE LNSEKHRLIN YHQEQAVNCL LGNVFYERLA GHGPKLGPVT
RKHPLVTRYF TFPFEEIDFS MEESMIHLPN KACFLMAHNG WVMGDDPLRN FAEPGSEVYL
RRELICWGDS VKLRYGNKPE DCPYLWAHMK KYTEITATYF QGVRLDNCHS TPLHVAEYML
DAARNLQPNL YVVAELFTGS EDLDNVFVTR LGISSLIREA MSAYNSHEEG RLVYRYGGEP
VGSFVQPCLR PLMPAIAHAL FMDITHDNEC PIVHRSAYDA LPSTTIVSMA CCASGSTRGY
DELVPHQISV VSEERFYTKW NPEALPSNTG EVNFQSGIIA ARCAISKLHQ ELGAKGFIQV

YVDQVDEDIV AVTRHSPSIH QSVVAVSRTA FRNPKTSFYS KEVPQMCIPG KIEEVVLEAR
TIERNTKPYR KDENSINGTP DITVEIREHI QLNESKIVKQ AGVATKGPNE YIQEIEFENL SPGSVIIFRV
SLDPHAQVAV GILRNHLTQF SPHFKSGSLA VDNADPILKI PFASLASRLT LAELNQILYR
CESEEKEDGG GCYDIPNWSA LKYAGLQGLM SVLAEIRPKN DLGHPFCNNL RSGDWMIDYV
SNRLISRSGT IAEVGKWLQA MFFYLKQIPR YLIPCYFDAI LIGAYTTLLD TAWKQMSSFV
QNGSTFVKHL SLGSVQLCGV GKFPSLPILS PALMDVPYRL NEITKEKEQC CVSLAAGLPH
FSSGIFRCWG RDTFIALRGI LLITGRYVEA RNIILAFAGT LRHGLIPNLL GEGIYARYNC
RDAVWWWLQC IQDYCKMVPN GLDILKCPVS RMYPTDDSAP LPAGTLDQPL FEVIQEAMQK
HMQGIQFRER NAGPQIDRNM KDEGFNITAG VDEETGFVYG GNRFNCGTWM DKMGESDRAR
NRGIPATPRD GSAVEIVGLS KSAVRWLLEL SKKNIFPYHE VTVKRHGKAI KVSYDEWNRK
IQDNFEKLFH VSEDPSDLNE KHPNLVHKRG IYKDSYGASS PWCDYQLRPN FTIAMVVAPE
LFTTEKAWKA LEIAEKKLLG PLGMKTLDPD DMVYCGIYDN ALDNDNYNLA KGFNYHQGPE
WLWPIGYFLR AKLYFSRLMG PETTAKTIVL VKNVLSRHYV HLERSPWKGL PELTNENAQY
CPFSCETOAW SIATILETLY DL

Sequence without tag. The proposed Strep-Tag is based on experience s with the expression system, a different complexity of the protein could make another tag necessary. In case you have a special request, please contact us.

#### Characteristics:

#### Key Benefits:

- Made in Germany from design to production by highly experienced protein experts.
- Protein expressed with ALiCE® and purified by multi-step, protein-specific process to ensure correct folding and modification.
- These proteins are normally active (enzymatically functional) as our customers have reported (not tested by us and not guaranteed).
- State-of-the-art algorithm used for plasmid design (Gene synthesis).

This protein is a **made-to-order protein** and will be made for the first time for your order. Our experts in the lab will ensure that you receive a correctly folded protein.

The big advantage of ordering our **made-to-order proteins** in comparison to ordering custom made proteins from other companies is that there is no financial obligation in case the protein cannot be expressed or purified.

## Expression System:

 ALiCE®, our Almost Living Cell-Free Expression System is based on a lysate obtained from Nicotiana tabacum c.v.. This contains all the protein expression machinery needed to produce even the most difficult-to-express proteins, including those that require posttranslational modifications. During lysate production, the cell wall and other cellular components that are not required for
protein production are removed, leaving only the protein production machinery and the
mitochondria to drive the reaction. During our lysate completion steps, the additional
components needed for protein production (amino acids, cofactors, etc.) are added to
produce something that functions like a cell, but without the constraints of a living system all that's needed is the DNA that codes for the desired protein!

#### Concentration:

- The concentration of our recombinant proteins is measured using the absorbance at 280nm.
- The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.
- · We use the Expasy's ProtParam tool to determine the absorption coefficient of each protein.

#### Purification:

Two step purification of proteins expressed in Almost Living Cell-Free Expression System (ALiCE®):

- 1. In a first purification step, the protein is purified from the cleared cell lysate using StrepTag capture material. Eluate fractions are analyzed by SDS-PAGE.
- Protein containing fractions of the best purification are subjected to second purification step through size exclusion chromatography. Eluate fractions are analyzed by SDS-PAGE and Western blot.

Purity:

>80 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.

Endotoxin Level:

Low Endotoxin less than 1 EU/mg (< 0.1 ng/mg)

Grade:

Crystallography grade

## **Target Details**

Target:	AGL
Alternative Name:	AGL (AGL Products)
Background:	Glycogen debranching enzyme (Glycogen debrancher) [Includes: 4-alpha-glucanotransferase
	(EC 2.4.1.25) (Oligo-1,4-1,4-glucantransferase), Amylo-alpha-1,6-glucosidase (Amylo-1,6-
	glucosidase) (EC 3.2.1.33) (Dextrin 6-alpha-D-glucosidase)],FUNCTION: Multifunctional enzyme
	acting as 1,4-alpha-D-glucan:1,4-alpha-D-glucan 4-alpha-D-glycosyltransferase and amylo-1,6-
	glucosidase in glycogen degradation.
Molecular Weight:	174.8 kDa
UniProt:	P35573

# Target Details

Expiry Date:

Target Details	
Pathways:	Cellular Glucan Metabolic Process
Application Details	
Application Notes:	In addition to the applications listed above we expect the protein to work for functional studies as well. As the protein has not been tested for functional studies yet we cannot offer a guarantee though.
Comment:	ALiCE®, our Almost Living Cell-Free Expression System is based on a lysate obtained from Nicotiana tabacum c.v This contains all the protein expression machinery needed to produce even the most difficult-to-express proteins, including those that require post-translational modifications.  During lysate production, the cell wall and other cellular components that are not required for protein production are removed, leaving only the protein production machinery and the mitochondria to drive the reaction. During our lysate completion steps, the additional components needed for protein production (amino acids, cofactors, etc.) are added to produce something that functions like a cell, but without the constraints of a living system - all that's needed is the DNA that codes for the desired protein!
Restrictions:	For Research Use only
Handling	
Format:	Liquid
Buffer:	The buffer composition is at the discretion of the manufacturer. If you have a special request, please contact us.
Handling Advice:	Avoid repeated freeze-thaw cycles.
Storage:	-80 °C
Storage Comment:	Store at -80°C.

Unlimited (if stored properly)



**Image 1.** "Crystallography Grade" protein due to multi-step, protein-specific purification process