antibodies .- online.com





GLG1 Protein (AA 30-1145) (His tag)



Image



Go to Product page

Overview

Quantity:	1 mg
Target:	GLG1
Protein Characteristics:	AA 30-1145
Origin:	Human
Source:	Insect Cells
Protein Type:	Recombinant
Purification tag / Conjugate:	This GLG1 protein is labelled with His tag.
Application:	SDS-PAGE (SDS), Western Blotting (WB), ELISA, Crystallization (Crys)

Product Details

Sequence:

LPGQGVHSQG QGPGANFVSF VGQAGGGGPA GQQLPQLPQS SQLQQQQQQ QQQQQPQPPQ
PPFPAGGPPA RRGGAGAGGG WKLAEEESCR EDVTRVCPKH TWSNNLAVLE CLQDVREPEN
EISSDCNHLL WNYKLNLTTD PKFESVAREV CKSTITEIKE CADEPVGKGY MVSCLVDHRG
NITEYQCHQY ITKMTAIIFS DYRLICGFMD DCKNDINILK CGSIRLGEKD AHSQGEVVSC
LEKGLVKEAE EREPKIQVSE LCKKAILRVA ELSSDDFHLD RHLYFACRDD RERFCENTQA
GEGRVYKCLF NHKFEESMSE KCREALTTRQ KLIAQDYKVS YSLAKSCKSD LKKYRCNVEN
LPRSREARLS YLLMCLESAV HRGRQVSSEC QGEMLDYRRM LMEDFSLSPE IILSCRGEIE
HHCSGLHRKG RTLHCLMKVV RGEKGNLGMN CQQALQTLIQ ETDPGADYRI DRALNEACES
VIQTACKHIR SGDPMILSCL MEHLYTEKMV EDCEHRLLEL QYFISRDWKL DPVLYRKCQG
DASRLCHTHG WNETSEFMPQ GAVFSCLYRH AYRTEEQGRR LSRECRAEVQ RILHQRAMDV
KLDPALQDKC LIDLGKWCSE KTETGQELEC LQDHLDDLVV ECRDIVGNLT ELESEDIQIE
ALLMRACEPI IQNFCHDVAD NQIDSGDLME CLIQNKHQKD MNEKCAIGVT HFQLVQMKDF

RFSYKFKMAC KEDVLKLCPN IKKKVDVVIC LSTTVRNDTL QEAKEHRVSL KCRRQLRVEE
LEMTEDIRLE PDLYEACKSD IKNFCSAVQY GNAQIIECLK ENKKQLSTRC HQKVFKLQET
EMMDPELDYT LMRVCKQMIK RFCPEADSKT MLQCLKQNKN SELMDPKCKQ MITKRQITQN
TDYRLNPMLR KACKADIPKF CHGILTKAKD DSELEGQVIS CLKLRYADQR LSSDCEDQIR
IIIQESALDY RLDPQLQLHC SDEISSLCAE EAAAQEQTGQ VEECLKVNLL KIKTELCKKE
VLNMLKESKA DIFVDPVLHT ACALDIKHHC AAITPGRGRQ MSCLMEALED KRVRLQPECK
KRLNDRIEMW SYAAKVAPAD GFSDLAMOVM TSPSKN

Sequence without tag. Tag location is at the discretion of the manufacturer. If you have a special request, please contact us.

Characteristics:

- Made in Germany from design to production by highly experienced protein experts.
- Human GLG1 Protein (raised in Insect Cells) purified by multi-step, protein-specific process to ensure crystallization grade.
- State-of-the-art algorithm used for plasmid design (Gene synthesis).

This protein is a made to order protein and will be made for the first time for your order. Our experts in the lab will ensure that you receive a correctly folded protein.

The big advantage of ordering our made-to-order proteins in comparison to ordering custom made proteins from other companies is that there is no financial obligation in case the protein cannot be expressed or purified.

In the unlikely event that the protein cannot be expressed or purified we do not charge anything (other companies might charge you for any performed steps in the expression process for custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression experiments or purification optimization).

When you order this made-to-order protein you will only pay upon receival of the correctly folded protein. With no financial risk on your end you can rest assured that our experienced protein experts will do everything to make sure that you receive the protein you ordered. The concentration of our recombinant proteins is measured using the absorbance at 280nm. The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.

The concentration of the protein is calculated using its specific absorption coefficient. We use the Expasy's protparam tool to determine the absorption coefficient of each protein.

Purification:

Two step purification of proteins expressed in baculovirus infected SF9 insect cells:

- 1. In a first purification step, the protein is purified from the cleared cell lysate using three different His-tag capture materials: high yield, EDTA resistant, or DTT resistant. Eluate fractions are analyzed by SDS-PAGE.
- 2. Protein containing fractions of the best purification are subjected to second purification step

Product Details

Sterility: 0.22 µm filtered Endotoxin Level: Protein is endotoxin free. Crystallography grade Target Details Target: GLG1 Alternative Name: GLG1 (SLG1 Products) Blackground: Binds fibroblast growth factor and E-selectin (cell-adhesion lectin on endothelial cells mediating the binding of neutrophils). Molecular Weight: 128.5 kDa Including tag. UniProt: Q97896 Application Details Application Notes: In addition to the applications listed above we expect the protein to work for functional studies as well. As the protein has not been tested for functional studies yet we cannot offer a gurantee though. Comment: In cases in which it is highly likely that the recombinant protein with the default tag will be insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible options with you in detail to assure that you receive your protein of interest. Restrictions: For Research Use only Handling Format: Liquid Buffer: 100 mM NaCL 20 mM Hepes, 10% glycerol, pH value is at the discretion of the manufacturer. Handling Advice: Avoid repeated freeze-thaw cycles.	1 Toddot Detaile	
Sterility: 0.22 µm filtered Endotoxin Level: Protein is endotoxin free. Crystallography grade Target Details Target: GLG1 Alternative Name: GLG1 (SLG1 Products) Blackground: Binds fibroblast growth factor and E-selectin (cell-adhesion lectin on endothelial cells mediating the binding of neutrophils). Molecular Weight: 128.5 kDa Including tag. UniProt: Q97896 Application Details Application Notes: In addition to the applications listed above we expect the protein to work for functional studies as well. As the protein has not been tested for functional studies yet we cannot offer a gurantee though. Comment: In cases in which it is highly likely that the recombinant protein with the default tag will be insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible options with you in detail to assure that you receive your protein of interest. Restrictions: For Research Use only Handling Format: Liquid Buffer: 100 mM NaCL 20 mM Hepes, 10% glycerol, pH value is at the discretion of the manufacturer. Handling Advice: Avoid repeated freeze-thaw cycles.		
Endotoxin Level: Protein is endotoxin free. Grade: Crystallography grade Target Details Target: GLG1 Alternative Name: GLG1 (GLG1 Products) Background: Binds fibroblast growth factor and E-selectin (cell-adhesion lectin on endothelial cells mediating the binding of neutrophilis). Molecular Weight: 128.5 kDa Including tag. UniProt: Q97896 Application Details Application Notes: In addition to the applications listed above we expect the protein to work for functional studies as well. As the protein has not been tested for functional studies yet we cannot offer a gurantee though. Comment: In cases in which it is highly likely that the recombinant protein with the default tag will be insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible options with you in detail to assure that you receive your protein of interest. Restrictions: For Research Use only Handling Format: Liquid Buffer: 100 mM NaCL, 20 mM Hepes, 10% glycerol, pH value is at the discretion of the manufacturer. Avoid repeated freeze-thaw cycles.	Purity:	>95 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.
Target Details Target: GLG1 Alternative Name: GLG1 (GLG1 Products) Background: Binds fibroblast growth factor and E-selectin (cell-adhesion lectin on endothelial cells mediating the binding of neutrophils). Molecular Weight: 128.5 kDa Including tag. UniProt: Q92896 Application Details Application Notes: In addition to the applications listed above we expect the protein to work for functional studies as well. As the protein has not been tested for functional studies yet we cannot offer a gurantee though. Comment: In cases in which it is highly likely that the recombinant protein with the default tag will be insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible options with you in detail to assure that you receive your protein of interest. Restrictions: For Research Use only Handling Format: Liquid Buffer: 100 mM NaCL, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer. Avoid repeated freeze-thaw cycles.	Sterility:	0.22 µm filtered
Target Details Target: GLG1 Alternative Name: GLG1 (GLG1 Products) Background: Binds fibroblast growth factor and E-selectin (cell-adhesion lectin on endothelial cells mediating the binding of neutrophils). Molecular Weight: 128.5 kDa Including tag. UniProt: Q92896 Application Details Application Notes: In addition to the applications listed above we expect the protein to work for functional studies as well. As the protein has not been tested for functional studies yet we cannot offer a guranter though. Comment: In cases in which it is highly likely that the recombinant protein with the default tag will be insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible options with you in detail to assure that you receive your protein of interest. Restrictions: For Research Use only Handling Format: Liquid Buffer: 100 mM NaCL, 20 mM Hepes, 10% glycerol, pH value is at the discretion of the manufacturer. Handling Advice: Avoid repeated freeze-thaw cycles.	Endotoxin Level:	Protein is endotoxin free.
Target: GLG1 Alternative Name: GLG1 (GLG1 Products) Background: Binds fibroblast growth factor and E-selectin (cell-adhesion lectin on endothelial cells mediating the binding of neutrophils). Molecular Weight: 128.5 kDa Including tag. UniProt: 092896 Application Details Application Notes: In addition to the applications listed above we expect the protein to work for functional studies as well. As the protein has not been tested for functional studies yet we cannot offer a guranter though. Comment: In cases in which it is highly likely that the recombinant protein with the default tag will be insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible options with you in detail to assure that you receive your protein of interest. Restrictions: For Research Use only Handling Format: Liquid Buffer: 100 mM NaCL, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer. Avoid repeated freeze-thaw cycles.	Grade:	Crystallography grade
Alternative Name: GLG1 (GLG1 Products) Binds fibroblast growth factor and E-selectin (cell-adhesion lectin on endothelial cells mediating the binding of neutrophils). Molecular Weight: 128.5 kDa Including tag. UniProt: O92896 Application Details Application Notes: In addition to the applications listed above we expect the protein to work for functional studies as well. As the protein has not been tested for functional studies yet we cannot offer a gurantee though. Comment: In cases in which it is highly likely that the recombinant protein with the default tag will be insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible options with you in detail to assure that you receive your protein of interest. Restrictions: For Research Use only Handling Format: Liquid 100 mM NaCL, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer. Avoid repeated freeze-thaw cycles.	Target Details	
Binds fibroblast growth factor and E-selectin (cell-adhesion lectin on endothelial cells mediating the binding of neutrophils). Molecular Weight: 128.5 kDa Including tag. UniProt: 092896 Application Details Application Notes: In addition to the applications listed above we expect the protein to work for functional studies as well. As the protein has not been tested for functional studies yet we cannot offer a gurantee though. Comment: In cases in which it is highly likely that the recombinant protein with the default tag will be insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible options with you in detail to assure that you receive your protein of interest. Restrictions: For Research Use only Handling Format: Liquid Buffer: 100 mM NaCL, 20 mM Hepes, 10% glycerol, pH value is at the discretion of the manufacturer. Handling Advice: Avoid repeated freeze-thaw cycles.	Target:	GLG1
the binding of neutrophils). Molecular Weight: 128.5 kDa Including tag. UniProt: 092896 Application Details Application Notes: In addition to the applications listed above we expect the protein to work for functional studies as well. As the protein has not been tested for functional studies yet we cannot offer a guranter though. Comment: In cases in which it is highly likely that the recombinant protein with the default tag will be insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible options with you in detail to assure that you receive your protein of interest. Restrictions: For Research Use only Handling Format: Liquid Buffer: 100 mM NaCL, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer. Handling Advice: Avoid repeated freeze-thaw cycles.	Alternative Name:	GLG1 (GLG1 Products)
Application Details Application Notes: In addition to the applications listed above we expect the protein to work for functional studies as well. As the protein has not been tested for functional studies yet we cannot offer a guranter though. Comment: In cases in which it is highly likely that the recombinant protein with the default tag will be insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible options with you in detail to assure that you receive your protein of interest. Restrictions: For Research Use only Handling Format: Liquid Buffer: 100 mM NaCL, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer. Handling Advice: Avoid repeated freeze-thaw cycles.	Background:	Binds fibroblast growth factor and E-selectin (cell-adhesion lectin on endothelial cells mediating the binding of neutrophils).
Application Details Application Notes: In addition to the applications listed above we expect the protein to work for functional studies as well. As the protein has not been tested for functional studies yet we cannot offer a gurantee though. Comment: In cases in which it is highly likely that the recombinant protein with the default tag will be insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible options with you in detail to assure that you receive your protein of interest. Restrictions: For Research Use only Handling Format: Liquid Buffer: 100 mM NaCL, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer. Handling Advice: Avoid repeated freeze-thaw cycles.	Molecular Weight:	128.5 kDa Including tag.
Application Notes: In addition to the applications listed above we expect the protein to work for functional studies as well. As the protein has not been tested for functional studies yet we cannot offer a gurantee though. Comment: In cases in which it is highly likely that the recombinant protein with the default tag will be insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible options with you in detail to assure that you receive your protein of interest. Restrictions: For Research Use only Handling Format: Liquid 100 mM NaCL, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer. Handling Advice: Avoid repeated freeze-thaw cycles.	UniProt:	Q92896
as well. As the protein has not been tested for functional studies yet we cannot offer a gurantee though. Comment: In cases in which it is highly likely that the recombinant protein with the default tag will be insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible options with you in detail to assure that you receive your protein of interest. Restrictions: For Research Use only Handling Format: Liquid Buffer: 100 mM NaCL, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer. Handling Advice: Avoid repeated freeze-thaw cycles.	Application Details	
insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible options with you in detail to assure that you receive your protein of interest. Restrictions: For Research Use only Handling Format: Liquid Buffer: 100 mM NaCL, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer. Handling Advice: Avoid repeated freeze-thaw cycles.	Application Notes:	as well. As the protein has not been tested for functional studies yet we cannot offer a gurantee
Handling Format: Liquid Buffer: 100 mM NaCL, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer. Handling Advice: Avoid repeated freeze-thaw cycles.	Comment:	insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible options with you in detail to assure that you
Format: Liquid Buffer: 100 mM NaCL, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer. Handling Advice: Avoid repeated freeze-thaw cycles.	Restrictions:	For Research Use only
Buffer: 100 mM NaCL, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer. Handling Advice: Avoid repeated freeze-thaw cycles.	Handling	
Handling Advice: Avoid repeated freeze-thaw cycles.	Format:	Liquid
	Buffer:	100 mM NaCL, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer.
Storage: -80 °C	Handling Advice:	Avoid repeated freeze-thaw cycles.
	Storage:	-80 °C

Handling

Storage Comment:	Store at -80°C.
Expiry Date:	Unlimited (if stored properly)

Images




Image 1. "Crystallography Grade" protein due to multi-step, protein-specific purification process