antibodies .- online.com





Importin 5 Protein (IPO5) (AA 2-1097) (His tag)



Overview

Quantity:	1 mg
Target:	Importin 5 (IPO5)
Protein Characteristics:	AA 2-1097
Origin:	Human
Source:	Insect Cells
Protein Type:	Recombinant
Purification tag / Conjugate:	This Importin 5 protein is labelled with His tag.
Application:	Western Blotting (WB), SDS-PAGE (SDS), ELISA, Crystallization (Crys)

Product Details

Sequence:

AAAAAEQQQF YLLLGNLLSP DNVVRKQAEE TYENIPGQSK ITFLLQAIRN TTAAEEARQM

AAVLLRRLLS SAFDEVYPAL PSDVQTAIKS ELLMIIQMET QSSMRKKVCD IAAELARNLI

DEDGNNQWPE GLKFLFDSVS SQNVGLREAA LHIFWNFPGI FGNQQQHYLD VIKRMLVQCM

QDQEHPSIRT LSARATAAFI LANEHNVALF KHFADLLPGF LQAVNDSCYQ NDDSVLKSLV

EIADTVPKYL RPHLEATLQL SLKLCGDTSL NNMQRQLALE VIVTLSETAA AMLRKHTNIV

AQTIPQMLAM MVDLEEDEDW ANADELEDDD FDSNAVAGES ALDRMACGLG GKLVLPMIKE

HIMQMLQNPD WKYRHAGLMA LSAIGEGCHQ QMEGILNEIV NFVLLFLQDP HPRVRYAACN

AVGQMATDFA PGFQKKFHEK VIAALLQTME DQGNQRVQAH AAAALINFTE DCPKSLLIPY

LDNLVKHLHS IMVLKLQELI QKGTKLVLEQ VVTSIASVAD TAEEKFVPYY DLFMPSLKHI

VENAVQKELR LLRGKTIECI SLIGLAVGKE KFMQDASDVM QLLLKTQTDF NDMEDDDPQI

SYMISAWARM CKILGKEFQQ YLPVVMGPLM KTASIKPEVA LLDTQDMENM SDDDGWEFVN

LGDQQSFGIK TAGLEEKSTA CQMLVCYAKE LKEGFVEYTE QVVKLMVPLL KFYFHDGVRV

AAAESMPLLL ECARVRGPEY LTQMWHFMCD ALIKAIGTEP DSDVLSEIMH SFAKCIEVMG DGCLNNEHFE ELGGILKAKL EEHFKNQELR QVKRQDEDYD EQVEESLQDE DDNDVYILTK VSDILHSIFS SYKEKVLPWF EQLLPLIVNL ICPHRPWPDR QWGLCIFDDV IEHCSPASFK YAEYFLRPML QYVCDNSPEV RQAAAYGLGV MAQYGGDNYR PFCTEALPLL VRVIQSADSK TKENVNATEN CISAVGKIMK FKPDCVNVEE VLPHWLSWLP LHEDKEEAVQ TFNYLCDLIE SNHPIVLGPN NTNLPKIFSI IAEGEMHEAI KHEDPCAKRL ANVVRQVQTS GGLWTECIAQ LSPEQQAAIQ ELLNSA

Sequence without tag. Tag location is at the discretion of the manufacturer. If you have a special request, please contact us.

Characteristics:

- Made in Germany from design to production by highly experienced protein experts.
- Human IPO5 Protein (raised in Insect Cells) purified by multi-step, protein-specific process to ensure crystallization grade.
- State-of-the-art algorithm used for plasmid design (Gene synthesis).

This protein is a made to order protein and will be made for the first time for your order. Our experts in the lab will ensure that you receive a correctly folded protein.

The big advantage of ordering our made-to-order proteins in comparison to ordering custom made proteins from other companies is that there is no financial obligation in case the protein cannot be expressed or purified.

In the unlikely event that the protein cannot be expressed or purified we do not charge anything (other companies might charge you for any performed steps in the expression process for custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression experiments or purification optimization).

When you order this made-to-order protein you will only pay upon receival of the correctly folded protein. With no financial risk on your end you can rest assured that our experienced protein experts will do everything to make sure that you receive the protein you ordered. The concentration of our recombinant proteins is measured using the absorbance at 280nm. The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.

The concentration of the protein is calculated using its specific absorption coefficient. We use the Expasy's protparam tool to determine the absorption coefficient of each protein.

Purification:

Two step purification of proteins expressed in baculovirus infected SF9 insect cells:

- 1. In a first purification step, the protein is purified from the cleared cell lysate using three different His-tag capture materials: high yield, EDTA resistant, or DTT resistant. Eluate fractions are analyzed by SDS-PAGE.
- 2. Protein containing fractions of the best purification are subjected to second purification step

	through size exclusion chromatography. Eluate fractions are analyzed by SDS-PAGE and Western blot.
Purity:	>95 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.
Sterility:	0.22 µm filtered
Endotoxin Level:	Protein is endotoxin free.
Grade:	Crystallography grade
Target Details	
Target:	Importin 5 (IPO5)
Alternative Name:	IPO5 (IPO5 Products)
Background:	Functions in nuclear protein import as nuclear transport receptor. Serves as receptor for nuclear localization signals (NLS) in cargo substrates. Is thought to mediate docking of the importin/substrate complex to the nuclear pore complex (NPC) through binding to nucleoporin and the complex is subsequently translocated through the pore by an energy requiring, Randependent mechanism. At the nucleoplasmic side of the NPC, Ran binds to the importin, the importin/substrate complex dissociates and importin is re-exported from the nucleus to the cytoplasm where GTP hydrolysis releases Ran. The directionality of nuclear import is thought to be conferred by an asymmetric distribution of the GTP- and GDP-bound forms of Ran between the cytoplasm and nucleus (By similarity). Mediates the nuclear import of ribosomal proteins RPL23A, RPS7 and RPL5. Binds to a beta-like import receptor binding (BIB) domain of RPL23A. In vitro, mediates nuclear import of H2A, H2B, H3 and H4 histones. In case of HIV-1 infection, binds and mediates the nuclear import of HIV-1 Rev. {ECO:0000250, ECO:0000269 PubMed:9687515}.
Molecular Weight:	124.5 kDa Including tag.
UniProt:	000410
Pathways:	Protein targeting to Nucleus
Application Details	
Application Notes:	In addition to the applications listed above we expect the protein to work for functional studies as well. As the protein has not been tested for functional studies yet we cannot offer a gurantee though.

Application Details

Comment:	In cases in which it is highly likely that the recombinant protein with the default tag will be insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible options with you in detail to assure that you receive your protein of interest.
Restrictions:	For Research Use only
Handling	
Format:	Liquid
Format: Buffer:	Liquid 100 mM NaCL, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer.
	·
Buffer:	100 mM NaCL, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer.
Buffer: Handling Advice:	100 mM NaCL, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer. Avoid repeated freeze-thaw cycles.