# antibodies .- online.com





## Retinoic Acid Early Transcript 1G (RAET1G) (AA 26-218) protein (rho-1D4 tag)



Go to Product page

( )	11/0	K\ /	iew	1
	$\cup$	ועוי	$\square \vee \vee$	ı

1 mg	
Retinoic Acid Early Transcript 1G (RAET1G)	
AA 26-218	
Human	
Insect Cells	
Recombinant	
rho-1D4 tag	
ELISA, Western Blotting (WB), Crystallization (Crys), SDS-PAGE (SDS)	
GLADPHSLCY DITVIPKFRP GPRWCAVQGQ VDEKTFLHYD CGSKTVTPVS PLGKKLNVTT	
AWKAQNPVLR EVVDILTEQL LDIQLENYIP KEPLTLQARM SCEQKAEGHG SGSWQLSFDG	
QIFLLFDSEN RMWTTVHPGA RKMKEKWEND KDMTMSFHYI SMGDCTGWLE DFLMGMDSTL	
EPSAGAPPTM SSG	
Sequence without tag. Tag location is at the discretion of the manufacturer. If you have a	
special request, please contact us.	
<ul> <li>Made in Germany - from design to production - by highly experienced protein experts.</li> <li>Human RAET1G Protein (raised in Insect Cells) purified by multi-step, protein-specific process to ensure crystallization grade.</li> </ul>	
State-of-the-art algorithm used for plasmid design (Gene synthesis).	
This protein is a made to order protein and will be made for the first time for your order. Our	
experts in the lab will ensure that you receive a correctly folded protein.	

The big advantage of ordering our made-to-order proteins in comparison to ordering custom made proteins from other companies is that there is no financial obligation in case the protein cannot be expressed or purified.

In the unlikely event that the protein cannot be expressed or purified we do not charge anything (other companies might charge you for any performed steps in the expression process for custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression experiments or purification optimization).

When you order this made-to-order protein you will only pay upon receival of the correctly folded protein. With no financial risk on your end you can rest assured that our experienced protein experts will do everything to make sure that you receive the protein you ordered.

The concentration of our recombinant proteins is measured using the absorbance at 280nm.

The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.

The concentration of the protein is calculated using its specific absorption coefficient. We use the Expasy's protparam tool to determine the absorption coefficient of each protein.

#### Purification:

Three step purification of membrane proteins expressed in baculovirus infected SF9 insect cells:

- 1. Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized with different detergents (detergent screen). Samples are analyzed by Western blot.
- 2. The best performing detergent is used for solubilization and the proteins are purified via their rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate fractions are analyzed by Western blot.
- Protein containing fractions of the best purification are subjected to second purification step through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and Western blot.

 Purity:
 >95 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.

 Sterility:
 0.22 μm filtered

 Endotoxin Level:
 Protein is endotoxin-free.

 Grade:
 Crystallography grade

### **Target Details**

Target:	Retinoic Acid Early Transcript 1G (RAET1G)	
Alternative Name:	RAET1G (RAET1G Products)	
Background: Acts as a ligand for the KLRK1 receptor and mediates NK cell cytotoxicity via the recep		

Storage Comment:

Expiry Date:

Target Details		
	{ECO:0000269 PubMed:15240696}.	
Molecular Weight:	22.8 kDa Including tag.	
UniProt:	Q6H3X3	
Application Details		
Application Notes:	In addition to the applications listed above we expect the protein to work for functional studies as well. As the protein has not been tested for functional studies yet we cannot offer a guarantee though.	
Comment:	In cases in which it is highly likely that the recombinant protein with the default tag will be insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible options with you in detail to assure that you receive your protein of interest.	
Restrictions:	For Research Use only	
Handling		
Format:	Liquid	
Buffer:	100 mM NaCL, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer.	
Handling Advice:	Avoid repeated freeze-thaw cycles.	
Storage:	-80 °C	

Store at -80°C.

Unlimited (if stored properly)