antibodies.com

Datasheet for ABIN3105228 ADTRP Protein (AA 1-230) (Strep Tag)



Overview

1 mg
ADTRP
AA 1-230
Human
Tobacco (Nicotiana tabacum)
Recombinant
This ADTRP protein is labelled with Strep Tag.
ELISA, SDS-PAGE (SDS), Western Blotting (WB)
MTKTSTCIYH FLVLSWYTFL NYYISQEGKD EVKPKILANG ARWKYMTLLN LLLQTIFYGV
TCLDDVLKRT KGGKDIKFLT AFRDLLFTTL AFPVSTFVFL AFWILFLYNR DLIYPKVLDT
VIPVWLNHAM HTFIFPITLA EVVLRPHSYP SKKTGLTLLA AASIAYISRI LWLYFETGTW
VYPVFAKLSL LGLAAFFSLS YVFIASIYLL GEKLNHWKWG DMRQPRKKRK
Sequence without tag. The proposed Strep-Tag is based on experience s with the expressior
system, a different complexity of the protein could make another tag necessary. In case you
have a special request, please contact us.
Key Benefits:
Made in Germany - from design to production - by highly experienced protein experts.
Protein expressed with ALiCE® and purified by multi-step, protein-specific process to ensure
correct folding and modification.These proteins are normally active (enzymatically functional) as our customers have

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/4 | Product datasheet for ABIN3105228 | 04/30/2024 | Copyright antibodies-online. All rights reserved. reported (not tested by us and not guaranteed).

• State-of-the-art algorithm used for plasmid design (Gene synthesis).

This protein is a **made-to-order protein** and will be made for the first time for your order. Our experts in the lab will ensure that you receive a correctly folded protein.

The big advantage of ordering our **made-to-order proteins** in comparison to ordering custom made proteins from other companies is that there is no financial obligation in case the protein cannot be expressed or purified.

Expression System:

• ALICE®, our Almost Living Cell-Free Expression System is based on a lysate obtained from Nicotiana tabacum c.v.. This contains all the protein expression machinery needed to produce even the most difficult-to-express proteins, including those that require post-translational modifications.

During lysate production, the cell wall and other cellular components that are not required for
protein production are removed, leaving only the protein production machinery and the
mitochondria to drive the reaction. During our lysate completion steps, the additional
components needed for protein production (amino acids, cofactors, etc.) are added to
produce something that functions like a cell, but without the constraints of a living system all that's needed is the DNA that codes for the desired protein!

Concentration:

- The concentration of our recombinant proteins is measured using the absorbance at 280nm.
- The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.
- We use the Expasy's ProtParam tool to determine the absorption coefficient of each protein.

Purification:	Two step purification of proteins expressed in Almost Living Cell-Free Expression System
	(ALICE®):
	1. In a first purification step, the protein is purified from the cleared cell lysate using StrepTag
	capture material. Eluate fractions are analyzed by SDS-PAGE.
	2. Protein containing fractions of the best purification are subjected to second purification step
	through size exclusion chromatography. Eluate fractions are analyzed by SDS-PAGE and
	Western blot.
Purity:	>80 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.
Endotoxin Level:	Low Endotoxin less than 1 EU/mg (< 0.1 ng/mg)

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 2/4 | Product datasheet for ABIN3105228 | 04/30/2024 | Copyright antibodies-online. All rights reserved.

Product Details

Grade:

Crystallography grade

Target Details

Target:	ADTRP
Alternative Name:	ADTRP (ADTRP Products)
Background:	Androgen-dependent TFPI-regulating protein (Fatty acid esters of hydroxy fatty acids hydrolase ADTRP) (FAHFA hydrolase ADTRP) (EC 3.1),FUNCTION: Hydrolyzes bioactive fatty-acid esters of hydroxy-fatty acids (FAHFAs), but not other major classes of lipids (PubMed:27018888). Show a preference for FAHFAs with branching distal from the carboxylate head group of the lipids (PubMed:27018888). Regulates the expression and the cell-associated anticoagulant activity of the inhibitor TFPI in endothelial cells (in vitro) (PubMed:21868574). {ECO:0000269 PubMed:21868574, ECO:0000269 PubMed:27018888}.
Molecular Weight:	26.8 kDa
UniProt:	Q96IZ2
Application Details	
Application Notes:	In addition to the applications listed above we expect the protein to work for functional studies as well. As the protein has not been tested for functional studies yet we cannot offer a guarantee though.
Comment:	 ALICE®, our Almost Living Cell-Free Expression System is based on a lysate obtained from Nicotiana tabacum c.v This contains all the protein expression machinery needed to produce even the most difficult-to-express proteins, including those that require post-translational modifications. During lysate production, the cell wall and other cellular components that are not required for protein production are removed, leaving only the protein production machinery and the mitochondria to drive the reaction. During our lysate completion steps, the additional components needed for protein production (amino acids, cofactors, etc.) are added to produce something that functions like a cell, but without the constraints of a living system - all that's needed is the DNA that codes for the desired protein!
Restrictions:	For Research Use only

Handling

Format:	Liquid
Buffer:	The buffer composition is at the discretion of the manufacturer. If you have a special request, please contact us.
Handling Advice:	Avoid repeated freeze-thaw cycles.
Storage:	-80 °C
Storage Comment:	Store at -80°C.
Expiry Date:	Unlimited (if stored properly)