antibodies

Datasheet for ABIN3111943 PLAP Protein (AA 23-506) (rho-1D4 tag)



Overview

Image

Quantity:	1 mg
Target:	PLAP (ALPP)
Protein Characteristics:	AA 23-506
Origin:	Human
Source:	Insect Cells
Protein Type:	Recombinant
Purification tag / Conjugate:	This PLAP protein is labelled with rho-1D4 tag.
Application:	Western Blotting (WB), SDS-PAGE (SDS), ELISA, Crystallization (Crys)

Product Details

Sequence:	IIPVEEENPD FWNREAAEAL GAAKKLQPAQ TAAKNLIIFL GDGMGVSTVT AARILKGQKK
	DKLGPEIPLA MDRFPYVALS KTYNVDKHVP DSGATATAYL CGVKGNFQTI GLSAAARFNQ
	CNTTRGNEVI SVMNRAKKAG KSVGVVTTTR VQHASPAGTY AHTVNRNWYS DADVPASARQ
	EGCQDIATQL ISNMDIDVIL GGGRKYMFRM GTPDPEYPDD YSQGGTRLDG KNLVQEWLAK
	RQGARYVWNR TELMQASLDP SVTHLMGLFE PGDMKYEIHR DSTLDPSLME MTEAALRLLS
	RNPRGFFLFV EGGRIDHGHH ESRAYRALTE TIMFDDAIER AGQLTSEEDT LSLVTADHSH
	VFSFGGYPLR GSSIFGLAPG KARDRKAYTV LLYGNGPGYV LKDGARPDVT ESESGSPEYR
	QQSAVPLDEE THAGEDVAVF ARGPQAHLVH GVQEQTFIAH VMAFAACLEP YTACDLAPPA GTTD
	Sequence without tag. Tag location is at the discretion of the manufacturer. If you have a
	special request, please contact us.
Characteristics:	 Made in Germany - from design to production - by highly experienced protein experts. Human ALPP Protein (raised in Insect Cells) purified by multi-step, protein-specific process to

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/4 | Product datasheet for ABIN3111943 | 09/11/2023 | Copyright antibodies-online. All rights reserved.

 Three step purification of membrane proteins expressed in baculovirus infected SF9 insect cells: 1. Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized with different detergents (detergent screen). Samples are analyzed by Western blot. 2. The best performing detergent is used for solubilization and the proteins are purified via their rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate fractions are analyzed by Western blot. 3. Protein containing fractions of the best purification are subjected to second purification step through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and Western blot. >95 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.
 Three step purification of membrane proteins expressed in baculovirus infected SF9 insect cells: 1. Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized with different detergents (detergent screen). Samples are analyzed by Western blot. 2. The best performing detergent is used for solubilization and the proteins are purified via their rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate fractions are analyzed by Western blot. 3. Protein containing fractions of the best purification are subjected to second purification step through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and
 Three step purification of membrane proteins expressed in baculovirus infected SF9 insect cells: 1. Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized with different detergents (detergent screen). Samples are analyzed by Western blot. 2. The best performing detergent is used for solubilization and the proteins are purified via their rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate
 Three step purification of membrane proteins expressed in baculovirus infected SF9 insect cells: 1. Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized with different detergents (detergent screen). Samples are analyzed by Western blot. 2. The best performing detergent is used for solubilization and the proteins are purified via their
Three step purification of membrane proteins expressed in baculovirus infected SF9 insect cells: 1. Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized with
Three step purification of membrane proteins expressed in baculovirus infected SF9 insect cells:
Three step purification of membrane proteins expressed in baculovirus infected SF9 insect
the Expasy's protparam tool to determine the absorption coefficient of each protein.
The concentration of the protein is calculated using its specific absorption coefficient. We use
specific reference buffer.
The protein's absorbance will be measured in several dilutions and is measured against its
The concentration of our recombinant proteins is measured using the absorbance at 280nm.
protein experts will do everything to make sure that you receive the protein you ordered.
folded protein. With no financial risk on your end you can rest assured that our experienced
When you order this made-to-order protein you will only pay upon receival of the correctly
experiments or purification optimization).
custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression
(other companies might charge you for any performed steps in the expression process for
In the unlikely event that the protein cannot be expressed or purified we do not charge anything
cannot be expressed or purified.
made proteins from other companies is that there is no financial obligation in case the protein
The big advantage of ordering our made-to-order proteins in comparison to ordering custom
This protein is a made to order protein and will be made for the first time for your order. Our experts in the lab will ensure that you receive a correctly folded protein.
ensure crystallization grade.State-of-the-art algorithm used for plasmid design (Gene synthesis).

Grade:

Endotoxin Level:

Product Details

Crystallography grade

Protein is endotoxin-free.

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 2/4 | Product datasheet for ABIN3111943 | 09/11/2023 | Copyright antibodies-online. All rights reserved.

Target Details	
Target:	PLAP (ALPP)
Alternative Name:	ALPP (ALPP Products)
Molecular Weight:	53.9 kDa Including tag.
UniProt:	P05187
Application Details	
Application Notes:	In addition to the applications listed above we expect the protein to work for functional studies as well. As the protein has not been tested for functional studies yet we cannot offer a guarantee though.
Comment:	In cases in which it is highly likely that the recombinant protein with the default tag will be insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible options with you in detail to assure that you receive your protein of interest.
Restrictions:	For Research Use only
Handling	
Format:	Liquid
Buffer:	100 mM NaCL, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer.
Handling Advice:	Avoid repeated freeze-thaw cycles.
Storage:	-80 °C
Storage Comment:	Store at -80°C.
Expiry Date:	Unlimited (if stored properly)



Image 1. "Crystallography Grade" protein due to multi-step, protein-specific purification process

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 4/4 | Product datasheet for ABIN3111943 | 09/11/2023 | Copyright antibodies-online. All rights reserved.