

# Datasheet for ABIN3113356 CD6 Protein (CD6) (AA 18-668) (rho-1D4 tag)



Overview

Quantity:	1 mg
Target:	CD6
Protein Characteristics:	AA 18-668
Origin:	Human
Source:	Insect Cells
Protein Type:	Recombinant
Purification tag / Conjugate:	This CD6 protein is labelled with rho-1D4 tag.
Application:	Western Blotting (WB), SDS-PAGE (SDS), ELISA, Crystallization (Crys)

### Product Details

Sequence:	HPSPAPPDQL NTSSAESELW EPGERLPVRL TNGSSSCSGT VEVRLEASWE PACGALWDSR
	AAEAVCRALG CGGAEAASQL APPTPELPPP PAAGNTSVAA NATLAGAPAL LCSGAEWRLC
	EVVEHACRSD GRRARVTCAE NRALRLVDGG GACAGRVEML EHGEWGSVCD DTWDLEDAHV
	VCRQLGCGWA VQALPGLHFT PGRGPIHRDQ VNCSGAEAYL WDCPGLPGQH YCGHKEDAGA
	VCSEHQSWRL TGGADRCEGQ VEVHFRGVWN TVCDSEWYPS EAKVLCQSLG CGTAVERPKG
	LPHSLSGRMY YSCNGEELTL SNCSWRFNNS NLCSQSLAAR VLCSASRSLH NLSTPEVPAS
	VQTVTIESSV TVKIENKESR ELMLLIPSIV LGILLLGSLI FIAFILLRIK GKYALPVMVN HQHLPTTIPA
	GSNSYQPVPI TIPKEVFMLP IQVQAPPPED SDSGSDSDYE HYDFSAQPPV ALTTFYNSQR
	HRVTDEEVQQ SRFQMPPLEE GLEELHASHI PTANPGHCIT DPPSLGPQYH PRSNSESSTS
	SGEDYCNSPK SKLPPWNPQV FSSERSSFLE QPPNLELAGT QPAFSAGPPA DDSSSTSSGE
	WYQNFQPPPQ PPSEEQFGCP GSPSPQPDST DNDDYDDISA A
	Sequence without tag. Tag location is at the discretion of the manufacturer. If you have a

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/4 | Product datasheet for ABIN3113356 | 08/02/2024 | Copyright antibodies-online. All rights reserved.

### Product Details

<ul> <li>ensure crystallization grade.</li> <li>State-of-the-art algorithm used for plasmid design (Gene synthesis).</li> <li>This protein is a made to order protein and will be made for the first time for your order. Our experts in the lab will ensure that you receive a correctly folded protein.</li> <li>The big advantage of ordering our made-to-order proteins in comparison to ordering custon made proteins from other companies is that there is no financial obligation in case the prote cannot be expressed or purified.</li> <li>In the unlikely event that the protein cannot be expression plasmid, the first expression receives a custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression experiments or purification optimization).</li> <li>When you order this made-to-order protein you will only pay upon receival of the correctly folded protein. With no financial risk on your end you can rest assured that our experienced protein experts will do everything to make sure that you receive the protein you ordered. The concentration of our recombinant proteins is measured using the absorbance at 280nn The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.</li> <li>The concentration of the protein is calculated using its specific absorption coefficient. We u the Expasy's protparam tool to determine the absorption coefficient of each protein.</li> <li>Purification:</li> <li>Three step purification of membrane proteins expressed in baculovirus infected SF9 insect cells:</li> <li>Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized wild different detergents (detergent screen). Samples are analyzed by Western blot.</li> <li>Protein containing fractions of the best purification are subjected to second purified via t through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and</li> </ul>		special request, please contact us.
experts in the lab will ensure that you receive a correctly folded protein.         The big advantage of ordering our made-to-order proteins in comparison to ordering custon made proteins from other companies is that there is no financial obligation in case the protein cannot be expressed or purified.         In the unlikely event that the protein cannot be expressed or purified we do not charge anyth (other companies might charge you for any performed steps in the expression process for custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression experiments or purification optimization).         When you order this made-to-order protein you will only pay upon receival of the correctly folded protein. With no financial risk on your end you can rest assured that our experienced protein experts will do everything to make sure that you receive the protein you ordered.         The concentration of our recombinant proteins is measured using the absorbance at 280nn The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.         The concentration of the protein is calculated using its specific absorption coefficient. We u the Expasy's protparam tool to determine the absorption coefficient of each protein.         Purification:       Three step purification of membrane proteins expressed in baculovirus infected SF9 insect cells:         1. Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized with different detergents (detergent screen). Samples are analyzed by Western blot.         2. The best performing detergent is used for solubilization and the proteins are purified via t no1D4 antibody columns: one DTT resistant, the other one not. Eluate fractio	Characteristics:	Human CD6 Protein (raised in Insect Cells) purified by multi-step, protein-specific process to ensure crystallization grade.
The big advantage of ordering our made-to-order proteins in comparison to ordering custon made proteins from other companies is that there is no financial obligation in case the protein cannot be expressed or purified.         In the unlikely event that the protein cannot be expressed or purified we do not charge anyth (other companies might charge you for any performed steps in the expression process for custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression experiments or purification optimization).         When you order this made-to-order protein you will only pay upon receival of the correctly folded protein. With no financial risk on your end you can rest assured that our experienced protein experts will do everything to make sure that you receive the protein you ordered. The concentration of our recombinant proteins is measured using the absorbance at 280nn The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.         Purification:       Three step purification of the protein is calculated using its specific absorption coefficient. We u the Expasy's protparam tool to determine the absorption coefficient of each protein.         Purification:       Three step purification of membrane proteins expressed in baculovirus infected SF9 insect cells:         1. Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized wit different detergents (detergent screen). Samples are analyzed by Western blot.         2. The best performing detergent is used for solubilization and the proteins are purified via t rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate fractions are analyzed by WD5-PAGE and through size exclusion chromatograph. Eluate fractions are		This protein is a made to order protein and will be made for the first time for your order. Our
made proteins from other companies is that there is no financial obligation in case the protein cannot be expressed or purified.         In the unlikely event that the protein cannot be expressed or purified we do not charge anyth (other companies might charge you for any performed steps in the expression process for custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression experiments or purification optimization).         When you order this made-to-order protein you will only pay upon receival of the correctly folded protein. With no financial risk on your end you can rest assured that our experienced protein experts will do everything to make sure that you receive the protein you ordered.         The concentration of our recombinant proteins is measured using the absorbance at 280nn The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.         The concentration of the protein is calculated using its specific absorption coefficient. We use the Expasy's protparam tool to determine the absorption coefficient of each protein.         Purification:       Three step purification of membrane proteins expressed in baculovirus infected SF9 insect cells:         1. Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized wild different detergents (detergent screen). Samples are analyzed by Western blot.         2. The best performing detergent is used for solubilization and the proteins are purified via t rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate fractions are analyzed by Western blot.         3. Protein containing fractions of the best purification are subjected to second purification s through size ex		experts in the lab will ensure that you receive a correctly folded protein.
<ul> <li>cannot be expressed or purified.</li> <li>In the unlikely event that the protein cannot be expressed or purified we do not charge anyth (other companies might charge you for any performed steps in the expression process for custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression experiments or purification optimization).</li> <li>When you order this made-to-order protein you will only pay upon receival of the correctly folded protein. With no financial risk on your end you can rest assured that our experienced protein experts will do everything to make sure that you receive the protein you ordered.</li> <li>The concentration of our recombinant proteins is measured using the absorbance at 280nn The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.</li> <li>The concentration of the protein is calculated using its specific absorption coefficient. We u the Expasy's protparam tool to determine the absorption coefficient of each protein.</li> <li>Purification:</li> <li>Thre estep purification of membrane proteins expressed in baculovirus infected SF9 insect cells:         <ul> <li>Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized wild different detergents (detergent screen). Samples are analyzed by Western blot.</li> <li>The best performing detergent is used for solubilization and the proteins are purified via t rho1D4 any via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate fractions are analyzed by Western blot.</li> <li>Protein containing fractions of the best purification are subjected to second purification s through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and</li> </ul> </li></ul>		The big advantage of ordering our made-to-order proteins in comparison to ordering custom
In the unlikely event that the protein cannot be expressed or purified we do not charge anyth (other companies might charge you for any performed steps in the expression process for custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression experiments or purification optimization). When you order this made-to-order protein you will only pay upon receival of the correctly folded protein. With no financial risk on your end you can rest assured that our experienced protein experts will do everything to make sure that you receive the protein you ordered. The concentration of our recombinant proteins is measured using the absorbance at 280nn The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer. The concentration of the protein is calculated using its specific absorption coefficient. We u the Expasy's protparam tool to determine the absorption coefficient of each protein. Purification: Three step purification of membrane proteins expressed in baculovirus infected SF9 insect cells: 1. Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized wit different detergents (detergent screen). Samples are analyzed by Western blot. 2. The best performing detergent is used for solubilization and the proteins are purified via t rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate fractions are analyzed by Western blot. 3. Protein containing fractions of the best purification are subjected to second purification s through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and		made proteins from other companies is that there is no financial obligation in case the protein
(other companies might charge you for any performed steps in the expression process for custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression experiments or purification optimization).         When you order this made-to-order protein you will only pay upon receival of the correctly folded protein. With no financial risk on your end you can rest assured that our experienced protein experts will do everything to make sure that you receive the protein you ordered. The concentration of our recombinant proteins is measured using the absorbance at 280nn The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.         The concentration of the protein is calculated using its specific absorption coefficient. We u the Expasy's protparam tool to determine the absorption coefficient of each protein.         Purification:       Three step purification of membrane proteins expressed in baculovirus infected SF9 insect cells:         1. Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized with different detergents (detergent screen). Samples are analyzed by Western blot.         2. The best performing detergent is used for solubilization and the proteins are purified via t no1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate fractions are analyzed by Western blot.         3. Protein containing fractions of the best purification are subjected to second purification s through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and		cannot be expressed or purified.
custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression experiments or purification optimization).         When you order this made-to-order protein you will only pay upon receival of the correctly folded protein. With no financial risk on your end you can rest assured that our experienced protein experts will do everything to make sure that you receive the protein you ordered.         The concentration of our recombinant proteins is measured using the absorbance at 280nn The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.         The concentration of the protein is calculated using its specific absorption coefficient. We u the Expasy's protparam tool to determine the absorption coefficient of each protein.         Purification:       Three step purification of membrane proteins expressed in baculovirus infected SF9 insect cells:         1. Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized with different detergents (detergent screen). Samples are analyzed by Western blot.         2. The best performing detergent is used for solubilization and the proteins are purified via t rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate fractions are analyzed by Western blot.         3. Protein containing fractions of the best purification are subjected to second purification s through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and		In the unlikely event that the protein cannot be expressed or purified we do not charge anything
<ul> <li>experiments or purification optimization).</li> <li>When you order this made-to-order protein you will only pay upon receival of the correctly folded protein. With no financial risk on your end you can rest assured that our experienced protein experts will do everything to make sure that you receive the protein you ordered. The concentration of our recombinant proteins is measured using the absorbance at 280nn The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.</li> <li>The concentration of the protein is calculated using its specific absorption coefficient. We use the Expasy's protparam tool to determine the absorption coefficient of each protein.</li> <li>Purification: Three step purification of membrane proteins expressed in baculovirus infected SF9 insect cells:         <ul> <li>Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized with different detergents (detergent screen). Samples are analyzed by Western blot.</li> <li>The best performing detergent is used for solubilization and the proteins are purified via t rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate fractions are analyzed by Western blot.</li> <li>Protein containing fractions of the best purification are subjected to second purification s through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and</li> </ul> </li> </ul>		(other companies might charge you for any performed steps in the expression process for
<ul> <li>When you order this made-to-order protein you will only pay upon receival of the correctly folded protein. With no financial risk on your end you can rest assured that our experienced protein experts will do everything to make sure that you receive the protein you ordered. The concentration of our recombinant proteins is measured using the absorbance at 280nn The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer. The concentration of the protein is calculated using its specific absorption coefficient. We ut the Expasy's protparam tool to determine the absorption coefficient of each protein.</li> <li>Purification: Three step purification of membrane proteins expressed in baculovirus infected SF9 insect cells:         <ol> <li>Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized with different detergents (detergent screen). Samples are analyzed by Western blot.</li> <li>The best performing detergent is used for solubilization and the proteins are purified via t rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate fractions are analyzed by Western blot.</li> </ol> </li> </ul>		custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression
folded protein. With no financial risk on your end you can rest assured that our experienced protein experts will do everything to make sure that you receive the protein you ordered. The concentration of our recombinant proteins is measured using the absorbance at 280nm The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.         The concentration of the protein is calculated using its specific absorption coefficient. We u the Expasy's protparam tool to determine the absorption coefficient of each protein.         Purification:       Three step purification of membrane proteins expressed in baculovirus infected SF9 insect cells:         1. Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized with different detergents (detergent screen). Samples are analyzed by Western blot.         2. The best performing detergent is used for solubilization and the proteins are purified via t rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate fractions are analyzed by Western blot.         3. Protein containing fractions of the best purification are subjected to second purification s through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and		experiments or purification optimization).
<ul> <li>protein experts will do everything to make sure that you receive the protein you ordered. The concentration of our recombinant proteins is measured using the absorbance at 280nm The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer. The concentration of the protein is calculated using its specific absorption coefficient. We us the Expasy's protparam tool to determine the absorption coefficient of each protein.</li> <li>Purification: Three step purification of membrane proteins expressed in baculovirus infected SF9 insect cells:         <ol> <li>Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized with different detergents (detergent screen). Samples are analyzed by Western blot.</li> <li>The best performing detergent is used for solubilization and the proteins are purified via t rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate fractions are analyzed by Western blot.</li> <li>Protein containing fractions of the best purification are subjected to second purification s through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and</li> </ol> </li> </ul>		When you order this made-to-order protein you will only pay upon receival of the correctly
The concentration of our recombinant proteins is measured using the absorbance at 280nm         The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.         The concentration of the protein is calculated using its specific absorption coefficient. We use the Expasy's protparam tool to determine the absorption coefficient of each protein.         Purification:       Three step purification of membrane proteins expressed in baculovirus infected SF9 insect cells:         1. Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized with different detergents (detergent screen). Samples are analyzed by Western blot.         2. The best performing detergent is used for solubilization and the proteins are purified via t rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate fractions are analyzed by Western blot.         3. Protein containing fractions of the best purification are subjected to second purification s through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and		folded protein. With no financial risk on your end you can rest assured that our experienced
The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.         The concentration of the protein is calculated using its specific absorption coefficient. We use the Expasy's protparam tool to determine the absorption coefficient of each protein.         Purification:       Three step purification of membrane proteins expressed in baculovirus infected SF9 insect cells:         1. Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized with different detergents (detergent screen). Samples are analyzed by Western blot.         2. The best performing detergent is used for solubilization and the proteins are purified via t rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate fractions are analyzed by Western blot.         3. Protein containing fractions of the best purification are subjected to second purification s through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and		protein experts will do everything to make sure that you receive the protein you ordered.
<ul> <li>specific reference buffer.</li> <li>The concentration of the protein is calculated using its specific absorption coefficient. We use the Expasy's protparam tool to determine the absorption coefficient of each protein.</li> <li>Purification:</li> <li>Three step purification of membrane proteins expressed in baculovirus infected SF9 insect cells:         <ol> <li>Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized with different detergents (detergent screen). Samples are analyzed by Western blot.</li> <li>The best performing detergent is used for solubilization and the proteins are purified via t rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate fractions are analyzed by Western blot.</li> </ol> </li> <li>Protein containing fractions of the best purification are subjected to second purification s through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and</li> </ul>		The concentration of our recombinant proteins is measured using the absorbance at 280nm.
The concentration of the protein is calculated using its specific absorption coefficient. We use the Expasy's protparam tool to determine the absorption coefficient of each protein.         Purification:       Three step purification of membrane proteins expressed in baculovirus infected SF9 insect cells:         1. Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized with different detergents (detergent screen). Samples are analyzed by Western blot.         2. The best performing detergent is used for solubilization and the proteins are purified via t rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate fractions are analyzed by Western blot.         3. Protein containing fractions of the best purification are subjected to second purification set through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and		The protein's absorbance will be measured in several dilutions and is measured against its
the Expasy's protparam tool to determine the absorption coefficient of each protein.         Purification:       Three step purification of membrane proteins expressed in baculovirus infected SF9 insect cells:         1. Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized with different detergents (detergent screen). Samples are analyzed by Western blot.         2. The best performing detergent is used for solubilization and the proteins are purified via t rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate fractions are analyzed by Western blot.         3. Protein containing fractions of the best purification are subjected to second purification s through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and		specific reference buffer.
<ul> <li>Purification: Three step purification of membrane proteins expressed in baculovirus infected SF9 insect cells:</li> <li>1. Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized with different detergents (detergent screen). Samples are analyzed by Western blot.</li> <li>2. The best performing detergent is used for solubilization and the proteins are purified via t rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate fractions are analyzed by Western blot.</li> <li>3. Protein containing fractions of the best purification are subjected to second purification s through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and</li> </ul>		The concentration of the protein is calculated using its specific absorption coefficient. We use
<ul> <li>cells:</li> <li>1. Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized with different detergents (detergent screen). Samples are analyzed by Western blot.</li> <li>2. The best performing detergent is used for solubilization and the proteins are purified via the rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate fractions are analyzed by Western blot.</li> <li>3. Protein containing fractions of the best purification are subjected to second purification s through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and</li> </ul>		the Expasy's protparam tool to determine the absorption coefficient of each protein.
<ol> <li>Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized with different detergents (detergent screen). Samples are analyzed by Western blot.</li> <li>The best performing detergent is used for solubilization and the proteins are purified via the rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate fractions are analyzed by Western blot.</li> <li>Protein containing fractions of the best purification are subjected to second purification s through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and</li> </ol>	Purification:	Three step purification of membrane proteins expressed in baculovirus infected SF9 insect
<ul> <li>different detergents (detergent screen). Samples are analyzed by Western blot.</li> <li>2. The best performing detergent is used for solubilization and the proteins are purified via t rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate fractions are analyzed by Western blot.</li> <li>3. Protein containing fractions of the best purification are subjected to second purification s through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and</li> </ul>		cells:
<ol> <li>The best performing detergent is used for solubilization and the proteins are purified via t rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate fractions are analyzed by Western blot.</li> <li>Protein containing fractions of the best purification are subjected to second purification s through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and</li> </ol>		1. Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized with
<ul> <li>rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate fractions are analyzed by Western blot.</li> <li>3. Protein containing fractions of the best purification are subjected to second purification s through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and</li> </ul>		
3. Protein containing fractions of the best purification are subjected to second purification s through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and		rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate
through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and		fractions are analyzed by Western blot.
		3. Protein containing fractions of the best purification are subjected to second purification step
		through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and Western blot.
western blot.		Western blot.
Purity: >95 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.	Purity:	>95 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.
Sterility: 0.22 µm filtered	Sterility:	0.22 µm filtered
Endotoxin Level: Protein is endotoxin-free.	Endotoxin Level:	Protein is endotoxin-free.

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 2/4 | Product datasheet for ABIN3113356 | 08/02/2024 | Copyright antibodies-online. All rights reserved.

#### Product Details

Grade:

Crystallography grade

## Target Details

Target:	CD6
Alternative Name:	CD6 (CD6 Products)
Background:	Cell adhesion molecule that mediates cell-cell contacts and regulates T-cell responses via its
	interaction with ALCAM/CD166 (PubMed:15048703, PubMed:15294938, PubMed:16352806,
	PubMed:16914752, PubMed:24945728, PubMed:24584089). Contributes to signaling cascades
	triggered by activation of the TCR/CD3 complex (PubMed:24584089). Functions as
	costimulatory molecule, promotes T-cell activation and proliferation (PubMed:15294938,
	PubMed:16352806, PubMed:16914752). Contributes to the formation and maturation of the
	immunological synapse (PubMed:15294938, PubMed:16352806). Functions as calcium-
	dependent pattern receptor that binds and aggregates both Gram-positive and Gram-negative
	bacteria. Binds both lipopolysaccharide (LPS) from Gram-negative bacteria and lipoteichoic
	acid from Gram-positive bacteria (PubMed:17601777). LPS binding leads to the activation of
	signaling cascades and down-stream MAP kinases (PubMed:17601777). Mediates activation of
	the inflammatory response and the secretion of pro-inflammatory cytokines in response to LPS
	(PubMed:17601777). {ECO:0000269 PubMed:15048703, ECO:0000269 PubMed:15294938,
	EC0:0000269 PubMed:16352806, EC0:0000269 PubMed:16914752,
	ECO:0000269 PubMed:17601777, ECO:0000269 PubMed:24584089,
	ECO:0000269 PubMed:24945728}.
Molecular Weight:	71.2 kDa Including tag.
UniProt:	P30203
Application Details	
Application Notes:	In addition to the applications listed above we expect the protein to work for functional studies
	as well. As the protein has not been tested for functional studies yet we cannot offer a gurantee
	though.
Comment:	In cases in which it is highly likely that the recombinant protein with the default tag will be
	insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to
	increase solubility. We will discuss all possible options with you in detail to assure that you
	receive your protein of interest.
Restrictions:	For Research Use only
	-

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 3/4 | Product datasheet for ABIN3113356 | 08/02/2024 | Copyright antibodies-online. All rights reserved.

## Handling

Format:	Liquid
Buffer:	100 mM NaCL, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer.
Handling Advice:	Avoid repeated freeze-thaw cycles.
Storage:	-80 °C
Storage Comment:	Store at -80°C.
Expiry Date:	Unlimited (if stored properly)