

Datasheet for ABIN3114048

## HERV-FRD Provirus Ancestral Env Polyprotein (Herv-frd) (AA 16-538) protein (rho-1D4 tag)



[Go to Product page](#)

### 1 Image

#### Overview

Quantity:	1 mg
Target:	HERV-FRD Provirus Ancestral Env Polyprotein (Herv-frd)
Protein Characteristics:	AA 16-538
Origin:	Human
Source:	Insect Cells
Protein Type:	Recombinant
Purification tag / Conjugate:	rho-1D4 tag
Application:	ELISA, Western Blotting (WB), Crystallization (Crys), SDS-PAGE (SDS)

#### Product Details

Sequence: AYRHPDFPLL EKAQQLQST GSPYSTNCWL CTSSSTETPG TAYPASPREW TSIEAELHIS  
YRWDPNLKGL MRPANLLST VKQDFDIRQ KPPIFGPIFT NINLMGIAP I CVMAKRKNGT  
NVGTL PSTVC NVTFTVDSNQ QTYQTYTHNQ FRHQPRFPKP PNITFPQGLT LDKSSRFCQG  
RPSSCSTRNF WFRPADYNQC LQISNLSSTA EWVLLDQTRN SLFWENKTKG ANQSQTPCVQ  
VLAGMTIATS YLGISAVSEF FGTSLTPLFH FHISTCLKTQ GAFYICGQSI HQCLPSNWTG  
TCTIGYVTPD IFIAPGNLSL PIPIYGN SPL PRVRRAIHFI PLLAGLGILA GTGTGIAGIT KASLTYSQLS  
KEIANNIDTM AKALTTM QEQ IDSLA AVLQ NRRGLDMLTA AQGGICLALD EKCCFWVNQS  
GKVQDNIRQL LNQASSLRER ATQGWLNWEG TWKWFSWVLP LTGPLVSLLL LLLFGPCLLN  
LITQFVSSRL QAIKLQTNLS AGRHPRNIQE SPF

**Sequence without tag. Tag location is at the discretion of the manufacturer. If you have a special request, please contact us.**

## Product Details

---

### Characteristics:

- Made in Germany - from design to production - by highly experienced protein experts.
- Human ERVFRD-1 Protein (raised in Insect Cells) purified by multi-step, protein-specific process to ensure crystallization grade.
- State-of-the-art algorithm used for plasmid design (Gene synthesis).

This protein is a made to order protein and will be made for the first time for your order. Our experts in the lab will ensure that you receive a correctly folded protein.

The big advantage of ordering our made-to-order proteins in comparison to ordering custom made proteins from other companies is that there is no financial obligation in case the protein cannot be expressed or purified.

In the unlikely event that the protein cannot be expressed or purified we do not charge anything (other companies might charge you for any performed steps in the expression process for custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression experiments or purification optimization).

When you order this made-to-order protein you will only pay upon receipt of the correctly folded protein. With no financial risk on your end you can rest assured that our experienced protein experts will do everything to make sure that you receive the protein you ordered.

The concentration of our recombinant proteins is measured using the absorbance at 280nm. The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.

The concentration of the protein is calculated using its specific absorption coefficient. We use the Expasy's protparam tool to determine the absorption coefficient of each protein.

---

### Purification:

Three step purification of membrane proteins expressed in baculovirus infected SF9 insect cells:

1. Membrane proteins are fractionated by ultracentrifugation and subsequently solubilized with different detergents (detergent screen). Samples are analyzed by Western blot.
2. The best performing detergent is used for solubilization and the proteins are purified via their rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate fractions are analyzed by Western blot.
3. Protein containing fractions of the best purification are subjected to second purification step through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and Western blot.

---

### Purity:

>95 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.

---

### Sterility:

0.22 µm filtered

---

### Endotoxin Level:

Protein is endotoxin-free.

---

### Grade:

Crystallography grade

## Target Details

---

Target: HERV-FRD Provirus Ancestral Env Polyprotein (Herv-frd)

Alternative Name: ERVFRD-1 ([Herv-frd Products](#))

Target Type: Viral Protein

Background: This endogenous retroviral envelope protein has retained its original fusogenic properties and participates in trophoblast fusion and the formation of a syncytium during placenta morphogenesis. The interaction with MFSD2A is apparently important for this process (PubMed:18988732). {ECO:0000269|PubMed:18988732}. Endogenous envelope proteins may have kept, lost or modified their original function during evolution but this one can still make pseudotypes with MLV, HIV-1 or SIV-1 virions and confer infectivity. Retroviral envelope proteins mediate receptor recognition and membrane fusion during early infection. The surface protein mediates receptor recognition, while the transmembrane protein anchors the envelope heterodimer to the viral membrane through one transmembrane domain. The other hydrophobic domain, called fusion peptide, mediates fusion of the viral membrane with the target cell membrane (PubMed:14694139). {ECO:0000269|PubMed:14694139}.

Molecular Weight: 59.1 kDa Including tag.

UniProt: [P60508](#)

## Application Details

---

Application Notes: In addition to the applications listed above we expect the protein to work for functional studies as well. As the protein has not been tested for functional studies yet we cannot offer a guarantee though.

Comment: In cases in which it is highly likely that the recombinant protein with the default tag will be insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible options with you in detail to assure that you receive your protein of interest.

Restrictions: For Research Use only

## Handling

---

Format: Liquid

Buffer: 100 mM NaCl, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer.

Handling Advice: Avoid repeated freeze-thaw cycles.

## Handling

---

Storage: -80 °C

Storage Comment: Store at -80°C.

Expiry Date: Unlimited (if stored properly)

## Images

---



**Image 1.** „Crystallography Grade“ protein due to multi-step, protein-specific purification process