antibodies

## Datasheet for ABIN3114189 FAM189B Protein (AA 1-668) (rho-1D4 tag)





Overview

| Quantity:                     | 1 mg   |
|-------------------------------|--|
| Target:                       | FAM189B  |
| Protein Characteristics:      | AA 1-668   |
| Origin:                       | Human  |
| Source:                       | Insect Cells   |
| Protein Type:                 | Recombinant  |
| Purification tag / Conjugate: | This FAM189B protein is labelled with rho-1D4 tag.                   |
| Application:                  | ELISA, Western Blotting (WB), Crystallization (Crys), SDS-PAGE (SDS) |

## Product Details

| Sequence: | MMPSPSDSSR SLTSRPSTRG LTHLRLHRPW LQALLTLGLV QVLLGILVVT FSMVASSVTT |
|-----------|---|
|           | TESIKRSCPS WAGFSLAFSG VVGIVSWKRP FTLVISFFSL LSVLCVMLSM AGSVLSCKNA |
|           | QLARDFQQCS LEGKVCVCCP SVPLLRPCPE SGQELKVAPN STCDEARGAL KNLLFSVCGL |
|           | TICAAIICTL SAIVCCIQIF SLDLVHTLAP ERSVSGPLGP LGCTSPPPAP LLHTMLDLEE |
|           | FVPPVPPPPY YPPEYTCSSE TDAQSITYNG SMDSPVPLYP TDCPPSYEAV MGLRGDSQAT |
|           | LFDPQLHDGS CICERVASIV DVSMDSGSLV LSAIGDLPGG SSPSEDSCLL ELQGSVRSVD |
|           | YVLFRSIQRS RAGYCLSLDC GLRGPFEESP LPRRPPRAAR SYSCSAPEAP PPLGAPTAAR |
|           | SCHRLEGWPP WVGPCFPELR RRVPRGGGRP AAAPPTRAPT RRFSDSSGSL TPPGHRPPHP |
|           | ASPPPLLLPR SHSDPGITTS SDTADFRDLY TKVLEEEAAS VSSADTGLCS EACLFRLARC |
|           | PSPKLLRARS AEKRRPVPTF QKVPLPSGPA PAHSLGDLKG SWPGRGLVTR FLQISRKAPD |
|           | PSGTGAHGHK QVPRSLWGRP GRESLHLRSC GDLSSSSSLR RLLSGRRLER GTRPHSLSLN |
|           | GGSRETGL  |

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/4 | Product datasheet for ABIN3114189 | 09/11/2023 | Copyright antibodies-online. All rights reserved.

| special request, please contact us.           Characteristics:         • Made in Germany - from design to production - by highly experienced protein experts.           • Human FAM1898 Protein (raised in Insect Cells) punfied by multi-step, protein-specific process to ensure crystallization grade.         • State-of-the-art algorithm used for plasmid design (Cene synthesis).           This protein is a made to order protein and will be made for the first time for your order. Our experts in the lab will ensure that you receive a correctly folded protein.         • The big advantage of ordering our made-to-order proteins in comparison to ordering custom made proteins from other companies is that there is no financial obligation in case the protein cannot be expressed or purified.           In the unlikely event that the protein cannot be expressed or purified we do not charge anything (other companies might charge you for any performed steps in the expression process for custom-made proteins, eq. fees might apply for the expression plasmid, the first expression experiments or purification optimization).           When you order this made-to-order protein you will only pay upon receival of the correctly folded protein. With no financial risk on your end you can rest assured that our experiments specific reference buffer.           The concentration of our recombinant proteins is measured using the absorbance at 280m.           The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.           The concentration of the protein is calculated using its specific absorption coefficient. We use the Expasy's protparam tool to determine the absorbance officient of each protein.           <   |                  | Sequence without tag. Tag location is at the discretion of the manufacturer. If you have a   |
|--|------------------|--|
| <ul> <li>Human FAM189B Protein (raised in Insect Cells) purified by multi-step, protein-specific process to ensure crystallization grade.</li> <li>State-of-the-art algorithm used for plasmid design (Gene synthesis).</li> <li>This protein is a made to order protein and will be made for the first time for your order. Our experts in the lab will ensure that you receive a correctly folded protein.</li> <li>The big advantage of ordering our made-to-order proteins in comparison to ordering custom made proteins from other companies is that there is no financial obligation in case the protein cannot be expressed or purified.</li> <li>In the unlikely event that the protein cannot be expressed or purified we do not charge anything (other companies might charge you for any performed steps in the expression process for custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression experiments or purification optimization).</li> <li>When you order this made-to-order protein you will only pay upon receival of the correctly folded protein. With no financial risk on your end you can rest assured that our experienced protein seperts will do everything to make sure that you receive the protein you ordered.</li> <li>The concentration of the protein is calculated using its specific absorption coefficient. We use the Expany's protparam tool to determine the absorption coefficient we use the Expany's protparam tool to determine the absorption coefficient of each protein.</li> <li>Purification:</li> <li>Three step purification of membrane proteins expressed in baculovirus infected SF9 insect cells:         <ul> <li>Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized with different detergents (detergent screen). Samples are analyzed by Western blot.</li> </ul> </li> <li>Purification: - The estep purification of rom theory of resisting the protein are purified via their rhatD4 are yoa wore not. Livate fractions</li></ul> |                  | special request, please contact us.  |
| experts in the lab will ensure that you receive a correctly folded protein.The big advantage of ordering our made-to-order proteins in companison to ordering custom<br>made proteins from other companies is that there is no financial obligation in case the protein<br>cannot be expressed or purified.In the unlikely event that the protein cannot be expressed or purified we do not charge anything<br>(other companies might charge you for any performed steps in the expression process for<br>custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression<br>experiments or purification optimization).When you order this made-to-order protein you will only pay upon receival of the correctly<br>folded protein. With no financial risk on your end you can rest assured that our experienced<br>protein experts will do everything to make sure that you receive the protein you ordered.<br>The concentration of our recombinant proteins is measured using the absorbance at 280nm.<br>The protein's absorbance will be measured in several dilutions and is measured against its<br>specific reference buffer.<br>The concentration of the protein is calculated using its specific absorption coefficient. We use<br>the Expasy's protparam tool to determine the absorption coefficient of each protein.Purification:Three step purification of membrane proteins expressed in baculovirus infected SF9 insect<br>cells:<br>1. Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized with<br>different detergents (detergent screen). Samples are analyzed by Western blot.2. The best performing detergent is used for solubilization and the proteins are purified via their<br>rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate<br>fractions are analyzed by Western blot.2. The best performing detergent biot.3. Protein con  | Characteristics: | Human FAM189B Protein (raised in Insect Cells) purified by multi-step, protein-specific process to ensure crystallization grade.   |
| made proteins from other companies is that there is no financial obligation in case the protein<br>cannot be expressed or purified.In the unlikely event that the protein cannot be expressed or purified we do not charge anything<br>(other companies might charge you for any performed steps in the expression process for<br>custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression<br>experiments or purification optimization).When you order this made-to-order protein you will only pay upon receival of the correctly<br>folded protein. With no financial risk on your end you can rest assured that our experienced<br>protein experts will do everything to make sure that you receive the protein you ordered.<br>The concentration of our recombinant proteins is measured using the absorbance at 280nm.<br>The protein's absorbance will be measured in several dilutions and is measured against its<br>specific reference buffer.<br>The concentration of the protein is calculated using its specific absorption coefficient. We use<br>the Expasy's protparam tool to determine the absorption coefficient of each protein.Purification:Three step purification of membrane proteins expressed in baculovirus infected SF9 insect<br>cells:<br>1. Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized with<br>different detergents (detergent screen). Samples are analyzed by Western blot.<br>2. The best performing detergent is used for solubilization and the proteins are purified via their<br>rho1D4 tag via two ho1D4 antibody columns: one DTT resistant, the other one not. Eluate<br>fractions are analyzed by Western blot.Purity:>95 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.  |                  |  |
| (other companies might charge you for any performed steps in the expression process for<br>custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression<br>experiments or purification optimization).When you order this made-to-order protein you will only pay upon receival of the correctly<br>folded protein. With no financial risk on your end you can rest assured that our experienced<br>protein experts will do everything to make sure that you receive the protein you ordered.<br>The concentration of our recombinant proteins is measured using the absorbance at 280nm.<br>The protein's absorbance will be measured in several dilutions and is measured against its<br>specific reference buffer.<br>The concentration of the protein is calculated using its specific absorption coefficient. We use<br>the Expasy's protparam tool to determine the absorption coefficient of each protein.Purification:Three step purification of membrane proteins expressed in baculovirus infected SF9 insect<br>cells:<br>1. Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized with<br>different detergents (detergent screen). Samples are analyzed by Western blot.<br>2. The best performing detergent is used for solubilization and the proteins are purified via their<br>rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate<br>fractions are analyzed by Western blot.Purity:>95 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.  |                  | made proteins from other companies is that there is no financial obligation in case the protein  |
| folded protein. With no financial risk on your end you can rest assured that our experienced protein experts will do everything to make sure that you receive the protein you ordered. The concentration of our recombinant proteins is measured using the absorbance at 280nm. The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.         The concentration of the protein is calculated using its specific absorption coefficient. We use the Expasy's protparam tool to determine the absorption coefficient of each protein.         Purification:       Three step purification of membrane proteins expressed in baculovirus infected SF9 insect cells: <ol> <li>Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized with different detergents (detergent screen). Samples are analyzed by Western blot.</li> <li>The best performing detergent is used for solubilization and the proteins are purified via their rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate fractions are analyzed by Western blot.</li> <li>Protein containing fractions of the best purification are subjected to second purification step through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and Western blot.</li> </ol>  |                  | (other companies might charge you for any performed steps in the expression process for custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression   |
| The protein's absorbance will be measured in several dilutions and is measured against its<br>specific reference buffer.<br>The concentration of the protein is calculated using its specific absorption coefficient. We use<br>the Expasy's protparam tool to determine the absorption coefficient of each protein.Purification:Three step purification of membrane proteins expressed in baculovirus infected SF9 insect<br>cells:<br>1. Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized with<br>different detergents (detergent screen). Samples are analyzed by Western blot.<br>2. The best performing detergent is used for solubilization and the proteins are purified via their<br>rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate<br>fractions are analyzed by Western blot.Purity:>95 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.  |                  | folded protein. With no financial risk on your end you can rest assured that our experienced protein experts will do everything to make sure that you receive the protein you ordered.   |
| the Expasy's protparam tool to determine the absorption coefficient of each protein.Purification:Three step purification of membrane proteins expressed in baculovirus infected SF9 insect<br>cells:<br>1. Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized with<br>different detergents (detergent screen). Samples are analyzed by Western blot.<br>2. The best performing detergent is used for solubilization and the proteins are purified via their<br>rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate<br>fractions are analyzed by Western blot.3. Protein containing fractions of the best purification are subjected to second purification step<br>through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and<br>Western blot.Purity:>95 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.  |                  | The protein's absorbance will be measured in several dilutions and is measured against its   |
| cells:1. Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized with<br>different detergents (detergent screen). Samples are analyzed by Western blot.2. The best performing detergent is used for solubilization and the proteins are purified via their<br>rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate<br>fractions are analyzed by Western blot.3. Protein containing fractions of the best purification are subjected to second purification step<br>through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and<br>Western blot.Purity:>95 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.  |                  |  |
| <ul> <li>different detergents (detergent screen). Samples are analyzed by Western blot.</li> <li>2. The best performing detergent is used for solubilization and the proteins are purified via their rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate fractions are analyzed by Western blot.</li> <li>3. Protein containing fractions of the best purification are subjected to second purification step through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and Western blot.</li> <li>Purity: &gt;95 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.</li> </ul>  | Purification:    |  |
|  |                  | <ul> <li>different detergents (detergent screen). Samples are analyzed by Western blot.</li> <li>2. The best performing detergent is used for solubilization and the proteins are purified via their rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate fractions are analyzed by Western blot.</li> <li>3. Protein containing fractions of the best purification are subjected to second purification step through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and</li> </ul> |
| Sterility: 0.22 µm filtered  | Purity:          | >95 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.   |
|  | Sterility:       | 0.22 µm filtered   |

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 2/4 | Product datasheet for ABIN3114189 | 09/11/2023 | Copyright antibodies-online. All rights reserved.

| Product Details     |  |
|---------------------|--|
| Endotoxin Level:    | Protein is endotoxin-free.   |
| Grade:              | Crystallography grade  |
| Target Details      |  |
| Target:             | FAM189B  |
| Alternative Name:   | FAM189B (FAM189B Products)   |
| Molecular Weight:   | 72.5 kDa Including tag.  |
| UniProt:            | P81408   |
| Application Details |  |
| Application Notes:  | In addition to the applications listed above we expect the protein to work for functional studies as well. As the protein has not been tested for functional studies yet we cannot offer a guarantee though.   |
| Comment:            | In cases in which it is highly likely that the recombinant protein with the default tag will be<br>insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to<br>increase solubility. We will discuss all possible options with you in detail to assure that you<br>receive your protein of interest. |
| Restrictions:       | For Research Use only  |
| Handling            |  |
| Format:             | Liquid   |
| Buffer:             | 100 mM NaCL, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer.   |
| Handling Advice:    | Avoid repeated freeze-thaw cycles.   |
| Storage:            | -80 °C   |
| Storage Comment:    | Store at -80°C.  |
| Expiry Date:        | Unlimited (if stored properly)   |

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 3/4 | Product datasheet for ABIN3114189 | 09/11/2023 | Copyright antibodies-online. All rights reserved.



**Image 1.** "Crystallography Grade" protein due to multi-step, protein-specific purification process

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 4/4 | Product datasheet for ABIN3114189 | 09/11/2023 | Copyright antibodies-online. All rights reserved.